

## Product Information

### SARS-CoV2 Main protease chromogenic substrate peptide

Catalog Number: **SAE0181**  
Storage Temperature  $-20\text{ }^{\circ}\text{C}$

Synonym: Thr-Ser-Ala-Val-Leu-Gln-pNA

#### Product Description

This hexapeptide is a substrate peptide for the Main protease (Mpro) from the SARS coronavirus<sup>1,2</sup> and also that of SARS-Cov2.<sup>3</sup> This substrate allows measurement of the activity of Mpro, using a standard plate reader or spectrophotometer. Proteolytic cleavage by Mpro will release the free pNA group from the peptide, with a maximum absorption peak at 405 nm ( $A_{405_{1\mu\text{M}}} = 0.00916$ ):



This product is lyophilized from 0.1% TFA in H<sub>2</sub>O.

#### Precautions and Disclaimer

For Research use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

#### Storage/Stability

The product is stable in lyophilized form for at least 5 years when stored at  $-20\text{ }^{\circ}\text{C}$ . This product has a limited lifetime in solution. Long-term storage in solution should be avoided.

#### Procedure

##### Solubilization:

The peptide is soluble in water at a concentration of 20 mg/mL, or in DMSO at a concentration of 20-50 mg/mL.

##### Working concentration:

Peptide: suggested working concentration of 200  $\mu\text{g/mL}$ .

Mpro: we recommend using Mpro (Catalog Number SAE0172) at a working concentration of 2-20  $\mu\text{g/mL}$ .

Assay buffer: Mpro protease is active under a wide variety of conditions and temperatures. However, we strongly advise against use of Tris-based buffers, because of their interactions with the given substrate.

Recommended assay buffer: 25 mM HEPES buffer, pH 7.0, 0.2% Tween-20.

Reaction temperature: 0 - 30  $^{\circ}\text{C}$

The reaction can be stopped by addition of acetic acid to a final concentration of 2%.

#### References

1. Liu, P. *et al.*, *Eur. J. Med. Chem.*, **206**, 112702 (2020).
2. Liu, Z. *et al.*, *J. Chem. Inf. Model.*, **45(1)**, 10-17 (2005).
3. Li, C. *et al.*, *J. Biol. Chem.*, **285(36)**, 28134-28140 (2010).

Tween is a registered trademark of Croda International PLC.

DT, EV, HE, LB, GCY, SBC-12/20-1