Instructions for Use

Eosin Y Solutions

Procedure No. HT110

Intended Use

Alcoholic and aqueous Eosin Y solutions are intended for use as general-purpose counterstains. Eosin Y Solutions are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual qualitative procedure is used for the detection of cytoplasm in routine histology or cytology human specimens. This data can be used as an aid to diagnosis to certain clinical conditions and pathophysiological states as Eosin Y stains cytoplasm of certain cells including muscle, collagen, and red blood cells. It should be reviewed in conjunction with other clinical diagnostic tests or information.

Eosin Y is commonly used as a cytoplasmic stain.¹ As early as 1885 it was advocated by List as a counterstain for methyl green. It is also employed in conjunction with basic blue dyes, and when used in combination with hematoxylin it is referred to as "H & E" stains.

Eosin Y is an acid stain which interacts with cellular proteins rich in basic amino acids. A dye protein complex is formed that is characterized by vivid pink cytoplasmic staining.

Reagents

Eosin Y Solution, Alcoholic (Cat. No. HT1101: HT110116-500ML; HT110132-1L; HT110180-2.5L; HT1101128-4L)

Certified eosin Y, 0.5% (w/v), C.I. 45380, in acidified ethanol, 90% (v/v).

Eosin Y Solution, Aqueous (Cat. No. HT1102: HT110216-500ML; HT110232-1L; HT110280-2.5L; HT1102128-4L)

Certified eosin Y, 0.5% (w/v), C.I. 45380

Eosin Y Solution, Alcoholic, with Phloxine (Cat. No. HT1103: HT110316-500ML; HT110332-1L; HT110380-2.5L; HT1103128-4L)

Certified eosin Y, 0.1% (w/v), C.I. 45380, and certified phloxine B, 0.1% (w/v) in acidified ethanol.

Hematoxylin Solution, Harris Modified (Cat. No. HHS: HHS16-500ML; HHS32-1L; HHS80-2.5L; HHS128-4L)

Certified hematoxylin, 7.0 g/L, C.I. 75290, sodium iodate, aluminum ammonium sulfate • 12 H₂O, and stabilizers.

Hematoxylin Solution, Gill No. 1 (Cat. No. GHS1: GHS116-500ML; GHS132-1L; GHS1128-4L)

Certified hematoxylin, 2 g/L, C.I. 75290, sodium iodate, 0.2 g/L, aluminum ammonium sulfate • 12 H₂O, 17.6 g/L, and stabilizers.

Hematoxylin Solution, Gill No. 2 (Cat. No. GHS2: GHS216-500ML; GHS232-1L; GHS280-2.5L; GHS2128-4L)

Certified hematoxylin, 4 g/L, C.I. 75290, sodium iodate, 0.4 g/L, aluminum ammonium sulfate • 12 H₂O, 35.2 g/L, and stabilizers.

Hematoxylin Solution, Gill No. 3 (Cat. No. GHS3: GHS316-500ML; GHS332-1L; GHS380-2.5L; GHS3128-4L)

Certified hematoxylin, 6 g/L, C.I. 75290, sodium iodate, 0.6 g/L, aluminum ammonium sulfate • 12 H₂O, 52.8 g/L, and stabilizers.

Scott's Tap Water Substitute Concentrate (Cat. No. S5134-6x100ML)

Magnesium sulfate, 200 g/L, sodium bicarbonate, 20 g/L, and preservative.

Special Materials Required but Not Provided

- Ethanol, 95% and 100% or Reagent Alcohol
- Xylene or xylene substitute
- Synthetic Mounting Media
- Differentiation Solution (Cat. No. A3179-1L or A3429-4L)

Storage and Stability

All reagents are stored at room temperature (18–26°C). Reagents are stable until expiration date shown on the label or until pink tint is no longer observed.

Preparation

Eosin Y Alcoholic Solutions; Eosin Y Solution, Alcoholic, with Phloxine; and Harris Modified Hematoxylin solutions should be filtered before use. Hematoxylin solutions, Gill Nos. 1, 2 and 3 are provided ready to use.

Acidified Eosin Y Solution, Aqueous is prepared by slowly adding up to 0.5 mL of glacial acetic acid per 100 mL of stain.

Scott's Tap Water Substitute is prepared by diluting 1 part of Scott's Tap Water Substitute Concentrate with 9 parts of deionized water.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Standard histology texts provide the necessary details.^{2,3}

Notes

- The times given in the insert are approximate. Personal preferences will vary, and the times can be adjusted to suit personal preferences. Stain solutions which are heavily used will lose their staining powers and the staining times should be lengthened or new solutions should be used.⁴
- Positive control slides should be included in each run.
- Some tap water supplies are acidic and unsuitable for use in the "blueing" portion of this procedure. If tap water is acidic, use a dilute alkaline solution such as Scott's Tap Water Substitute.
- Purple or red-brown nuclei are indicative of inadequate "blueing".
- If eosin staining is excessive, nuclear staining may be masked. Proper eosin staining will demonstrate a 3-tone effect. To increase differentiation of eosin, extend time in alcohols or use a first alcohol with a higher water content. The times in the alcohols may be adjusted to obtain the proper degree of eosin staining.
- Addition of new stock to depleted working solutions of hematoxylin or eosin is not recommended.
- Avoid excessive water carry-over into alcoholic eosin solutions.

Procedures

Eosin Y Solutions, Alcoholic

- Deparaffinize to water or fix and hydrate frozen sections.
- Stain in hematoxylin.
- Rinse slides in running tap water.
- Differentiate if regressive hematoxylin is used. Rinse in running tap water.
- Blue in Scott's Tap Water Substitute.
- Rinse in running tap water.
- Rinse slides in 95% ethanol or reagent alcohol for 30 seconds.
- Counterstain in Alcoholic Eosin Y Solution (30 seconds – 3 minutes).
- Dehydrate, clear and mount.

Eosin Y Solution, Aqueous

- Deparaffinize to water or fix and hydrate frozen sections.
- Stain in hematoxylin.
- Rinse slides in running tap water.
- Differentiate if regressive hematoxylin is used. Rinse in running tap water.
- Blue in Scott's Tap Water Substitute.
- Rinse in running tap water.
- Aqueous Eosin Y Solution may be acidified by adding up to 0.5 mL of glacial acetic acid per 100 mL of stain.
- Counterstain in Acidified Aqueous Eosin Y Solution (30 seconds – 3 minutes).
- Dehydrate, clear and mount.

Performance Characteristics

Expected Results

Cytoplasm should be pink to red. Nuclei should be blue to blue-black, depending on the hematoxylin used.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT1101	EOSIN Y SOLUTION, ALCOHOLIC	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
HT1102	EOSIN Y SOLUTION, AQUEOUS	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
HT1103	EOSIN Y SOLUTION, ALCOHOLIC, WITH PHLOXINE B	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT110116, HT110132, HT110180, HT1101128:

H225: Highly flammable liquid and vapor.

H317: May cause an allergic skin reaction.

H319: Causes serious eye irritation.

H371: May cause damage to organs (Eyes).

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P242: Use non-sparking tools.

P243: Take action to prevent static discharges.

P260: Do not breathe mist or vapors.

P264: Wash skin thoroughly after handling.

P270: Do not eat, drink or smoke when using this product.

P272: Contaminated work clothing must not be allowed out of the workplace.

P280: Wear protective gloves/ protective clothing/ eye protection/face protection.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P311: IF exposed or concerned: Call a POISON CENTER/ doctor.

P333 + P313: If skin irritation or rash occurs: Get medical advice/ attention.

P337 + P313: If eye irritation persists: Get medical advice/ attention.

P370 + P378: In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

P403 + P235: Store in a well-ventilated place. Keep cool.

P405: Store locked up.

P501: Dispose of contents/ container to an approved waste disposal plant.

HT110216, HT110232, HT110280, HT1102128:

Not a hazardous substance or mixture.

HT110316, HT110332, HT110380, HT1103128:

H225: Highly flammable liquid and vapor.

H319: Causes serious eye irritation.

H371: May cause damage to organs (Eyes).

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P311: IF exposed or concerned: Call a POISON CENTER/ doctor.

HHS16, HHS32, HHS80, HHS128:

Not a hazardous substance or mixture.

GHS116, GHS132, GHS1128:

H302: Harmful if swallowed.

H373: May cause damage to organs (Kidney) through prolonged or repeated exposure if swallowed.

P260: Do not breathe dust/ fume/ gas/ mist/ vapors/ spray.

P264: Wash skin thoroughly after handling.

P270: Do not eat, drink or smoke when using this product.

P301 + P312: IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell.

P314: Get medical advice/ attention if you feel unwell.

P501: Dispose of contents/ container to an approved waste disposal plant.

GHS216, GHS232, GHS280, GHS2128:

H318: Causes serious eye damage.

P280: Wear eye protection/ face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

GHS316, GHS332, GHS3128:

H290: May be corrosive to metals.

H318: Causes serious eye damage.

H373: May cause damage to organs (Kidney) through prolonged or repeated exposure if swallowed.

P234: Keep only in original packaging.

P260: Do not breathe dust/ fume/ gas/ mist/ vapors/ spray.

P280: Wear eye protection/ face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P314: Get medical advice/ attention if you feel unwell.

P390: Absorb spillage to prevent material damage.

GHS380:

H226: Flammable liquid and vapor.

H290: May be corrosive to metals.

H302: Harmful if swallowed.

H318: Causes serious eye damage.

H373: May cause damage to organs (Kidney) through prolonged or repeated exposure if swallowed.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280: Wear eye protection/ face protection.

P301 + P312 + P330: IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Rinse mouth.

P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.

P314: Get medical advice/ attention if you feel unwell.

S5134:

H315: Causes skin irritation.

H317: May cause an allergic skin reaction.

H318: Causes serious eye damage.

P261: Avoid breathing mist or vapors.

P264: Wash skin thoroughly after handling.

P272: Contaminated work clothing should not be allowed out of the workplace.

P280: Wear protective gloves/ eye protection/ face protection.

P302 + P352: IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 232
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan and BB Hrapchak, Editors, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, NY, 2002, p 129

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 7.0	2022
Rev. 8.0	2022
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.
Rev. 9.0	2022
	Updated Warnings and Hazards section.
Rev. 10.0	2025
	Updated Warnings and Hazards section. Added CH-REP and importer information

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, S. 232
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan and BB Hrapchak, Editors, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, herausgegeben von Bancroft JD und Gamble, M, Churchill Livingstone, New York, 2002, S. 129

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Références

- Conn's Biological Stains, 10th ed., R.W. Horobin and J.A. Kiernan, Editors, Taylor & Francis, NY, 2002, p 232
- Theory and Practice of Histotechnology, 2nd ed., D.C. Sheehan and B.B. Hrapchak, Editors, C.V. Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., L.G. Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft J.D. and Gamble, M., Churchill Livingstone, New York, NY, 2002, p 129

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Referanslar

- Conn's Biological Stains, 10. baskı, RW Horobin and JA Kiernan, Editörler, Taylor & Francis, NY, 2002, s. 232
- Theory and Practice of Histotechnology, 2. baskı, DC Sheehan and BB Hrapchak, Editörler, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3. baskı, LG Luna, Editör, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Editörler: Bancroft JD ve Gamble, M, Churchill Livingstone, New York, NY, 2002, s. 129

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Referințe

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 232
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan and BB Hrapchak, Editors, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, NY, 2002, p129

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Reference

- Conn's Biological Stains, 10. izd., RW Horobin i JA Kiernan, urednici, Taylor & Francis, NY, 2002., str. 232
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan and BB Hrapchak, Editors, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, NY, 2002, p 129

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Atsauces

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, 232. lpp.
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan and BB Hrapchak, Editors, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, NY, 2002, 129. lpp.

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Odkazy

- Conn's Histological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 232
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan and BB Hrapchak, Editors, CV Mosby, St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, NY, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, NY, 2002, p 129

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Ak chcete urobiť objednávku, navštívte našu webovú lokalitu na [SigmaAldrich.com](https://www.sigmaaldrich.com). Informácie o technickom servise nájdete na stránke technického servisu na našej webovej lokalite na [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

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MDSS GmbH
MDSS CH GmbH
Chemie GmbH