



**FlowCollect™ Bcl-2 Activation Dual Detection Kit**  
**25 Tests**

**Cat. No. FCCS025108**

**FOR RESEARCH USE ONLY**  
**Not for use in diagnostic procedures.**

This page left blank intentionally.

---

## Application

The signaling pathways that regulate apoptosis are important in understanding tissue homeostasis, the development of the immune system, and therapies for cancer. Bcl-2 is a member of the Bcl-2 family that regulates apoptosis by regulating the proteins that control mitochondrial membrane permeability (1,2).

The Bcl-2 family consists of two protein sub-groups, the anti-apoptotic proteins and the pro-apoptotic proteins (3). Bcl-2 is an anti-apoptotic family member that binds to pro-apoptotic family members, such as Bax to sequester them and inhibit insertion into the mitochondrial membrane (4,5). When cells receive an apoptotic signal, Bcl-2 releases Bax, allowing Bax to form a complex on the mitochondrial membrane that releases cytochrome C into the cytoplasm leading to caspase activation and cell death (4).

Bcl-2 has been shown to be regulated at the level of phosphorylation by JNK at a multi-site phosphorylation loop which includes serine 70 (6). Phosphorylation of Bcl-2 at Serine 70 has been shown to be important in both apoptosis and in autophagy (6-9).

---

## Test Principle

Millipore's FlowCelect™ Bcl-2 Activation Dual Detection kit includes two directly conjugated antibodies: Anti-Bcl-2-Alexa Fluor® 488 and Anti-pBcl-2 (Serine 70)-PE, along with optimized fixation, permeabilization, wash, and assay buffers necessary for cell preparation and analysis. Millipore's FlowCelect™ Activation Dual Detection Kit contains sufficient reagents for 25 2-color samples. Detailed assay instructions are included to assist in your analysis and to ensure that the correct cell concentration is obtained during acquisition of sample data.

Millipore's FlowCelect™ Bcl-2 Activation Dual Detection Kit is designed to enable a researcher a quick and easy way to detect total Bcl-2 expression as well as the phosphorylation of Bcl-2 at serine 70 simultaneously. Millipore's FlowCelect™ Activation Dual Detection Kit was developed and tested using Jurkat cells treated with staurosporine, but the kit can be used with other human cell lines to determine the effect of other chemical reagents on Bcl-2 expression and phosphorylation.

---

## Kit Components

- 20X Anti-Bcl-2 Alexa Fluor<sup>®</sup> 488 Antibody (Part No. CS205126) One vial containing 150  $\mu$ L antibody.
- 20X Anti-pBcl-2(Ser70) PE Antibody (Part No. CS205125) One vial containing 150  $\mu$ L antibody.
- Fixation Buffer: (Part No. CS202122) One bottle containing 13 mL buffer.
- 10X Wash Buffer: (Part No. CS202123) One bottle containing 13 mL buffer.
- 5X Assay Buffer: (Part No. CS202124) One bottle containing 55 mL buffer.
- 1X Permeabilization Buffer: (Part No. CS203284) Two bottles containing 14 mL buffer.

---

## Materials Not Supplied

1. Flow Cytometer
2. Cells of interest and/or Jurkat cells as positive control
3. Media for cell line of interest
4. Pipettors with corresponding tips capable of accurately measuring 1 – 1000  $\mu$ L
5. Tabletop centrifuge capable of 2500 rpm.
6. Sample tubes capable of holding 1 mL
7. Staurosporine or other activating reagent

---

## Precautions

- The instructions provided have been designed to optimize the kit's performance. Deviation from the kit's instructions may result in suboptimal performance and may produce inaccurate data.
- Some assay components included in the kit may be harmful. The kit includes a fixation solution containing formaldehyde. Please refer to the MSDS sheet which can be found at [www.millipore.com](http://www.millipore.com) for specific information on hazardous materials.
- All fluorochrome conjugated antibodies are light sensitive and must be stored in the dark at 2-8°C.
- During storage and shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For maximum recovery of product, centrifuge vial briefly prior to removing cap.
- Do not use reagents beyond 4 months from date of receipt.

---

## Storage

Upon receipt, all antibodies and buffers should be stored at 2-8°C. **Caution:** *Fluorochrome conjugated antibodies should always be stored at 2-8°C. Do not freeze fluorescent antibodies. Any deviation in temperature for long periods of time may compromise the performance of the antibodies.*

---

## Preparation of Buffers

- Make 4 mL of 1X Fixation Buffer (Enough for 8 samples): Mix 1 mL of Fixation Buffer with 0.4 mL of 10X Wash Buffer and 8.6 mL of deionized water.
- Make 100 mL of 1X Assay Buffer: Mix 20 mL of 5X Assay Buffer with 80 mL of deionized water.
- Make 20 mL of 1X Wash Buffer: Mix 2 mL of 10X Wash Buffer with 18 mL of deionized water.
- Place Permeabilization Buffer on ice. *Permeabilization buffer must be ice cold for optimal results.*

**Note:** Prepared 1X Fixation Buffer is stable up to one month if stored at 2-8°C. Prepared 1X Wash and Assay Buffers are stable up to four months if stored at 2-8°C, and can be prepared when you receive the kit or when you start the first assay.

---

## Example Cell Treatment Protocol

**Note:** *This assay protocol has been optimized for Jurkat cells treated with staurosporine. However, this kit is suitable for flow cytometric analysis of a variety of human cell types and treatments but treatments may need to be optimized.*

### Cell Preparation for Staurosporine Treatment

1. Seed cells into 2 T-75 flasks ( $10 \times 10^6$  cells in 20 mL of media per flask).
2. Label one flask **Untreated** and label the other flask **Treated**.
3. Label the other flask **Treated**.
4. To the **Treated** flask and add 20  $\mu$ L of 1 mM staurosporine (in DMSO) for a final concentration of 1  $\mu$ M. To the **Untreated** flask and add 20  $\mu$ L of DMSO.
5. Incubate the flasks in a 37°C, 5% CO<sub>2</sub> incubator for 6-16 hours.
6. Count cells.

### Example Cell Staining Protocol:

- **Note:** *You will need  $0.5 \times 10^6$  cells for one test on Guava cytometer and  $1 \times 10^6$  cells for one test on other flow cytometers. Permeabilization Buffer must be ice cold prior to assay.*
1. Aliquot  $0.5 \times 10^6$  cells in media into a microfuge tube for each test to be performed. If using the assay for the first time, you will need to set up additional tubes to use for adjusting instrument settings as follows: **Untreated** cells will be labeled tube 3. **Treated** cells will be labeled tube 1, 2, and 4 (see table below).
  2. Spin down cells at 2500 rpm (600Xg) for 3 minutes and discard media.
  3. Add 0.5 mL of 1X Wash Buffer by gently pipetting up and down.

4. Spin down cells at 2500 rpm (600Xg) for 3 minutes and discard buffer.
5. Add 0.5 mL of 1X Fixation Buffer by gently pipetting up and down.
6. Incubate at room temperature for 20 minutes.
7. Spin down cells at 2500 rpm (600Xg) for 3 minutes and discard buffer.
8. Resuspend cells in 0.5 mL 1X Assay Buffer by gently pipetting up and down.
9. Spin down cells at 2500 rpm (600Xg) for 3 minutes and discard buffer.
10. Resuspend cells in 0.5 mL of **ice cold** 1X Permeabilization Buffer by gently pipetting up and down.
11. Incubate samples on ice for 10 minutes. (If necessary, cells can be stored in 1X Permeabilization Buffer overnight at 4 degrees and the assay can be continued the next day.)
12. Spin down cells at 2500 rpm (600Xg) for 3 minutes and discard buffer.
13. Resuspend cells in 0.5 mL 1X Assay Buffer by gently pipetting up and down.
14. Spin down cells at 2500 rpm (600Xg) for 3 minutes and discard buffer.
15. Resuspend cells in tubes 1 and 2 in 95  $\mu$ L 1X Assay Buffer by pipetting up and down. Resuspend cells in tubes 3 and 4 in 90  $\mu$ L 1X Assay Buffer by pipetting up and down.
16. Add 5  $\mu$ L of 20X antibody to tubes using the following chart.

<b>Tube #</b>	<b>Description</b>
1	Treated cells stained with Bcl-2 AlexaFluor <sup>®</sup> 488
2	Treated cells stained with pBcl-2 (Ser70)-PE
3	Untreated cells stained with both antibodies
4	Treated cells stained with both antibodies

17. Incubate tubes in the dark on ice for 60 minutes.
18. Add 900  $\mu$ L of 1X Assay Buffer to each tube then pellet the cells at 2500 rpm for 3 minutes and discard buffer.
19. If using a Guava cytometer, resuspend cells in 500  $\mu$ L of 1X Assay Buffer. If using another cytometer, resuspend cells in 1 mL of 1X Assay Buffer.

20. Transfer 200  $\mu$ L of cells to a 96-well plate or transfer 500  $\mu$ L- 1 mL to a sample tube and analyze on a Flow Cytometer.

**Note:** Samples 1 and 2 are used for adjusting compensation. Once the initial compensation for a given cell type/treatment has been performed the settings file can be saved and used for subsequent data collection without the need for samples 1 and 2.

## Sample Data




Figure 1: A) Single color histogram analysis of Bcl-2 expression. Jurkat cells are stained with anti-Bcl-2-AlexaFlour<sup>®</sup>488 (green shaded histogram) or with isotype control (grey shaded histogram). B) Single color histogram analysis of Bcl-2 phosphorylation at serine 70. Jurkat cells treated with staurosporine for 6 hours (yellow shaded histogram) or untreated Jurkat cells (grey shaded histogram) and stained with anti-phospho Bcl-2(S70)-PE.




Figure 2: Multi-color dot blot analysis of Bcl-2 expression and Bcl-2 phosphorylation at serine 70. Bcl-2 expression remains constant but pBcl-2(S70) phosphorylation decreases upon staurosporine treatment. Untreated Jurkat cells (A) and Jurkat cells treated with staurosporine for 6 hours (B) were stained with anti-phospho Bcl-2(S70)-PE and anti-Bcl-2-AlexaFlour<sup>®</sup>488.

---

## Technical Hints

- Permeabilization Buffer must be ice cold prior to assay
- For cellular staining and analysis to be most effective, make sure that test cells have good viability prior to use.
- Certain cell cultures cell pellets may become hazy or transparent following the fixation step, making it difficult to see. After the washing steps they will become easier to visualize.
- If sampling a small collection of cells for flow analysis, it is recommended that all steps be performed in a smaller collection tube (e.g. centrifuge tube)
- Do not mix or interchange reagents from various kit lots.
- If minor precipitate is detected in the 10X Wash Buffer place the bottle in a warm water bath for 30 minutes, followed by mixing the contents on a mechanical vortex.

---

## Troubleshooting

Potential Problem	Experimental Suggestions
Acquisition rate decreases dramatically Instrument clogging Too many cells	<ul style="list-style-type: none"><li>• Cell concentration too high - Decrease the number of cells per microliter by diluting sample to 300 – 500 cells per microliter. Guava cytometers give the most accurate data when the flow rate is less 500 cells per microliter.</li><li>• Run three Quick Cleans (for Guava cytometers) to rinse out the flow cell. This procedure can be performed during or after an assay. This will wash away any material forming within the glass capillary walls.</li></ul>
Too few cells	<ul style="list-style-type: none"><li>• Spin down cells and resuspend in a smaller volume. The assay instructions are optimized to give you a range of cells between 100-500 cells/<math>\mu</math>L in the final sample volume. However, cell loss is common during washing steps in the assay procedure. A substantial decrease in cell numbers can lead to difficulty in adjusting settings. Make sure to leave the cell pellet intact when discarding buffer. If the cells are not generating a compact pellet after centrifugation, increase the time to 5 minutes and/or increase the speed by 500 rpm until a compact and visible cell pellet forms.</li></ul>
Background staining and/or non-specific staining of cells	<ul style="list-style-type: none"><li>• This assay was optimized using Jurkat cells. Therefore, further antibody titrations may be necessary for other cell types and conditions to capture the ideal staining concentration. Non-specific staining and background may indicate that less antibody will need to be used during the staining procedure.</li></ul>



Variability in day to day experiments	<ul style="list-style-type: none"> <li>• Monitor experimental cell cultures to ensure that cell viability and cell numbers being analyzed are consistent. Any changes in culture conditions or viability can influence experimental results.</li> <li>• Make sure that a quality check on the flow cytometer to be used (e.g. calibration) is performed on a daily basis prior to use.</li> </ul>
Staining is weak	<ul style="list-style-type: none"> <li>• Some cell lines may require higher concentrations of fixation buffer to stain correctly. Use twice as much Fixation Buffer.</li> </ul>
Staining is weak	<ul style="list-style-type: none"> <li>• Some cell lines may require higher concentrations of conjugated antibodies. Try using twice as much antibody for your staining.</li> </ul>
Staining is too bright	<ul style="list-style-type: none"> <li>• Some cell lines may require lower concentrations of conjugated antibodies. Try using half as much antibody for your staining.</li> </ul>

*\*For further support, please contact Millipore's Technical services at +1(800) 437-7500*

---

## References

1. Cory S, Adams JM. The Bcl-2 family: regulators of the cellular life-or-death switch. *Nat Rev Can.* 2002;2:647–656.
2. Adams JM, Cory S. The Bcl-2 protein family: arbiters of cell survival. *Science (New York, NY.* 28. 1998;281(5381):1322–6.
3. Huang DC, Strasser A. BH3-Only proteins-essential initiators of apoptotic cell death. *Cell.* 8. 2000;103(6):839–42.
4. Gross A, Jockel J, Wei MC, Korsmeyer SJ. Enforced dimerization of BAX results in its translocation, mitochondrial dysfunction and apoptosis. *The EMBO journal.* 15. 1998;17(14):3878–85.
5. Yang E, Zha J, Jockel J, Boise LH, Thompson CB, Korsmeyer SJ. Bad, a heterodimeric partner for Bcl-XL and Bcl-2, displaces Bax and promotes cell death. *Cell.* 27. 1995;80(2):285–91.
6. Wei Y, Sinha S, Levine B. Dual role of JNK-1 mediated phosphorylation of Bcl-2 in autophagy and apoptosis regulation. *Autophagy*, 2008;4(7):949-51
7. Ito, T., Deng, X., Carr, B., May, W. S. Bcl-2 phosphorylation required for anti-apoptosis function. *J. Biol. Chem.* 2008;272,11671-11673
8. Haldar, S., Jena, N., Croce, C. M. (1995) Inactivation of Bcl-2 by phosphorylation. *Proc. Natl. Acad. Sci. USA* 92,4507-4511
9. Haldar, S., Chintapalli, J., Croce, C. M. (1996) Taxol induces bcl-2 phosphorylation and death of prostate cancer cells. *Cancer Res* 56,1253-1255.

---

## Warranty

**Millipore Corporation** ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the Company all monies paid for such applicable Product. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY COMPANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

(c) 2010: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing

**Cat No. FCCS025108**

05/2010  
Revision A: FCCS025108MAN