

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# **ProductInformation**

## Chlorpromazine hydrochloride

Product Number C8138 Store at Room Temperature

### **Product Description**

Molecular Formula: C<sub>17</sub>H<sub>19</sub>ClN<sub>2</sub>S · HCl

Molecular Weight: 355.3 CAS Number: 69-09-0

Melting Point: 179-180 °C (with decomposition, capillary method); 194-196 °C (microblock method)<sup>1</sup>

 $\lambda_{max}$ : 253 nm (0.12 M HCl) Extinction Coefficient:  $E^{mM} = 31.5$  (0.12 M HCl)

Synonyms: 2-chloro-N,N-dimethyl-10H-phenothiazine-10-propanamine:

2-chloro-10-(3-dimethylaminopropyl)phenothiazine; N-(3-dimethylaminopropyl)-3-chlorophenothiazine<sup>1</sup>

Chlorpromazine is a compound that is used in cell signaling and neuroscience research. It is a dopamine inhibitor, notably at the D2-dopaminergic receptor. It also has α-adrenergic blocking activity and stimulates prolactin release. 1,2 The antifungal properties of chlorpromazine and other compounds has been reviewed.<sup>3</sup> The role of chlorpromazine and other tricyclic compounds as inhibitors of neuronal voltagegated K channels has been discussed.4

The use of chlorpromazine (0 - 0.2 mM) to induce stress in *E. coli* to study regulation of the RcsCB His-Asp phosphorelay system has been described.5 The effect of chlorpromazine (1 µM) and other dopamine-active drugs on the production of oxytocin and vasopressin in rat neurohypophyseal tissue cultures has been studied. 6 Chlorpromazine has been used to probe long-term potentiation of corticostriatal synaptic transmission in rat. Chlorpromazine (12 or 25 µM) has been shown to limit lipid accumulation in the human breast cell line HBL-100, as initially induced by tetraphenylphosphonium chloride.

An analytical method for chlorpromazine that couples high performance liquid chromatography to inductively coupled plasma mass spectroscopy has been described.

#### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

#### **Preparation Instructions**

This product is soluble in water (50 mg/ml), yielding a clear, faint yellow to very faint gray solution.

#### References

- 1. The Merck Index, 12th ed., Entry# 2238.
- Martindale The Extra Pharmacopoeia, 31st ed., Reynolds, J. E. F., ed., Royal Pharmaceutical Society (London, UK: 1996), pp. 687-695.
- Afeltra, J., and Verweij, P. E., Antifungal activity of nonantifungal drugs. Eur. J. Clin. Microbiol. Infect. Dis., **22(7)**, 397-407 (2003).
- Mathie, A., et al., Voltage-activated potassium channels in mammalian neurons and their block by novel pharmacological agents. Gen. Pharmacol., **30(1)**, 13-24 (1998).
- Conter, A., et al., The RcsCB His-Asp phosphorelay system is essential to overcome chlorpromazine-induced stress in Escherichia coli. J. Bacteriol., 184(10), 2850-2853 (2002).
- Galfi, M., et al., Effects of dopamine and dopamine-active compounds on oxytocin and vasopressin production in rat neurohypophyseal tissue cultures. Regul. Pept., 98(1-2), 49-54 (2001).
- 7. Zhang, C., et al., Regeneration of dopaminergic function in 6-hydroxydopamine-lesioned rats by neuroimmunophilin ligand treatment. J. Neurosci., 21(15), RC156, 1-6 (2001).
- Sathasivam, N., et al., Inhibition of tetraphenylphosphonium-induced NMR-visible lipid accumulation in human breast cells by chlorpromazine. Biochim. Biophys. Acta, 1633(3), 149-160 (2003).
- 9. Duckett, C. J., et al., Quantitation in gradient high performance liquid chromatography/inductively coupled mass spectrometry investigated using diclofenac and chlorpromazine. Rapid Commun. Mass Spectrom., 16(4), 245-247 (2002).

GCY/JRC 5/06