

User Guide

Millicell® Standing Cell Culture Inserts

MILlicell® (MCE)

PIHA01250, PIHA03050

MILlicell® (PTFE)

PICM01250, PICM03050

MILlicell® (PC)

PIHP01250, PIHP03050, PITP01250, PI8P01250, PIXP01250

Introduction


Millicell® cell culture inserts are sterile, general purpose devices for the growth and differentiation of various cell types.

Millicell® standing cell culture inserts are single-use products consisting of a filter sealed to a cylindrical polystyrene holder. The 12 mm inserts are used in 12 or 24 well culture plates. The 30 mm inserts are used in 6-well culture plates.

NOTE: When working with Millicell® inserts, use appropriate tissue culture media, reagents, and growth environments as well as aseptic technique.

For information on ECM coating (required for attachment dependent cells on the Millicell® insert with Biopore™ membrane), staining, or other Millicell® related procedures, contact Technical Service at SigmaAldrich.com/techservice.

Diagram of Millicell® Standing Cell Culture Insert



Directions

Perform the following steps in a bacteriological LAF (laminar flow) hood or equivalent controlled environment.

CAUTION: Do not use at temperatures above 50 °C to ensure accurate results. Do not reuse inserts.

1. Prepare culture plate wells before opening the Millicell® inserts. Add tissue culture medium to each well as follows:

- 600 µL for 12 mm inserts, or
- 1.5 to 2.0 mL for 30 mm inserts.

Prepare enough wells for experiments and controls.

CAUTION: Do not use with strong acids or bases, or organic solvents incompatible with the membrane, and/or polystyrene.

2. Peel off the cover sheet from the blister package containing the Millicell® insert.

CAUTION: Do not use inserts that have damaged membrane, such as cracks or holes.

3. Use sterilized forceps to remove the insert from the package. Do not touch the membrane.

4. Repeat steps 2 and 3 until a Millicell® insert is placed into each prepared well. Make sure the inserts are positioned correctly, with the membrane section closest to the well bottom.

5. Allow several minutes for the membrane in each insert to become moistened with the tissue culture medium that was added in step 1. If using:

- **12 mm inserts**, the recommended inside volume is 400 µL and the outside volume is 600 µL.
- **30 mm inserts**, the recommended inside and outside volumes are 1.5 to 2.0 mL each.

- Inoculate the cells onto the inside of the insert above the membrane.
- Follow standard tissue culture incubation and feeding procedures for cell growth and monolayer formation.

NOTE: Be sure not to puncture the membrane or disturb cultured cells during medium addition or removal.

Specifications

Dimensions	12 mm	30 mm
Height (including feet):	1.0 cm	1.3 cm
Height of Feet:	1–2 mm	1–2 mm
Outer Diameter:	1.2 cm	3.0 cm
Inner Diameter:	1.0 cm	2.7 cm
Membrane Area (effective):	0.6 cm ²	4.2 cm ²
Temperature Limit:	50 °C	50 °C
Membrane Pore Sizes:	0.4, 0.45, 3.0, 8.0, and 12.0 µm	
Solvent Compatibility:	Incompatible with strong acids, strong bases, and certain organic solvents.	
Sterility:	Inserts are shipped sterile and ready-for-use. If the package is open or damaged before use, sterility cannot be guaranteed.	

Materials of Construction

Available Membrane Types	
MF-Millipore™:	Mixed cellulose esters (MCE)
Biopore™:	Hydrophilic Polytetrafluoroethylene (PTFE)
Isopore™:	Tissue culture-treated, track-etched PVP-free polycarbonate (PC)
Plastic Holder	Polystyrene

Product Ordering

Purchase products online at SigmaAldrich.com/products.

Millicell® Membrane	Pore Size	Insert Size	Cat. No. (50/box)
MF-Millipore™ (MCE)	0.45 µm	12 mm 30 mm	PIHA01250 PIHA03050
Biopore™ (PTFE)	0.4 µm	12 mm 30 mm	PICM01250 PICM03050
Isopore™ (PC)	0.4 µm	12 mm 30 mm	PIHP01250 PIHP03050
	3.0 µm	12 mm	PITP01250
	8.0 µm	12 mm	PI8P01250
	12.0 µm	12 mm	PIXP01250

Organotype (Low height: 5 mm)

Biopore™ (PTFE)	0.4 µm	30 mm	PICM0RG50
-----------------	--------	-------	-----------

Related Products

Millicell®-ERS Electrical Resistance System (measures membrane potential and resistance of epithelial cells in culture)	MERS00002
---	-----------

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Merck, Millipore, Biopore, Isopore, Millicell, MF-Millipore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.