

For life science research only.
Not for use in diagnostic procedures.



TriPure Isolation Reagent

 **Version: 09**

Content Version: April 2021

Monophasic solution of phenol and guanidine thiocyanate for RNA, DNA, and protein isolation.

Cat. No. 11 667 157 001 50 ml

Cat. No. 11 667 165 001 200 ml

Store the product at +2 to +25°C.

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	3
	Storage Conditions (Product)	3
1.3.	Additional Equipment and Reagent required	3
1.4.	Application	4
2.	How to Use this Product	5
2.1.	Before you Begin	5
	General Considerations	5
	Prevention of RNase contamination	5
	Safety Information	5
	Precautions	5
2.2.	Protocols	5
	Protocol 1 Extraction of RNA, DNA, and protein from tissues or cells	5
	Protocol 2 Phase separation	7
	Protocol 3 Isolation of RNA	7
	Protocol 4 Isolation of DNA	8
	Protocol 5 Isolation of protein	10
3.	Results	11
4.	Troubleshooting	12
5.	Additional Information on this Product	13
5.1.	Test Principle	13
	Isolation procedures	13
6.	Supplementary Information	14
6.1.	Conventions	14
6.2.	Changes to previous version	14
6.3.	Ordering Information	14
6.4.	Trademarks	15
6.5.	License Disclaimer	15
6.6.	Regulatory Disclaimer	15
6.7.	Safety Data Sheet	15
6.8.	Contact and Support	15

1. General Information

1.1. Contents

Vial / bottle	Label	Function / description	Catalog number	Content
1	TriPure Isolation Reagent	Clear, red ready-to-use solution.	11 667 157 001	1 bottle, 50 ml
			11 667 165 001	1 bottle, 200 ml

1.2. Storage and Stability

Storage Conditions (Product)

When stored at +2 to +25°C, the product is stable through the expiry date printed on the label.

Vial / bottle	Label	Storage
1	TriPure Isolation Reagent	Store at +2 to +25°C. ⚠ Keep in the original polypropylene bottle and do not expose to light for long periods of time (days).

1.3. Additional Equipment and Reagent required

For extraction and phase separation

- Sterile, disposable polypropylene tubes that can withstand 12,000 × *g* in the presence of TriPure Isolation Reagent and chloroform
- Homogenization apparatus (tissue and certain cells only)
- Reagents for density gradients (white blood cells only)
- Chloroform (free of all additives, such as isoamyl alcohol)
- Glycogen* (for tissue samples < 10 mg)

For RNA isolation

- Isopropanol
- 75% ethanol
- Diethylpyrocarbonate (DEPC)-treated RNase-free water or DEPC-treated 0.5% SDS*

For DNA isolation

- 100% ethanol
- 75% ethanol
- 8 mM NaOH
- 0.1 M sodium citrate in 10% ethanol

For protein isolation

- Isopropanol
- 1% SDS*
- 100% ethanol
- 0.3 M guanidine hydrochloride in 95% ethanol

1.4. Application

TriPure Isolation Reagent is used for the isolation of:

- Total RNA, DNA, and protein from the same sample in a single-step liquid-phase separation. This procedure is an improvement of the single-step RNA isolation method.
- This reagent performs well with both small and large quantities of tissue or cells from human, animal, plant, or bacterial origin.

 *RNA is undegraded and free of protein or DNA contamination.*

2. How to Use this Product

2.1. Before you Begin

General Considerations

Prevention of RNase contamination

- Wear clean disposable gloves at all times.
- Use sterile disposable plasticware and pipettes/tips reserved for RNA work only.

Safety Information

Precautions

- TriPure Isolation Reagent contains phenol, a poison and guanidine thiocyanate, an irritant. This reagent is toxic. This reagent can also cause burns if it touches the skin.
- Do not allow guanidine thiocyanate to contact sodium hypochlorite (bleach) solution or acids. These mixtures produce a highly toxic gas.
- When using TriPure Isolation Reagent, work under a fume hood.
- Wear gloves, lab coat, and appropriate eye protection, such as a shield or safety goggles.
- Avoid breathing vapor.
- If you get the reagent on your skin, wash skin immediately with soap or mild detergent, and flush with large amounts of water for 15 to 30 minutes, until no evidence of chemical remains.
- If you get the reagent in your eyes, flush eyes with large amounts of water for 15 to 30 minutes, until no evidence of chemical remains.

2.2. Protocols

Protocol 1 Extraction of RNA, DNA, and protein from tissues or cells

Starting with tissue

i Keep track of the amount of TriPure Isolation Reagent that you use for each sample in this procedure. The amount of TriPure used will determine the amount of reagents required in later procedures.

⚠ **TriPure Isolation Reagent is a hazardous reagent. Observe all safety precautions, see section, Safety Information, Precautions.**

- 1 For each 50 to 100 mg of tissue to be processed, add 1 ml TriPure Isolation Reagent to a polypropylene centrifuge tube at +15 to +25°C.
 - The capacity of the tube should be at least twice the volume of the tissue and TriPure Isolation Reagent combined.

i For small (<10 mg) tissue samples, use 0.8 ml TriPure Isolation Reagent plus 5 to 10 µg Glycogen*, used as a carrier to assist RNA precipitation.

⚠ **The volume of the sample should not be greater than 10% of the volume of TriPure Isolation Reagent.**

- 2 Add fresh or frozen tissue samples to the tube.

- 3 Homogenize the tissue with a glass-Teflon or power homogenizer (Polytron or equivalent).

i The molecular weight of DNA isolated in Protocol 4 is dependent on the shearing forces applied during homogenization. Use a loosely fitting homogenizer, and avoid use of a Polytron homogenizer if you want high molecular weight DNA. The use of a hand-held homogenizer may decrease overall yields but will minimize shearing of high molecular weight DNA.
-

2. How to Use this Product

- Optional step: only necessary for tissues that contain large amounts of proteins, fat, polysaccharides, or extracellular material, such as muscle tissue, fatty tissue, or tuberous parts of plants.
 - Clarify the homogenate by centrifuging it at $12,000 \times g$ for 10 minutes at +2 to +8°C.
 - Extracellular membranes, polysaccharides, and high molecular weight DNA pellet during the centrifugation, while excess fat collects as a layer on top of the supernatant.
 - Remove any fatty layer, then transfer the supernatant which contains the RNA, most DNA, and protein to a fresh polypropylene centrifuge tube.
-

Starting with adherent cells grown in a monolayer option A

- Decant culture medium from cells.
 - At +15 to +25°C, add directly to the culture dish or flask, 1 ml TriPure Isolation Reagent for each 10 cm² area covered by cells, regardless of cell number.
 - Adding too little TriPure Isolation Reagent may lead to contamination of isolated RNA with DNA.*
 - Pass the cell lysate through a pipette several times.
 - Transfer the cell lysate to a polypropylene centrifuge tube.
 - The capacity of the tube should be at least twice the volume of the cells and TriPure Isolation Reagent combined.
-

Starting with adherent cells grown in a monolayer option B

- Decant culture medium from cells.
 - At +15 to +25°C, add directly to the culture dish or flask, 3 ml TriPure Isolation Reagent for each 5 to 10×10^6 cells.
 - Scrape cells from the walls of the flask or dish into the reagent to lyse them.
 - Pass the cell lysate through a pipette several times.
 - Transfer the cell lysate to a polypropylene centrifuge tube.
 - The capacity of the tube should be at least twice the volume of the cells and TriPure Isolation Reagent combined.
-

Starting with cells grown in suspension

- Pellet cells by centrifugation in a polypropylene centrifuge tube.
 - The capacity of the tube should be at least twice the volume of the pelleted cells and TriPure Isolation Reagent combined.
 - Remove the supernatant.
 - ⚠ Do not wash cells before adding TriPure Isolation Reagent.**
 - Add TriPure Isolation Reagent directly to the tube containing the cell pellet at +15 to +25°C.
 - Use 1 ml reagent for each 5 to 10×10^6 animal, plant, or yeast cells; or 1 ml reagent for each 1×10^7 bacterial cells.
 - For small amounts (<10⁶) of cells, use 0.8 ml TriPure Isolation Reagent plus 5 to 10 µg Glycogen*, used as a carrier to assist RNA precipitation.*
 - Lyse cells by repetitive pipetting or homogenization.
 - The molecular weight of the isolated DNA is dependent on the shearing forces applied during homogenization. Use a loosely fitting homogenizer and avoid use of a Polytron homogenizer if you want high molecular weight DNA.*
-

- Transfer the cell lysate to a polypropylene centrifuge tube.
-

Starting with human blood

- Isolate buffy coat or lymphocytes from human whole blood using low-speed centrifugation or density gradients.
 - For each 5 to 10×10^6 white blood cells, add 1 ml TriPure Isolation Reagent directly to the centrifuge tube.
 - Lyse cells by repetitive pipetting or homogenization.
 - The molecular weight of the isolated DNA is dependent on the shearing forces applied during homogenization. Use a loosely fitting homogenizer and avoid use of a Polytron homogenizer if you want high molecular weight DNA.*
-

Protocol 2 Phase separation

- The sample homogenates obtained in Protocol 1 (Extraction) can be stored at -60°C or below for at least 1 month before they are used in phase separation.*
- Incubate each homogenized sample from Protocol 1 for 5 minutes at $+15$ to $+25^\circ\text{C}$ to ensure the full dissociation of nucleoprotein complexes.
 - Add chloroform to each sample.
 - Use 0.2 ml chloroform for each 1 ml TriPure Isolation Reagent required in the initial homogenization (Protocol 1).
 - Cap tube securely, and shake it vigorously for 15 seconds.
 - Incubate tube at $+15$ to $+25^\circ\text{C}$ for 2 to 15 minutes.
 - To separate the solution into three phases, centrifuge tube at $12,000 \times g$ for 15 minutes at $+2$ to $+8^\circ\text{C}$.
 - ⚠ Do not exceed $12,000 \times g$ during centrifugation.**
 - After centrifugation, use the three phases of each solution as follows:
 - From the colorless upper aqueous phase (which is approximately 60% of the volume of TriPure Isolation Reagent used in Protocol 1), isolate RNA according to Protocol 3.
 - From the interphase and lower red organic phase, isolate DNA according to Protocol 4, and protein according to Protocol 5.
-

Protocol 3 Isolation of RNA

- Transfer the colorless upper aqueous phase obtained in Protocol 2 to a new polypropylene centrifuge tube.
 - Save the interphase and lower red organic phase obtained in Protocol 2 from each sample, and use these phases to isolate DNA according to Protocol 4, and protein according to Protocol 5.*
 - Precipitate the RNA from the colorless aqueous phase of each sample by performing the following steps:
 - Add isopropanol to the aqueous phase; use 0.5 ml isopropanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization (Protocol 1).
 - Cap the tube, then invert it several times to mix thoroughly.
 - Incubate sample for 5 to 10 minutes at $+15$ to $+25^\circ\text{C}$ to allow the RNA precipitate to form.
 - Centrifuge the sample at $12,000 \times g$ for 10 minutes at $+2$ to $+8^\circ\text{C}$.
 - Discard the supernatant.
-

2. How to Use this Product

- 3 Add 75% ethanol to each centrifuge tube.
 - Use at least 1 ml of 75% ethanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization procedure.
 - i* The RNA precipitate can be stored in 75% ethanol for at least 1 week at +2 to +8°C, or at least 1 year at –15 to –25°C.

- 4 Wash the RNA pellet in the ethanol by vortexing.

- 5 Centrifuge the sample at 7,500 × *g* for 5 minutes at +2 to +8°C.
 - Discard the supernatant.

- 6 Remove the excess ethanol from the RNA pellet by air-drying or placing the sample under vacuum for 5 to 10 minutes.
 - ⚠ Do not dry the RNA pellet by centrifugation under vacuum. Do not let the RNA pellet dry completely as a dry pellet will be much less soluble.**

- 7 Resuspend the RNA pellet in diethylpyrocarbonate (DEPC)-treated RNase-free water or DEPC-treated 0.5% SDS.

- 8 Dissolve the RNA pellet by passing the solution through a pipette tip several times, then incubating the solution for 10 to 15 minutes at +55°C to +60°C.

Protocol 4 Isolation of DNA

- i* The interphase and red organic phase obtained in Protocol 2 can be stored at +2 to +8°C overnight before they are used in Protocol 4.
- 1 Carefully remove all the remaining colorless upper aqueous phase from each sample obtained in Protocol 2, and discard.
 - i* You must get rid of all the aqueous phase in order to obtain high quality DNA.

- 2 Precipitate the DNA from the interphase and red organic phase of each sample by performing the following steps:
 - Add 100% ethanol to the interphase and organic phase; use 0.3 ml 100% ethanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization (Protocol 1).
 - Cap the sample, then invert it several times to mix it thoroughly.
 - Incubate sample for 2 to 3 minutes at +15 to +25°C to allow the DNA precipitate to form.
 - Centrifuge the sample at 2,000 × *g* for 5 minutes at +2 to +8°C.

- 3 Remove the supernatant containing phenol, ethanol, and protein from each sample, and save at +2 to +8°C.
 - Use this supernatant to isolate protein according to Protocol 5.

- 4 To remove any phenol present in the DNA from each sample, perform the following steps:
 - Add 0.1 M sodium citrate in 10% ethanol to the pellet remaining in the centrifuge tube; use 1 ml sodium citrate/ethanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization procedure.
 - Incubate the sample with occasional mixing for 30 minutes at +15 to +25°C.
 - Centrifuge the sample at 2,000 × *g* for 5 minutes at +2 to +8°C.
 - Discard the supernatant.

- 5 Repeat Step 4 twice for a total of 3 sodium citrate/ethanol washes.

- 6 After the 3 sodium citrate/ethanol washes, perform the following steps for each DNA precipitate:
 - Wash each DNA pellet in 75% ethanol; use 1.5 to 2.0 ml ethanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization procedure.
 - i* DNA suspended in 75% ethanol can be stored at +2 to +8°C for several months.
 - Incubate sample with occasional mixing at +15 to +25°C for 10 to 20 minutes.
 - Centrifuge the sample at 2,000 × *g* for 5 minutes at +2 to +8°C.
 - Discard the supernatant.

- 7 Remove the excess ethanol from each DNA pellet by air-drying or placing the sample under vacuum for 5 to 10 minutes.
-
- 8 Dissolve the DNA pellet by performing the following steps:
- To the pellet, add enough 8 mM NaOH to approach a DNA concentration of 0.2 to 0.3 µg/µl. For example, use 0.3 to 0.6 ml 8 mM NaOH to resuspend DNA isolated from 50 to 70 mg tissue or 10⁷ cells.
 - Pass the sample through a pipette.
- i* DNA will fully dissolve in the mild alkaline solution. However, the dissolved DNA (especially from tissue samples) may still contain insoluble gel-like material (fragments of membranes, etc.).
- (optional) If the dissolved DNA contains insoluble material, pellet the insoluble material by centrifugation at 12,000 × *g* for 10 minutes at +2 to +8°C, then transfer the clear supernatant (DNA) to a new tube before storing the DNA.
-
- 9 For long-term storage, adjust the pH of the DNA solution to pH 7 to 8 with HEPES, and adjust the EDTA concentration to 1 mM.
- i* Samples solubilized in 8 mM NaOH can be stored overnight at +2 to +8°C without neutralization.
-

Guidelines for using DNA isolated with TriPure Isolation Reagent

Calculation of cell number

To calculate cell number from DNA content, assume the amount of DNA in 10⁶ diploid cells of human, rat, and mouse origin equals 7.1 µg, 6.5 µg, and 5.8 µg, respectively.

Digestion of DNA by restriction endonucleases

- 1 Adjust the pH of the DNA solution isolated in Protocol 4 to the required value with HEPES according to the table, see **Amplification of DNA by PCR, Step 1**.
-
- 2 Alternatively, dialyze samples against 1 mM EDTA, pH 7.0 to 8.0.
-
- 3 Digest the DNA for 3 to 24 hours, with 3 to 5 units of enzyme per µg DNA under conditions optimal for the specific restriction enzyme. In a typical restriction assay, 80 to 90% of the DNA preparation is digestible.
-

Amplification of DNA by PCR

- 1 Adjust the pH of the DNA isolated in Protocol 4 to 8.4 using 0.1 M HEPES as shown in the following table.

For a final pH of...	add this amount of 0.1 M HEPES [µl]	or this amount of 1.0 M HEPES [µl]
8.4	86	-
8.2	93	-
8.0	101	-
7.8	117	-
7.5	159	-
7.2	-	23
7.0	-	32

- 2 Add 0.1 to 1.0 µg of the DNA sample to a standard PCR.
-

2. How to Use this Product

Protocol 5 Isolation of protein

i Isolate proteins from the phenol-ethanol supernatant obtained in Protocol 4, Step 3.

Standard Protocol 5

- 1 Precipitate proteins from the phenol-ethanol supernatant of each sample by performing the following steps:
 - Add isopropanol to the phenol-ethanol supernatant; use 1.5 ml isopropanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization (Protocol 1).
 - Cap the tube, then invert it several times to mix it thoroughly.
 - Incubate sample for a minimum of 10 minutes at +15 to +25°C to allow the protein precipitate to form.
 - Centrifuge the sample at 12,000 × *g* for 5 minutes at +2 to +8°C.
 - Discard the supernatant.

- 2 To wash the precipitated protein in each sample, perform the following steps:
 - Resuspend each protein pellet with 0.3 M guanidine hydrochloride in 95% ethanol. Use 2 ml guanidine hydrochloride/ethanol for each 1 ml TriPure Isolation Reagent required in the initial homogenization procedure.

i The protein pellet suspended in 0.3 M guanidine hydrochloride/95% ethanol can be stored for at least 1 month at +2 to +8°C or for at least one year at –15 to –25°C.

- Mix the sample, and incubate at +15 to +25°C for 20 minutes.
- Centrifuge the sample at 7,500 × *g* for 5 minutes at +2 to +8°C.
- Discard the supernatant.

- 3 Repeat Step 2 twice for a total of 3 guanidine hydrochloride/ethanol washes.

- 4 After the 3 washes, perform the following steps on each protein pellet:
 - Add 2 ml 100% ethanol to each protein pellet.
 - Vortex to wash the pellet.
 - Incubate at +15 to +25°C for 20 minutes.
 - Centrifuge at 7,500 × *g* for 5 minutes at +2 to +8°C.
 - Discard supernatant.

- 5 Remove the excess ethanol from each protein pellet by air-drying or placing the sample under vacuum for 5 to 10 minutes.

- 6 Dissolve the protein pellet by adding 1% SDS to it and repeatedly passing the solution through a pipette.

i You may need to incubate the SDS-pellet mixture at +50°C to completely solubilize the protein.

- 7 Sediment any insoluble material by centrifuging the resuspended protein at 10,000 × *g* for 10 minutes at +2 to +8°C.

- 8 Transfer the supernatant containing the protein to a new tube.

i The sample may be used immediately or may be stored at –15 to –25°C. The protein solution may be analyzed directly on a western blot.

Alternate Protocol 5

To recover proteins more efficiently, perform the following steps:

- 1 Dialyze the phenol-ethanol supernatant against 3 changes of 0.1% SDS at +2 to +8°C.

- 2 Centrifuge the dialyzed material at 10,000 × *g* for 10 minutes at +2 to +8°C.

- 3 Transfer the supernatant containing the protein to a new tube.

i The sample may be used immediately or may be stored at –15 to –25°C. The protein solution may be analyzed directly on a western blot.

3. Results

Expected yield of RNA

- **Tissues:** liver or spleen, 6 to 10 µg/mg tissue; kidney, 3 to 4 µg/mg tissue; skeletal muscle or brain, 1 to 1.5 µg/mg tissue; placenta, 1 to 4 µg/mg tissue.
- **Cultured cells:** epithelial cells, 8 to 15 µg/10⁶ cells; fibroblasts, 5 to 7 µg/10⁶ cells.

Expected yield of DNA

- **Tissues:** liver or kidney, 3 to 4 µg/mg tissue; skeletal muscle, brain, or placenta, 2 to 3 µg/mg tissue.
- **Cultured cells:** human, rat, or mouse cells, 5 to 7 µg/10⁶ cells.

A_{260}/A_{280} ratio

- **RNA:** Total RNA isolated with TriPure Isolation Reagent is undegraded and free of protein and DNA. The isolated RNA has an A_{260}/A_{280} ratio of 1.6 to 2.0. It is suitable for northern blot analysis, poly(A)+ selection, *in vitro* translation, RNase protection assays, cloning, or RT-PCR.
- **DNA:** DNA isolated with TriPure Isolation Reagent is free of RNA and proteins with an A_{260}/A_{280} ratio of >1.7. It is suitable for PCR, restriction digests, or Southern blots.
- **Protein:** Protein isolated with TriPure Isolation Reagent is suitable for western blots.

4. Troubleshooting

Protocol	Observation	Possible cause
RNA isolation	Low yield, A_{260}/A_{280} ratio < 1.65	Incomplete sample homogenization or lysis.
		Incomplete solubilization of final RNA pellet.
		Too little TriPure Isolation Reagent used for sample homogenization.
		Following homogenization, sample was not stored for 5 minutes at +15 to +25°C.
		Contamination of the aqueous phase with phenol phase.
	RNA is degraded.	Incomplete solubilization of final RNA pellet.
		Tissues were not immediately processed or frozen after removal from animal.
		Samples used for the isolation procedure were stored at –15 to –25°C instead of at –60°C or below.
		Cells were dispersed by trypsin digestion.
		Aqueous solutions or tubes not RNase-free.
DNA contamination	Too little TriPure Isolation Reagent used for sample homogenization.	
	Sample used for the procedure contained organic solvents, such as ethanol, DMSO, or strong buffers, or had an alkaline pH.	
DNA isolation	Low yield	Incomplete sample homogenization or lysis.
		Incomplete solubilization of final DNA pellet.
	A_{260}/A_{280} ratio < 1.7	Incomplete removal of phenol from the DNA preparation with sodium citrate/ethanol.
	DNA is degraded.	Tissues were not immediately processed or frozen after removal from animal.
		Samples used for the isolation procedure were stored at –15 to –25°C instead of at –60°C or below.
		Sample was homogenized with a Polytron or other high speed homogenizer.
RNA contamination	Too much aqueous phase remained with the interphase and organic phase during Protocol 4.	
Protein isolation	Low yield	Inadequate wash of the DNA pellet with 10% ethanol/0.1 M sodium citrate solution.
		Incomplete sample homogenization or lysis.
	Protein degradation	Incomplete solubilization of final protein pellet.
		Tissues were not immediately processed or frozen after removal from animal.
		Band deformation in PAGE.

5. Additional Information on this Product

5.1. Test Principle

Isolation procedures

- ① During sample homogenization or lysis, the TriPure Isolation Reagent disrupts cells and denatures endogenous nucleases, thus preserving the integrity of RNA and DNA in the sample.

- ② After chloroform is added to the extract, the entire mix is centrifuged.
 - The solution contains three phases: a colorless aqueous upper phase, a white interphase, and a red organic lower phase.

- ③ The upper phase is placed in a separate tube from the other two phases.

- ④ RNA is recovered from the colorless aqueous phase by isopropanol precipitation.

- ⑤ DNA and protein are sequentially isolated from the white interphase and the red organic phase by alcohol precipitation steps.

6. Supplementary Information

6.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols

i *Information Note: Additional information about the current topic or procedure.*

⚠ Important Note: Information critical to the success of the current procedure or use of the product.

① ② ③ etc. Stages in a process that usually occur in the order listed.

① ② ③ etc. Steps in a procedure that must be performed in the order listed.

* (Asterisk) The Asterisk denotes a product available from Roche Diagnostics.

6.2. Changes to previous version

Layout changes.

Editorial changes.

6.3. Ordering Information

Product	Pack Size	Cat. No.
Reagents, kits		
Glycogen, Molecular Biology Grade	custom fill	10 899 232 103
Sodium Dodecyl Sulfate (SDS)	1 kg	11 667 289 001

6.4. Trademarks

All product names and trademarks are the property of their respective owners.

6.5. License Disclaimer

For patent license limitations for individual products please refer to:

List of biochemical reagent products.

6.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

6.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

6.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site.**

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.

