

comprehensive solutions for your vaccine platforms

Process Development & Manufacturing

Millipore®

Expert Pharm/BioPharm Products & CTDMO Services

SAFC

Pharma & Biopharma Raw Material Solutions

BioReliance®

Pharma & Biopharma Manufacturing & **Testing Services**



MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

Our Mission Vaccines, empowered

You want to scale and deliver your vaccine to the world quickly. Our collaborative global vaccine capabilities will take your innovation from pre-clinical to full-scale GMP-manufacturing efficiently, safely, and cost-effectively.

Discover how expertise—empowered by collaboration—can refine and overcome manufacturing challenges in all vaccine modalities/platforms.

Our Portfolio of Brands

MilliporeSigma has brought together the world's leading Life Science brands, so whatever your life science problem, you can benefit from our expert products and services.

Millipore_®

The Millipore® portfolio from MilliporeSigma offers an ecosystem of industry-leading products and services, spanning preparation, separation, filtration and monitoring – all of which are deeply rooted in quality, reliability and time-tested processes. Our proven products, regulatory and application expertise are a strong foundation you can rely on to consistently perform at the highest level.

SAFC®

The SAFC® portfolio from MilliporeSigma offers customized and ready-to-use raw material solutions, backed by deep regulatory expertise. Our high-quality products and services are supported by an experienced and responsive team of raw material and regulatory experts who are committed to understand your requirements and provide tailored solutions that meet your exact needs.

BioReliance_®

The BioReliance® portfolio from MilliporeSigma encompasses biopharmaceutical characterization, safety testing and process development, as well as clinical and commercial biomanufacturing. Our experienced teams and operational expertise make us the partner who supports you all the way and always has your vital goal in mind.



Comprehensive Products & Services Portfolio

Key Capabilities to Increase Efficiency and Safety of Vaccines

This document describes a range of products most commonly used for vaccine process development and manufacturing. Because each vaccine workflow is unique, however, selection of products most appropriate for your specific workflow requires testing and confirmation through MSAT support.

steps

MAKE

Produce antigen

PURIFY

Remove cell debris, virus and other impurities

FORMULATE

Final drug product























Media preparation

Cell culture & propagation

(Inactivation)

DNA digestion

UF/DF Chromatography **UF/DF** Formulation

Sterile **Filtration**

Final Fill

example products



Durapore® filters



Mobius® mixers



Mobius® bioreactors



Emprove® chemicals



Millistak+® depth filters



Benzonase® endonuclease



Pellicon® cassettes Pellicon® capsules





Eshmuno® and Fractogel® resins Natrix[®] Q

membranes



Millipak[®] final fill filters



Mobius® single use assemblies

Testing Services

CDTMO Services Regulatory **Expertise**

Inactivated/Live-attenuated Vaccines

An **inactivated vaccine** (or killed vaccine) consists of **virus** particles that have been grown in culture and lose disease producing capacity but remain capable of raising immunogenicity.

2

The virus is either in the supernatant or cell lysis is required, followed by Benzonase® endonuclease. Formaldehyde/beta propiono lactone is used for inactivation, followed by purification of the virus using chromatography and TFF.

Process Challenges:

- 1. Virus purity
- 2. HCP, DNA clearance
- 3. Virus recovery




Examples of live-attenuated vaccine: VARIVAX® (MSD); examples of inactivated vaccines including: BBIBP-CorV (Sinopharm), Fluzone® (Sanofi).

Virus propagation is cell-culture based and common cell lines used are VERO, Per.c6, MDCK.


Production & Purification of Virus Based Vaccines

Inactivated/Live-attenuated Production Platform




Assurance & Regulatory Compliance

Virus-Based Vaccines (Cell based) - 1000 L Single-Use



Virus-Based Vaccines (Cell based) – 1000 L Hybrid*



Subunit Vaccines

Subunit vaccines contain fragments of protein and/or polysaccharide from the pathogen, which are likely to produce a strong and effective immune response. By restricting the immune system's access to the pathogen in this way,

Depending on the characteristics of the protein to be expressed, cell lines such as mammalian, yeast, bacteria or insect cells can be used.

the risk of side effects is minimized.

2


Following production of the antigen, purification includes clarification, IEX and/or HIC chromatography and TFF followed by sterile filtration.

3

Examples of subunit vaccines include:
Acellular pertussis, Meningitis, Shingles,
cell-based Influenza, and many COVID
vaccines under development by some top
sellers including FluBlok® (Sanofi), Shingrix®
(GSK), Nuvaxovid (Novavax), ZifiVax™
(Longcom), etc.

Production & Purification of Protein Subunit Vaccines


Insect Cell Expression System



Assurance & Regulatory Compliance


Production & Purification of Protein Subunit Vaccines

CHO Cell Expression System




Assurance & Regulatory Compliance

Protein Subunit Vaccines (CHO Cell-Based) - 1000 L Single-Use



Protein Subunit Vaccines (CHO Cell-Based) – 1000 L Hybrid*



Millipore_®

Virus-Like-Particle Vaccines (VLP)

VLPs contain repetitive, high density displays of viral surface proteins that present conformational viral epitopes able to elicit strong cell and B cell immune responses.

VLPs self assemble into structures morphologically resembling viruses.

Since VLPs cannot replicate, they provide a safer vaccine methodology and can be produced in a variety of cell culture systems (yeast, insect, bacteria) against different strains.

2

After VLP production, cells are harvested and separated into a pellet prior to being disrupted/ lysed for the release of VLPs. The lysate is subsequently clarified. Further purification involves the use of UF/DF, chromatography, inactivation/removal VLP, final formulation and final sterilizing filtration.



Examples of VLP vaccines include: Engerix® (Hepatitis B from GSK), GARDASIL®9 (HPV from MSD)

Production & Purification of VLP-Based Vaccines

Insect Cell Expression System



Media and Inoculum **Preparation**



Cell growth in **Bioreactor and Virus Inoculation**



Nucleic Acid Removal

Prostak™

 $0.65 \, \mu m$

• Pellicon® 2

0.1-0.65 µm

cassettes with

Microfiltration

Modules (0.1 -



Clarification



• Pellicon® 2

Biomax® /

300 kDa

Pellicon[®]

Ultracel®

membranes

capsules with

100-300 kDa

cassettes with

Ultracel® 100-

Bioburden Reduction



UF/DF

Purification Chromatography



Baculovirus



Polishina Inactivation Chromatography



UF/DF



Formulation and Sterile **Filtration**

- Durapore® 0.1 or 0.22 µm : Mobius® Mixing
- Millipore Express® SHR 0.1 μm
- Viresolve® barrier
- Single-use and filter validation services

- Sf-RVN® Insect Cell Line
- EX-CELL® CD Insect Cell Media
- Upstream Chemicals and Supplements
- Raw materials testing
- Cell bank & virus seed production, testing & storage
- In-process characterization &

- systems
- NovaSeptum[®] GO Sterile Sampling Systems

Benzonase[®]

Benzonase®

ELISA kit II

Benzonase[®]

In-process

characteriza-

tion & testing

Residual

- Durapore® membranes
 - Polygard® CN, CR
 - Clarigard® filters Clarisolve[®] depth
 - filters
 - Millistak+® CE, HC/HC Pro
 - Single-use and filter validation services
- endonuclease

 - Biosafety lot of bulk harvest
 - Biosafety testing of production
 - In-process characterization &

- Millipore Express[®] SHC and Millipore Express® SHF
- Durapore® 0.22 μm or 0.45 μm filters
- Milligard® PES filters

In-process

characteriza-

tion & testing

- membranes Single-use and filter validation services
- Emprove® Chemicals: Buffers, Salts, Aminoacids, Isotonicity Agents, Stabilizers
- Cleaning-In-Place solutions
- In-process ization &

- Fractogel® TMAE/DEAE/ DMAÉ Eshmuno[®] O
- resins Natrix® Q chromatography
- membranes
- Emprove® Chemicals: Buffers, Salts, Aminoacids, Isotonicity Agents, Stabilizers
- Cleaning-In-Place solutions
- Viral Clearance
- characterization & testing

• Mobius® mixing systems

• Emprove®

Viral

Clearance

characteriza-

Studies

tion & t

esting

Chemicals:

Surfactants

(e.g. PS-80,

- NovaSeptum® GO Sterile Eshmuno[®] O Sampling resins Systems
 - Natrix[®] Q chromatography membranes

• Fractogel®

DMAÉ

TMAE/DEAE/

- Emprove®
 - Chemicals: Buffers, Salts, Aminoacids, Isotonicity Agents, Stabilizers
 - Cleaning-In-Place solutions
 - In-process characterization & testing


- Pellicon® 2 cassettes with Biomax® / Ultracel® 100-300 kDa membranes
- Pellicon[®] capsules with Ultracel® 100-300 kDa membranes
- Single-use and filter validation services
- Emprove® Chemicals: Buffers, Salts, Aminoacids, Isotonicity Agents, Stabilizers
- Cleaning-In-Place solutions
- Biosafety, Analyticál and Stability testing of substance
- In-process characterization & testing

- Durapore® filters
- Millipore Express® SHC filters
- Millipak® Final Fill filters
- Single-use and filter validation services
- Emprove® Chemicals: Buffers, Stabilizers, pH, Surfactant, Isotonicity Agents, Adjuvants
- Analytical testing of final fill drug product
- In-process characterization & testing


Assurance & Regulatory Compliance

Millipore SAFC BioReliance

Virus-Like Particle Vaccines (Insect Cell-Based) – 1000 L Single-Use



Virus-Like Particle Vaccines (Insect Cell-Based) – 1000 L Hybrid*



Millipore_®

Viral Vector Vaccines – Adenovirus Vaccines

1

A live vector vaccine is a vaccine that uses a weakened or harmless microorganism to transport pieces of antigen to stimulate an immune response. Common viral vectors are adenovirus, canarypox, lentivirus, and alphavirus.

The vectors deliver genetic materials into human cells, and instruct them to make antigen, which then trigger an immune response. Viral vectors are genetically modified to have non-replicating viral vectors (e.g. adenovirus) and replicating viral vectors (e.g. weakened Measles).

2


Typical viral vectors purification process included cell harvest followed by cell lysis and nucleic acid treatment. Subsequent steps including UF/DF, chromatography, final sterilizing filtration may or may not be used, depending if the size of vector is too large for filtration. In this case, closed or aseptic manufacturing practices must be followed for large vectors.

3

Viral vector vaccines include: ERVEBO® Ebola Vaccine (MSD), Vaxzevria™ (Oxford/AstraZeneca)


Production & Purification of Viral Vector Vaccines

Adenovirus Production Platform




Assurance & Regulatory Compliance

Adenovirus Based Vectored Vaccines (Cell Culture-Based) – 2000 L Single-Use



Millipore SAFC

Adenovirus Based Vectored Vaccines (Cell Culture-Based) – 2000 L Hybrid*



*Mix of Single-Use & Multi-Use systems

Plasmid DNA (pDNA) Vaccines

Plasmid DNA vaccines are based on purified plasmid preparations containing one or more DNA sequences capable of inducing and/or promoting an immune response against a pathogen.

E. coli-based, the purified plasmid is delivered by lipofection or gene gun. pDNA enters the nucleus of transfected local cells and plasmid-encoded genes are expressed. Foreign antigens are generated and trigger an immune response by vaccinated host.

2


After harvest, cells are lysed and clarified. Downstream purification includes TFF, chromatography, and final sterile filtration.



Examples of pDNA vaccines include: COVID vaccines from ZyCoV-D[®].


Production & Purification of pDNA Vaccines

pDNA Production Platform




Assurance & Regulatory Compliance

Plasmid DNA Vaccines (*E. coli*) – 500 L Single-Use (excl. Fermenter)



Plasmid DNA Vaccines (E. coli) – 500 L Hybrid*



*Mix of Single-Use & Multi-Use systems

mRNA Vaccines

Delivery of an mRNA into the cytosol of a cell can induce the production of a target protein which can function as an antigen to trigger an immune response for vaccination purposes

E. coli-based, the target protein expressing gene is inserted into the plasmid DNA and the plasmid DNA acts as a precursor to the mRNA.

2


The DNA is linearized, purified and mRNA is produced through *in-vitro* transcription using enzymes. The mRNA is purified, end capped and formulated into lipids.

3

Examples of mRNA vaccines include: COVID vaccines from Spikevax $^{\text{TM}}$ (Moderna) and Comirnaty $^{\text{®}}$ (Pfizer-BioNTech).


Production & Purification of mRNA Vaccines

mRNA production platform




Assurance & Regulatory Compliance

mRNA Based Vaccines - 200 L Single-Use



mRNA Based Vaccines – 200 L Hybrid*



BioReliance® Biosafety Testing Services

Patient safety and compliance

For vaccine and viral therapies, BioReliance® biosafety testing services provide a comprehensive range of assays and services to support each stage of product development. From cell banking and cell line characterization to product characterization and lot release testing, our GMP-compliant testing services and regulatory expertise can help progress your vaccine from discovery to commercialization.


BioReliance® upstream services include:

- MCB/WCB bank manufacturing
- MCB/WCB bank characterization
- Biorepository services
- · Virus bank manufacturing
- Virus bank characterization
- Next generation sequencing
- Raw material testing

BioReliance® downstream services include:

- Analytical services for biologics
- Viral clearance studies
- Bulk lot release testing
- Final product release testing

Our Integrated mRNA-LNP CDMO offering





mRNA

- Preclinical Custom mRNA Manufacturing
- GMP Clinical & Commercial Custom mRNA Manufacturina*

Lipids

· Portfolio and Custom GMP Lipids

LNP/Liposomes

- Development
- cGMP Production
- Final Fill



^{*} Availability H2 2023. Timeline Subject to change - CoE: Center of Excellence

Emprove® Chemicals

Ensure the successful process and formulation of your product


Specifically developed for high-risk applications, our buffers, salts and stabilizers are low in bioburden and endotoxins. Supported by our Emprove® program, our excipients have extensive documentation, helping minimize regulatory and quality-associated risks in your manufacturing. All this simplifies the complexity of your supplier qualification and speeds up processes, thus reducing total cost of ownership.

Our portfolio of more of more than 400 raw and starting material comes with comprehensive and thorough documentation not only covering the latest regulatory requirements, but also anticipating industry expectations not yet covered by regulation.

- Elemental impurity information according to ICH Q3D
- Multi-compendial raw materials facilitate the international drug registration procedures
- Paper-free packaging of our excipients to minimize risk


Emprove® Chemicals include:

- Buffering Agents
- Solvents & Surfactants
- Emprove® Salts & Process Specialities
- Emprove® Stabilizers & Amino-acids





M Lab™ Collaboration Centers



Expertise

M Lab™ Collaboration Centers uniquely provide the pharma/biopharma industry access to scientific, technical, and engineering expertise in state-of-the-art, fully equipped facilities so they can explore ideas, learn innovative techniques, and collaborate – whether in person or remotely – to find solutions to their toughest development and manufacturing problems.

The M Lab[™] Collaboration Centers global network is comprised of **300+technical experts** ranging from process development scientists to bioprocess engineers supporting a variety of modalities.

These highly qualified experts are at your disposal to **troubleshoot your process, identify efficiencies**, and **adopt innovative techniques** to help you bring life-enhancing drug therapies to market faster and more efficiently.

Our technical experts have helped biopharma/pharma manufacturers solve over 5,000 problems and save 13,000+ hours annually in process troubleshooting and deviation investigations.

Industry innovators collaborate with our experts to troubleshoot existing processes, overcome barriers to single-use implementation, receive guidance for process development, learn best practices for adopting biopharma 4.0 technologies, explore applications for novel modalities, and much more.

M Lab™ Collaboration Centers are beacons of scientific leadership, producing 150+ peer-reviewed articles, technical presentations, and patent filings each year.





vaccines, empowered

To place an order or receive technical assistance

In the U.S. and Canada, call toll-free 1-800-645-5476

For other countries across Europe and the world, please visit: **EMDMillipore.com/offices**

For Technical Service, please visit: **EMDMillipore.com/techservice**

MilliporeSigma 400 Summit Drive Burlington, MA 01803





