



New Generation of Organic CRM



 ISO/IEC 17025 accreditation of HYDRANAL® service lab

SIGMA-ALDRICH®

New Generation of Organic CRM



Michael Weber Manager R&D/Innovation Europe

Dear Colleague,

It has now been sixty years since Fluka® was founded in 1950 in the middle of St. Gallen, Switzerland. After the addition of new facilities in Buchs in the Rhine Valley, Fluka established itself as one of the leading producers and suppliers of research chemicals, reagents and standards. In 1989, Fluka joined the Sigma-Aldrich Corporation and eventually changed its name to Sigma-Aldrich Switzerland. Today, within the impressive Sigma-Aldrich product portfolio, Fluka is the brand representing high-quality reagents and standards for analytical chemistry.

Hence, it is a pleasure to present a new service tool as well as a new class of certified standards which are now launched under the Fluka brand. With the introduction of the new webbased dynamic Custom Standards Platform (*sigma-aldrich.com/csp*) we provide the capability to compose your individual inorganic custom standard within minutes. You can actually choose between multi-elemental standards for spectroscopy as well as anion and cation mixtures for ion chromatography. The whole process from specifying the composition of your individual standard to the submission of your request is a fast and easy online experience. And of course, all these custom standards are produced and certified under double accreditation following both ISO/IEC 17025 and ISO Guide 34. Try it!

As a second highlight, we are introducing the first series of a new generation of organic Certified Reference Materials for calibration and validation of chromatography or for other analytical applications. After the implementation of high-precision quantification with 1H-NMR in Buchs, our lab received ISO/IEC 17025 accreditation for qNMR analysis. As a relative primary method, qNMR offers outstanding opportunities for the exact determination of the purity of organic substances. The new class of organic *Trace***CERT**[®] standards is traceable to NIST references and will be characterized not only by qNMR, but also by additional analytical techniques including HPLC, GC and other methods. Amino acids and PAHs have been selected to be the first CRM which are certified by this multi-step approach, but many more substances of high relevance to the analytical community will follow soon.

Happy birthday, Fluka brand, with 60 years of exceptional service, meeting the needs of the world-wide scientific community!

Tidoel Welles-

Michael Weber Manager R&D/Innovation Europe michael.weber@sial.com

2

Analytix is published five times per year by Sigma-Aldrich Chemie GmbH, MarCom Europe, Industriestrasse 25, CH-9471 Buchs SG, Switzerland Publisher: Sigma-Aldrich Marketing Communications Europe Publication director: Ingo Haag, PhD Editor: Daniel Vogler

Feature Article

Launch of a New Generation of Organic
 Certified Reference Materials
 First Series of Products Includes 21 Amino Acids

Microbiology

7 Overview of Chromogenic Media Chromogenic media offer a range of benefits for the enumeration, detection, and identification of microorganisms.

Standards

- 10 New Dynamic Platform for Inorganic Custom Standards
- 12 Certified Reference Materials for Spectroscopy

Chromatography

- 13 Enantiomeric Purities of Pharmaceuticals Using Carbohydrate-based Isothiocyanates
- 18 Analysis of Bioethanol Detection of Chloride and Sulfate with Dication and Trication Solutions

Titration

20 Increasing Reliability and Value of Analytical Results

ISO/IEC 17025 accreditation of HYDRANAL® service lab

Monthly Savings Program

22 35% Off – Volumetric Titration Reagents

New Product Corner

- 23 Analytical Standards of Alkylresorcinols
- 23 New Standards for Regulated Food Dyes

Launch of a New Generation of Organic Certified Reference Materials

First Series of Products Includes 21 Amino Acids

Alexander Rück, Senior Scientist R&D Europe alexander.rueck@sial.com Christine Hellriegel, Senior Scientist R&D Europe christine.hellriegel@sial.com Matthias Nold, Product Manager Analytical Standards matthias.nold@sial.com

We are pleased to present a first series of the new product line of organic *Trace***CERT**[®] certified standards. These certified reference materials (CRMs) are produced and certified in a double accredited laboratory fulfilling both ISO/IEC 17025 and ISO Guide 34, using quantitative NMR for content determination with direct traceability to NIST references and the SI base unit kg. The double accreditation represents the highest achievable level of reliability and therefore earns the label of the "gold standard" among CRM producers.

The new product class of organic CRMs for use as internal standards for qNMR has been presented in the previous issue of Analytix. Here, we introduce our first product line of organic CRMs for use as chromatography standards.

Since the use of certified reference materials is becoming increasingly important, especially in the areas of food & beverage analysis as well as in environmental analysis, the first two series of CRMs include amino acids (which are presented here) and polyaromatic hydrocarbons (PAH) (which we will present in the next issue of Analytix).

Why certified reference materials?

To ensure the reliability of analytical results, the quality of the reference material used for calibration of the instruments is crucial. Certified reference materials are most trustworthy since they fulfil several requirements:

- They must be traceable to an internationally accepted reference (e.g. from NIST, BAM or SI)
- The certified value must have a properly calculated uncertainty
- They must be shown to be sufficiently homogeneous
- The stability of the material must be carefully determined

If a manufacturer is accredited to ISO Guide 34 ("General requirements for the competence of reference material producers") and also according to ISO/IEC 17025 ("General requirements for the competence of testing and calibration laboratories"), then the fulfillment of all these requirements is guaranteed.

For more than two years now, the Buchs site of Sigma-Aldrich in Switzerland has held this "gold standard" for the production of inorganic *TraceCERT* standard solutions (*sigma-aldrich.com/inorganiccrm*) for AAS, ICP and IC. While for inorganic materials internationally accepted reference materials are readily available (e.g. from NIST and other metrological institutes), it is a significant challenge to establish traceability for organic reference materials.

Conventionally, the purity of a substance is determined by means of GC or HPLC and the content is calculated by subtracting water content, residue solvents and inorganic impurities from the measured chromatographic purity. However, since every chemical substance has its own characteristic UV-absorption behaviour, direct traceability is only possible if an internationally accepted reference standard is available for the compound in question. This is not necessary if a relative primary analytical method such as quantitative NMR (qNMR) is used. With qNMR technology, traceability between completely different organic compounds can be achieved since the NMR signal intensity is in direct relation to the numbers of protons involved.



Double accreditation for CRM production using qNMR

In the R&D lab of Sigma-Aldrich Buchs (Switzerland) a new Bruker Avance III 600 MHz Ultrashield NMR spectrometer was recently installed and the instrument performance was optimized for the high-resolution quantification of organic substances. By the end of 2009, the new NMR lab was fully accredited under both ISO/IEC 17025 and ISO Guide 34 for the certification of organic reference materials using ¹H qNMR. Quantitative NMR is a non-destructive technique that demonstrates many advantages over other analytical methods: it requires easy sample preparation; analysis is fast; it provides full structural information; data about impurities or isomers are available; and it is less time consuming than many other analytical techniques, since it requires no equilibration time.



¹H qNMR is marked by its linearity, specificity and robustness regarding most parameter settings. The most outstanding attribute of ¹H qNMR is that it may be considered a primary ratio method because it does not need a standard reference of the same material. The NMR spectrum of a solution shows all soluble organic substances present, with the peak area being directly proportional to the amount of nuclei being measured.

Hence, the signal ratio of two different protons can be measured with tremendous precision and the only significant contribution to the measurement uncertainty

is the integration of the signals. In other words, the direct response of a qNMR experiment is of the highest trueness. In the previous issue of Analytix, we already presented CRMs which have been tailor-made as qNMR calibrants. In the following, we will describe the certification process, exemplified by amino acids, our first series of organic CRMs.

Certification of amino acids by qNMR

In the certification process, the following tests have been performed for any candidate substance: hygroscopy/volatility; choice of the appropriate deuterated solvent; chemical stability; 2D-NMR (H-H COSY) to detect potential underlying impurity signals and determination of the ¹H relaxation time. Once having selected a suitable starting material, a defined number of subsamples are taken from different positions within the bulk material. From those samples, some are stored at room temperature for periodical stability tests in order to get data on aging of the material. Other samples are stored at higher temperature for stress tests, helping to ensure stability at even extreme transportation conditions. Further samples go into quality control for chromatography, LC-MS or other standard techniques if required. However, 10 subsamples are taken for homogeneity assessment and content determination by gNMR. Therefore, it is necessary that a reference standard has certain required properties in order to be suitable for the quantification. Since the reference standard and the amino acid are weighed together into the same vial, both substances have to be soluble in the same deuterated solvent, but must not react with each other or with the solvent. Each of the substances must show at least one ¹H-NMR-signal that is not overlapped, and the content of the reference standard must be traceable to NIST or SI. The weighing procedure is performed on an ultramicro balance in such amounts that equal NMR signal intensities will result. Taking into account air buoyancy, the integrals of 10 subsamples contribute to the resulting content value. This is followed by a comprehensive uncertainty evaluation, resulting finally in a certificate associated with each amino acid which conforms with ISO Guide 31. A complete series of the 20 proteinogenic amino acids, along with Cystin, have been guantified by gNMR using potassium hydrogen phthalate (KHP) or maleic acid, respectively, as the internal reference standard. Figure 1 shows the gNMR-spectrum of L-Serine with the internal standard KHP as an example. An illustration of the certification process is shown in Figure 2.



(continued on page 6)



Figure 2 Multi-component certification approach

Table 1 lists the organic *Trace***CERT**® products for the proteinogenic amino acids and cystin. The products are supplied with a printed certificate containing certified value, comprehensive documentation, proper uncertainty calculation, lot-specific values, expiration date and storage information. On our website, customers can download complete example certificates as pdfs.

The amino acid CRMs described in this article are the first products launched by this innovative approach to highaccuracy qNMR quantification, with an exciting expanded portfolio to follow. As a next group we will launch CRMs of polyaromatic hydrocarbons (PAHs). Please visit our website at *sigma-aldrich.com/organicCRM* for an up-to-date product listing.

Cat. No.	Description	Pack Size
44526	L-Alanine	1 g
90538	L-Arginine hydrochloride	1 g
51363	L-Asparagine	1 g
51572	L-Aspartic Acid	1 g
95437	L-Cysteine	1 g
95436	L-Glutamic acid	1 g
76523	L-Glutamine	1 g
76524	Glycine	1 g
73767	L-Histidine	1 g
56241	L-Isoleucine	1 g
76526	L-Leucine	1 g
67448	L-Lysine hydrochloride	1 g
39496	L-Methionine	1 g
40541	L-Phenylalanine	1 g
93693	L-Proline	1 g
54763	L-Serine	1 g
61506	L-Threonine	1 g
51145	L-Tryptophan	1 g
91515	L-Tyrosin	1 g
50848	L-Valine	1 g
49603	L-Cystine	1 g

Table 1 Fluka labelled certified amino acid TraceCERT standards

Microbiology

Overview of Chromogenic Media

Chromogenic media offer a range of benefits for the enumeration, detection, and identification of microorganisms.

Jvo Siegrist, Product Manager Microbiology ivo.siegrist@sial.com

The use of traditional versus improved media formulation containing chromogenic substrate is currently an important topic in the field of microbiology. The focus behind such developments was to produce media that would make the detection and identification of microorganisms more rapid and more reliable. Chromogenic substrates such as ONPG, X-Gal, or X-Glu, together with a specified selectivity of the medium, is the simple principle behind chromogenic media. The target organisms are characterized by enzyme systems that metabolize the substrates to release the chromogen. The chromogen can then be visually detected by direct observation of a distinct colour change in the medium. Direct confirmation of the target organism without further testing is sometimes possible.

Did you know?

Sigma-Aldrich provides thirty-seven different chromogenic media.



Over the past 15 years, we have continued to develop and add new media to our ever expanding product range designed to meet customer needs

Figure 1 HiCrome[™] Rapid Coliform Broth

5			
Species	Enzyme	Substrate	Selective Agents
Bacillus cereus	β-glucosidase, Phosphatidylinositol-specific Phospholipase C	indoxyl-β-glucopyranoside, indoxyl-myo-inositol-1-phosphate	polymyxin B
Campylobacter	na	na	deoxycholate, cefoperazone, amphotericin B
Candida	β-acetylgalactosaminidase, alkaline phosphatase	indoxyl-N-acetyl-β-D-glucosaminide, indoxyl-phosphate	Chloramphenicol, Gentamicin
Clostridium perfingens	β-glucosidase (plus sucrose fermentation)	indoxyl-β-D-glucoside	D-cycloserine, polymixin B
Coliforms/E. coli	β-glucuronidase, β-galactosidase	indoxyl-β-glucuronide, Indoxyl-β-galactoside	bile salts, tergitol 7°, SDS, novobiocin, cefsulodin
Cronobacter (E. sakazakii)	a-glucosidase	indoxyl-a-D-glucoside	deoxycholate, crystal violet, sodium thiosulfate
E. coli O157	β -glucosidase, α -galactosidase	indoxyl-β-D-glucuronide, indoxyl-α-galatoside	bile salts, SDS, crystal violet, potassium tellurite, novobiocin, cefixime
Enterococci	β-D-glucosidase	indoxyl-β-glucoside	sodium azide, polysorbate 80
Extended Spectrum &-Lactamase Enterobacteria (ESBL)	β-D-glucosidase	indoxyl-β-glucoside	cefpodoxime, cefotaxime, ceftazidime
Klebsiella	β-D-ribofuranosidase, β-D-glucosidase	indoxyl-β-D-ribofuranoside, indoxyl-β-D-glucoside	bile salts, SDS, carbenicillin
Listeria spp.	β- glucosidase	indoxyl-β-glucoside	lithium chloride, ceftazidime, amphotericin B, nalidixic acid, polymyxin B
L. monoctogenes	Phosphatidylinositol-specific Phospholipase C, β-glucosidase,	indoxyl-β-glucoside, indoxyl-myo-inositol-1-phosphate	lithium chloride, ceftazidime, amphotericin B, nalidixic acid, polymyxin B
Pseudomonas	β-Alanyl arylamidase	7-Amido-1-pentyl-phenoxazin-3-one	cetrimide
Salmonella	α-galactosidase, lipase	indoxyl-α-galactoside, indoxyl-fatty acid ester	sodium deoxycholate
MRSA (Methicillin Resistant Staphylococcus aureus)	α-glucosidase	indoxyl-a-D-glucopyranoside	methicillin, high concentration of sodium chloride
Staphylococcus aureus	α-glucosidase, phosphatase, deoxyribonuclease	indoxyl-α,D-glucoside, phenolphthalein phosphate, indoxyl-phosphate	tellurite, lithium chloride
Streptococci	β-glucuronidase	indoxyl-β-glucuronide	sodium azide
UTI (Urinary Tract Infections)	β -glucosidase, β -galactosidase	indoxyl-β-glucopyranoside, indoxyl-β-galactoside	-
Vibrio	β -glucosidase, β -galactosidase	indoxyl-β-glucopside, indoxyl-β-galactoside, indoxyl-β-galactoside	high concentration of sodium chloride, sodium thiosulphate, sodium citrate, sodium cholate
VRE (Vancomycin Resistant Enterococci)	α-glucosidase, β-glucosidase, β-galactosidase	indoxyl-α-glucopyranoside, indoxyl-β-glucopyranoside, indoxyl-β-galactoside	vancomycin
Yeasts and Moulds	R-N-acetylgalactosaminidase	indoxyl-N-acetyl-R-D-alucosaminide	ovytetracycline

β-xylosidase indoxyl-β-D-xyloside

Table 1 Summary of possible enzyme activities, chromogenic substrates and selectivity system for microorganisms

(continued on page 8)

Advantage of chromogenic media:

- Faster results (compared to traditional method)
- Reliable visual detection (often no further testing required)
- Additional testing possible directly from the media

Within recent years, great strides have been taken in the sector of chromogenic media. Initial research concentrated on the use of synthetic substrates for the detection of enzymatic microbial activities. Nitrophenol and nitroanilline compounds were used at this time, producing a yellow colouration. The colour of nitrophenol, however, is influenced by a pH-change, making it difficult to use reliably in microbiology. Later developments included the use of naphtol or naphtylamine. Today, while diverse modern chromogenic substrates are available, most of the modern substrates are based on the indoxylsubstrate. The use of different chromophore and metabolite derivatives then makes it possible to detect diverse enzyme activities all in one assay. The color of the indoxyl-substrates can be as follows: blue (5-bromo-4-chloro-3-indoxyl- = X, 3-indoxyl- =Y), magenta (5-bromo-6-chloro-3-indoxyl-), salmon (6-chloro-3-indoxyl-), purple (5-iodo-3-indoxyl-) and green (N-methylindoxyl-). One of the major advantages of the indoxyl-substrate and these chromophores is that they remain in the cell, making the characterization of a single cell possible (no diffusion into the media).

Additional advancements in the knowledge about enzyme and species specificity have also occurred within the past year. These recent gains in the development of selective agents and diverse chromogenic substrates have led toward an impressive range of chromogenic media available to meet our customers' unique analytical emphases (See **Tables 1** and **2**).

Organisms	Media		
Bacillus cereus	Fluka®	92325	HiCrome™ Bacillus Agar
Candida albicans	Fluka	94382	Candida Ident Agar, modified
Cl. perfringens	Fluka	12398	CP ChromoSelect Agar
	Fluka	75605	m-CP Agar
E. coli	Fluka	70722	HiCrome™ E. coli Agar B
	Fluka	09142	HiCrome [™] ECD Agar with MUG
	Fluka	92435	TBX Agar
E. coli & Coliforms	Fluka	81938	HiCrome™ Coliform Agar
	Fluka	73009	HiCrome™ ECC Agar
	Fluka	85927	HiCrome™ ECC Selective Agar
	Fluka	51489	HiCrome™ Rapid Coliform Broth
	Fluka	39734	Membrane Lactose Glucuronide Agar
E. coli 0157:H7	Fluka	39894	HiCrome™ EC O157 Agar
	Fluka	72557	HiCrome™ EC 0157:H7 Selective Agar, Base
	Fluka	80330	HiCrome [™] Enrichment Broth Base for EC O157:H7
	Fluka	83339	HiCrome™ Mac Conkey Sorbitol Agar
Thermotolerant E. coli	Fluka	90924	HiCrome™ m-TEC Agar
Enterobacter sakazakii (Cronobacter spp.)	Fluka	92324	HiCrome™ Cronobacter spp. Agar
	Fluka	14703	HiCrome™ Cronobacter spp. Agar, modified
Enterococci	Fluka	52441	HiCrome™ Enterococci Broth
	Fluka	51759	HiCrome™ Rapid Enterococci Agar
Enterococcus faecium	Fluka	90919	HiCrome™ Enterococcus faecium Agar Base
Klebsiella	Fluka	90925	HiCrome™ Klebsiella Selective Agar Base
Listeria	Fluka	53707	HiCrome™ Listeria Agar Base, modified
	Fluka	77408	Listeria mono Differential Agar (Base)
Proteus, enteropathogenic gram-positive organisms	Fluka	16636	HiCrome™ UTI Agar, modified
Salmonella	Fluka	00563	HiCrome™ MM Agar
	Fluka	90918	HiCrome™ RajHans Medium, Modified
	Fluka	78419	HiCrome™ Salmonella Agar
	Fluka	05538	HiCrome™ Salmonella Agar, Improved
	Fluka	84369	Salmonella Chromogen Agar
Staphylococcus aureus	Fluka	05662	HiCrome™ Aureus Agar Base
	Fluka	68879	Phenolphthalein Phosphate Agar
Methicillin-resistant Staph. aureus	Fluka	90923	HiCrome™ MeReSa Agar Base
Vibrio	Fluka	92323	HiCrome™ Vibrio Agar
Yeasts and fungi	Fluka	66481	HiCrome™ OGYE Agar Base

Table 2 Sigma-Aldrich's product line of chromogenic media according to organisms detected (Fluorogenic media are not listed. Complete product listings are available at *sigma-aldrich.com/microbiology*)





Figure 2 HiCrome[™] ECC Agar (Fluka 73009)

Figure 3 HiCrome[™] Coliform Agar (Fluka 81938)

The NEW CP ChromoSelect Agar

The European Directive on drinking water quality recommends mCP agar in the reference method for recovering C. perfringens. In the present study, three media (mCP, TSCF and CP ChromoSelect Agar) were evaluated for recovery of C. perfringens in different water samples. Out of 139 water samples tested, using a membrane filtration technique, 131 (94.2%) of the samples analyzed were found to be presumptively positive for C. perfringens on at least one of the culture media.

Green-colored colonies on CP *ChromoSelect* agar (CCP agar) were counted as presumptive C. perfringens isolates. Out of 483 green colonies on CCP agar, 96.9% (465 colonies, indole negative) were identified as C. perfringens, 15 colonies (3.1%) were indole positive and were identified as C. sordelli, C. bifermentans or C. tetani. Only 3 strains (0.6%) gave false positive results and were identified as C. fallax, C. botulinum, and C. tertium. Variance analysis of the data showed no statistically significant differences in the counts obtained between various media employed in this work.

The mCP method is very onerous for routine screening, and bacterial colonies could not be used for further biochemical testing. Conversely, the colonies on CCP and TSC were easy to count and subculture for confirmation tests. TSCF detects all sulfite-reducing clostridia, not only C. perfringens; however, in some cases, excessive blackening of the agar frustrated counting of the lower dilutions. If the contamination was too high, TSCF did not consistently produce black colonies and, as a consequence, the colonies' white color provided false negative results.

The identification of typical and atypical colonies isolated from all media demonstrated that CCP agar was the most specific medium for C. perfringens recovery in water samples.



Figure 4 CP ChromoSelect Agar (Cl. Perfringens appears as green colonies, Fluka 12398)

New Dynamic Platform for Inorganic Custom Standards

Matthias Nold, Product Manager Analytical Standards matthias.nold@sial.com Juerg Wuethrich, Senior Scientist R&D Europe juerg.wuethrich@sial.com

Under the label *Trace***CERT**[®] Sigma-Aldrich offers a wide range of Certified Reference Materials which are tailor-made for many common procedures and analytical techniques. However, every laboratory has its unique expertise, leading to many different applications. Customers who cannot find a suitable standard in our product portfolio can now use our new dynamic internet-platform for customized standard solutions.

Reliable quality and documentation

All certified custom standards mixes are certified according to **ISO/IEC 17025** and **ISO Guide 34**. This double accreditation is the highest attestation of quality and is also called the "Gold Standard for CRM Producers". The certified values

are traceable to internationally accepted references (SI unit, NIST SRM or other). Customers will get a printed lot-specific certificate according to ISO Guide 31. An example of the packaging and documentation is shown here. With our special packaging using pre-



cleaned HDPE bottles in coated aluminum bags, we can ensure maximum shelf life and low uncertainties of 0.2% to 0.6% relative, depending on element composition and concentration levels. The shelf life is normally between 6 months and 3 years.

Customer benefits and workflow

One of the key benefits of customized standard solutions is time savings: making individual multi analyte blends from single stock solutions or other raw materials is time consuming and costly.

A major advantage that our new dynamic web service provides for our customers is that, while clicking on the desired elements, the customer can directly see if this combination is chemically possible or if the analyte is stable in the chosen matrix.

After examining the technical feasibility and calculation of the production, customers will get a quotation within a short time, usually within two days. The target delivery time for a custom standard mix is less than four weeks.

The order volume is generally free from minimum requirements. Nevertheless, we recommend ordering more than one bottle of a custom standard in order to take advantage of the lower manufacturing cost per bottle. Especially for mixtures with longer shelf life, this can be a good opportunity to save money. From our long-standing expertise with the HDPE-bottles (100 mL, 250 mL or 500 mL), we can ensure the printed expiration date for each lot, since we know exactly the leaching and storage behavior of these containers.

There are sometimes limitations

Some aspects can lead to limitations regarding the analyte concentrations in certain mixtures.

When high concentrations of some of the analytes are required (more than 1000 mg/L per analyte), their solubility in the mixture can become a limitation. On the other hand, very low concentration levels can also cause difficulties. Particularly at the ppb level, many analytes are not stable in solution for an extended time.

Another serious issue often occurs when a very high concentration of one of the analytes is required in a mixture while all the other analytes should be present in a very low concentration. In this case, the impurity profile of the highly concentrated analyte imposes certain limits. For example, consider when a mix with 5000 mg/L calcium is required together with 10 mg/L strontium and 10 mg/L sodium. Assuming that the starting material is calcium carbonate and the maximum impurities for strontium and sodium shall not exceed 1% of the target concentration (10 mg/L), the maximum allowable strontium and sodium impurity levels in the calcium carbonate are then 8 ppm. Often it is not possible to get a calcium carbonate with such a high purity. As a consequence, the realization of custom standards with high variances of analyte concentrations in one mixture is sometimes not possible and/or the analyte concentrations have high uncertainties.

Using the platform is easy

On the landing page *sigma-aldrich.com/csp* we offer the possibility to choose between elemental standards (for AAS, ICP) and ion standards (for IC). Clicking on the corresponding links leads to an interactive periodic table (in case of the elemental standards) or to a table of anions and cations (in case of the ions standards).

Step 1: Defining matrix

Before the elements of the standard can be chosen, the matrix has to be defined by choosing between HCl and HNO₃ (**Figure 1**). It is necessary to choose the matrix first, because some analytes are not compatible with both matrices.

Standards



Figure 2 Picking components for the multi-element standards from the periodic table of elements.



Figure 1 Customers can select from two different matrixes for multi-element standards

Step 2: Defining elements or ions

The allowed elements for a specific matrix will then be marked green and the components for the standard can be picked by clicking on the green squares. The selected elements are then immediately marked in blue and a box for the definition of the concentrations appears on the right side of the periodic table (**Figure 2**).

For the multi-ion standards it works in the same way as for elements. One exception is that no pre-definition of the matrix is required since the matrix is automatically defined by the chosen analytes. For the anion standards, the matrix is water, whereas for cation standards, nitric acid (<0.1%) is used. With the first component defined, it is automatically defined whether a multi-anion or a multi-cation standard is composed. Thus, all cations will be disabled if you choose an anion and vice versa (**Figure 3**).

Bromide	Chloride	Chromate	Cyanide	Fluoride
lodide	Nitrate	Nitrite	Phosphate	Sulfate
Chlorite	Chlorate	Perchlorate	Bromate	Thiocyanate
Thiosulfate	Sulfite			
Ammonium	Barium	Calcium	Cadmium	Cobalt
Copper	Lead	Lithium	Magnesium	Manganese
Nickel	Potassium	Sodium	Strontium	Zinc
Ceslum				

Figure 3 List of anions and cations (in green) that can be chosen. This portfolio will be expanded continuously.

Step 3: Defining concentrations and package size

After all components are chosen for the custom standard, the concentrations for all of the components must be entered before clicking the "continue" button. After this, the user can specify whether the custom standard should be produced in mass per mass (mg/kg) or mass per volume (mg/L) units. Finally, the number of required bottles and/or the total volume can be entered. The choices for container volume are 100 mL, 250 mL, or 500 mL HDPL packages, or a freely selected bulk quantity (**Figure 4**).



Figure 4 Required information for a custom standard request

After all parameters have been defined, the request can be submitted online and a summary of the request can be printed out before sending. With the new Custom Standards Platform, Sigma-Aldrich is proud to offer an easy and fast way for the design of certified custom standards. Try it!

Matthias Nold, Product Manager Analytical Standards matthias.nold@sial.com

Optical Emission (OES) and X-Ray Fluorescence (XRF) Spectroscopy are efficient and fast techniques to determine the composition of metals, including trace elements as well as the major alloying ingredients. These techniques are therefore widely used for compositional control in the metal industry. The principle of OES is based on spark-discharge excitation of the solid sample, with measurement of the emitted light; this allows for quantification of the elements present. The principle for XRF is similar, although the excitation is provided by X-rays.

To ensure accurate measurement values, the instruments should be calibrated and monitored using one or more appropriate reference materials.

Sigma-Aldrich is launching a new product line of certified reference materials (CRMs) for spectroscopy. Certified values are given for the major alloy ingredients and the most rele-

vant impurities for each corresponding metal or alloy. The CRM samples (in the form of solid discs) are produced in the United Kingdom by MBH Analytical Ltd, a recognized manufacturer in the field of reference materials for spectroscopic analysis.

MBH CRM discs are made by casting metals that have been prepared to special formulations, usually including a range of (undesirable) trace elements. Products are checked for homogeneity, and then traceable analysis is performed by an international panel of analysts mostly operating within the terms of ISO/IEC 17025. All stages of processing, analysis and certification are in accordance with ISO guidelines for the preparation of CRMs, and each disc is supplied with a certificate showing its composition and details of the manufacturing and analytical techniques involved in its preparation.

Cat. No.	Description	Cert	tified	Value	s %																		
97593	Steel (low Alloy)	C	Si	S	Р	Mn	Ni	Cr	Мо	Cu	Sn	Al	V	W	As	Co	Nb	Ti	Pb	Sb	В	Та	Ν
		0.261	0.753	0.0476	0.043	1.19	0.886	2.09	0.788	0.130	0.047	1.136	0.296	0.358		0.470	0.322	0.260				0.26	
05602	Duplex Steel	C	Si	S	Р	Mn	Ni	Cr	Мо	Cu	AI	V	W	Co	Nb	N							
		1.492	1.89	0.0098	10.97	10.97	3.03	24.71	0.086	0.114	(0.028)	0.126			2.47	0.478			_	_	_	_	_
518/0	Copper-Nickel Alloy	C	SI	Mn	20.40	Fe	Cr	Mo	0.020	(0.06.2)	AI	Mg	5	P	Pb	Nb 0.156	Cđ	NI					
E1470	Loaded Prace	Ph	0.295 Sn	70	50.40	1.24 Ni	0.073	0.025 Si	0.029	(0.003) Bi	0.080	0.035	0.0009	0.0034	50	Cd	 R	00.38	_	-	-		
514/9	Leaded Diass	1 274	0.0912	35.34	0.299	0.282	0 384	0.022	0.098	0.0257	0.104	0.042		0.0053	0.0010	0.0030	0.0014	62.10	_	_	_	_	_
56912	Aluminium Bronze	Sn	Pb	Zn	Fe	Ni	Al	Si	As	Mn	P	Cr	Co	Ma	Aq	C.00000	Cu	02.10					
50512	Alaminani biolize	0.103	0.266	0.247	4.05	4.53	9.60	0.295	0.0066	0.056	0.048	0.0031		0.0035			80.70	_	_	_	_		
97625	Phosphor Bronze	Sn	Pb	Zn	Fe	Ni	Al	Si	As	Mn	Bi	Sb	Р	S	Со	Mg	С	Se	Cu				
		3.39	1.02	1.48	0.555	0.90	(0.001)	0.097	0.199	0.195	0.0033	0.470	0.885	0.023	0.095	0.0035			90.6				
74808	Gunmetal	Sn	Pb	Zn	Fe	Ni	Al	Si	As	Mn	Bi	Sb	Р	S	Cr	Co	Ag	Cu					
		5.22	5.13	4.19	0.136	1.008	(0.007)	(0.001)	0.059	0.0018	0.0225	0.093	0.0112	0.061		0.0662	0.0099	83.98					
66955	Copper (Impurities)	Sn	Pb	Zn	Fe	Ni	Al	Si	As	Mn	Bi	Sb	Р	S	Co	Cr	Cd	Mg	Ag	Au	Se	Те	In
		(37)	460	(24)	15	211	2		196	23	202	200	(20)	56	2		23		323	129	120	240	9
96687	Galvanizing Alloy	Pb	Mg	Al	Cd	Fe	Sn	Cu	Ni	Mn	Bi	Sb	As	Cr	Co	v							
		0.0062	0.0034	0.514	0.0006	0.0028	0.0024	0.0321	0.0441	0.0089	0.0061	0.0287	(0.0003) 0.0007	0.0037								
95653	Zinc (Impurities)	Pb	Mg	AI	Cd	Fe	Sn	Cu	Ni	Mn	Bi	Sb	Ti	TI	In	Si	Hg						
40000	7	0.0026	0.0001	0.0008	0.001/	0.00//	0.0016	0.0018	0.0013	0.0005	0.0005	0.0006		0.0012	0.0006	 R-	0.0003					_	
49892	Alloy	0.0106	0.0828	3.89	0.0053	0.0081	0.0099	1.004	0.0024	0.0086	0.0022	0.0068	0.0008	0.0063	0.011								
43998	Aluminium - Magnesium	Cu	Mg	Si	Fe	Mn	Ni	Zn	Pb	Sn	Ti	Cr	Co	V	Be	Cd	Ga	Zr	Na	В	Li	As	Мо
	Alloy	0.093	3.97	0.21	0.41	0.40	0.090	0.092	0.107	0.096	0.25	0.095			0.0074								
91927	Aluminium-Silicon-	Cu	Mg	Si	Fe	Mn	Ni	Zn	Pb	Sn	Ti	Cr	Co	V	Ве	Bi	Sb	Ag	Ca	Р	Cd	Zr	
	Copper Alloy	1.48	0.304	10.19	1.001	0.547	0.234	1.688	0.292	0.368	1.57	0.0106	0.0179	0.0179	0.0025		0.0042				0.0032		
67226	Tin (Impurities)	As	Bi	Sb	Pb	Cu	Fe	Cd	Zn	Ni	Al	Ag	s	In	Au	Te	Se	Ga	Hg				
		0.0008	(0.0025) 0.0034	0.0113	0.0013	0.0008	0.0015	0.0034	0.0004	0.0197	0.0018		0.0047	0.0007	0.0009	0.0005	0.0027	0.0044				
12242	Tin-Base Solder	As	Bi	Sb	Pb	Cu	Fe	Cd	Zn	Ni	AI	Ag	Р	Se	Ge								
	(Lead free)	0.045	0.175	1.06	0.119	2.97	(0.016)	0.0072	(0.002)	0.0203	(0.008)	0.562	0.018	0.0115									
80901	Lead (Impurities)	Sn	Sb	Bi	Cu	As	Ag	Fe	Zn	Cd	In	Ni	Te	Se	Au	TI	Na	Hg	Pt				
		0.603	0.259	1.186	0.066	0.155	0.505			0.211	0.224	0.0020	0.0011	0.0064	0.0082		0.0073	0.089					
30978	Battery Alloy with	Sn	Sb	Bi	Cu	As	Ag	Zn	Cd	Ni	AI	Te	Ca	Hg									
	Calcium	0.899	(0.002)	0.0174	0.0009	0.0007	0.0082	0.0048	0.0019	(0.0004)	0.0240	(0.0005	0.0960										
87610	Lead (high Alloy)	Sn	Sb	Bi	Cu	As	Ag	Fe	Zn	Cd	In	Ni	Te	Se	Au	TI	S	Hg					
		1.45	6.05	0.0194	0.0291	0.217	0.0071	(0.010)	(0.0007) 0.0046		0.0062	0.0071	0.0149			(0.0036)						

Table 1 Spectroscopy Products with Certified Values

Standards

Enantiomeric Purities of Pharmaceuticals Using Carbohydrate-based Isothiocyanates

Manfred P. Schneider, FB C – Bergische Universität Wuppertal, D-42097 Wuppertal, Germany schneid@uni-wuppertal.de

Introduction

In the first two papers of this series (Analytix Nr. 1/2010 and 2/2010), the efficient enantiomeric analysis of (a) alkyl oxiranes, and (b) amino acids was described using inexpensive reversed phase columns (as an alternative to high cost so-called "chiral columns") in combination with monosaccharide based iso-thiocyanates as derivatization reagents (**Figure 1**). While alkyloxiranes have to be converted first into the corresponding β-amino alcohols by reaction with isopropyl amine, all amino acids (proteinogenic, non-proteinogenic and non-natural) can be converted directly (straight out of the bottle or the body fluid, respectively) into the corresponding diastereomeric thioureas, which can be injected directly into the HPLC. Base-line separations were observed in almost all cases, establishing a highly efficient and generally applicable method for the analysis of these classes of molecules.

Pharmaceuticals

The carbohydrate-based isothiocyanates have also been shown to be highly suitable for the enantiomeric analysis of neurotransmitters (e.g. adrenaline and related molecules), numerous pharmaceuticals carrying functional amino groups, such as ß-adrenergic blockers, various pharmaceuticals such as penicillamine and mexiletine, and fine chemicals such as 1-phenyl-2-aminoethanol. Representative examples of these classes of molecules are shown in **Figure 1**.

The close relationship between the biological (physiological) activities of these molecules and their absolute configurations is well established. Frequently, only one enantiomer (*Eutomer*) has the desired pharmacological activity, while its antipode is inactive (*Distomer*), shows undesired side-effects, or is even toxic. Inactive enantiomers are also frequently referred to as *Xenobiotics* in the sense of pollutants. Several of these compounds are known for their illicit use in doping, as narcotics or psychotropic agents, and for their illegal use in food and feed. It is well established that in so-called ß-adrenergic blockers, the pharmacological activity resides in the (S)-enantiomers, while the (R)-enantiomer of penicillamine is highly toxic. On the other hand, the neurotransmitter activity of adrenaline resides largely in the (R)-enantiomer. Many more similar examples can be found in the literature. Clearly, in view of the fact that novel pharmaceuticals of this kind are increasingly used in the form of single enantiomers, the determination of enantiomeric ratios/purities is of increasing importance. Examples include the monitor of (a) enantioselective syntheses, (b) guality control in manufacturing, (c) stability and metabolic rate in biological systems (e.g. serum), and (d) analysis of illicit drugs and narcotics in body fluids (e.g. doping/body building/athletic sports). Advantageously, all of these compounds can be analyzed without any prior manipulation ("straight out of the bottle" or the reaction medium, e.g. biological fluids). They react under mild conditions and at a rapid rate (at room temperature) with the mono-saccharide isothiocyanates, leading to the corresponding diastereomeric thioureas, as exemplified in Figure 2. These, in turn, can be injected - without the need for further purification - directly into the HPLC.

As derivatives of natural mono-saccharides (**Figure 1**), all of the employed reagents are enantiomerically pure by definition, and the ratios of thus produced diastereomers directly reflect the enantiomeric composition of the chiral amino compound in question. This requires, of course, that both enantiomers of a racemic mixture react rapidly and quantitatively, and with the same rate in order to avoid a diastereoselectivity during the derivatization process. For new target molecules, this must be ascertained in every case by calibration with the corresponding racemate. The described strategy frequently has distinct advantages over



amphetamine

(S)mexiletine



ecstasy



(1R,2S) - (-)-ephedrine

and antiarrythmica





adrenaline (epinephrine)





H₂N *atenolol*

propranolol

Figure 1 Pharmacologically active molecules with functional amino groups

(R)clenbuterol

sotalol

13

14

<u>Chromatograph</u>





Figure 2 Pharmaceuticals (here a *beta-blocker*) with functionalized amino groups: formation of diastereomeric thioureas (schematic) [* denotes center of chirality]

the so-called *direct* method employing chiral stationary phases in that (a) the separation of diastereomers is usually more simple to perform and often provides better resolutions, (b) the choice of chromatographic conditions is much greater and thus can be more easily optimized, and (c) the reagents contain chromophores (fluorophores) for convenient UV- or fluorescence detection. In view of the range of novel derivatization reagents which recently became available (PGITC, PGaIITC, NGaIITC) [5], the method is an interesting alternative to so-called *chiral columns*.

In principle, all of the above reagents can be employed for the above pharmaceuticals. Thus, Nimura et al. [1] achieved base-line separations in the analysis of adrenaline (epinephrine) and noradrenaline (norepinephrine) using GITC (Aldrich T5783) and AITC (Fluka 90245). Adrenaline is only present in minute quantities in lidocain local anesthetics; nevertheless this method allowed the guantitative determination of the enantiomeric ratio in more than 250 commercially available anesthetics [2]. Using the same reagents, a series of differently substituted amphetamines were analyzed [3]. The method also works well for the analysis of ephedrine, pseudoephedrine and norephedrine (4). Further examples include the enantiomeric analysis of the antiarrhythmic agent mexiletine in human plasma using GITC [5] and the corresponding analysis of a whole series of beta-adrenergic antagonists (*ß-blockers*) such as propranolol, atenolol and sotalol (Figure 2) [6]. The introduction of benzoyl groups (BGITC Aldrich 335622) [7] and naphthoyl groups (NGalITC Fluka 04669) [8] considerably enhanced the UV- and fluorescence detection of the corresponding thioureas by factors of 6 (BIGTC) to 40 (NGalITC). Furthermore, the introduction of these residues frequently resulted in a considerable improvement in the separation of these diastereomers, as did the incorporation of extremely bulky pivaloyl groups such as in PGITC (Fluka 44891) and PGaIITC (Fluka 88102). It should also be noted that the described methodology is not limited to the referenced pharmaceuticals, but can be potentially extended to all chiral compounds carrying functional amino groups such as a wide variety of amino alcohols.

PGalITC	t ₁ (min)	t ₂ (min)	D t _{1/2} (min)
norleucin	9.29	10.81	1.52
sotalol	20.68	22.10	1.42
atenolol	20.49	21.90	1.41
penicillamine	10.93	12.74	1.81
DOPA	20.32	24.41	4.09
2-amino-1-phenylethanol	21.01	22.31	1.30

 Table 1 Mobile phase: acetonitrile: 0.1% TFA/H₂0 = 70:30

NGalITC	t ₁ (min)	t ₂ (min)	D t _{1/2} (min)
atenolol	24.53	27.01	2.48
sotalol	23.20	25.28	2.08
penicillamine	25.50	28.41	2.91
clenbuterol	23.70	27.45	3.75

Table 2 Mobile phase: acetonitrile : 0.1% TFA/H2O = 60:40

BGITC	k _s	k _R	а
propranolol	2.65	3.65	1.38
pindolol	16.88	18.95	1.12
atenolol	12.89	16.81	1.30
sotalol	11.94	13.90	1.16
penicillamine	7.37	10.05	1.33

Table 3 Mobile phase: MeOH: H₂O = 80:20 to 90:10

While simple RP-18 columns are generally employed, the separation conditions can be varied widely in order to achieve the best separating conditions. Various different mobile phases have been used ranging from MeOH : phosphate buffer (pH 2.8) [1] over MeOH : H_2O : phosphate buffer (pH 7) to acetonitrile : water : 0.1% trifluoroacetic acid in order to optimize the separation conditions. In certain cases the reagent may interfere with the separation, having the same or similar retention time. The addition of small amounts of ethanolamine or hydrazine is sufficient to destroy excess reagent by formation of the corresponding thioureas, which elute at different retention times.

Summary

The method described above allows the rapid, efficient and inexpensive determination of enantiomeric purities in a wide variety of structurally varied pharmaceuticals and fine chemicals. By using the suitable derivatization reagent, base-line separations are observed in nearly all cases. The procedure is quite general and applicable to (a) detecting enantiomeric ratios of pharmaceuticals, in addition to biological samples; (b) determining racemizations and differences in metabolic degradation; (c) monitoring asymmetric syntheses; and (d) detecting molecules in illicit drug abuse and doping. The method is clearly adaptable to automation using reaction batteries and auto-samplers. The technique is applicable both on a laboratory scale and in on-line quality control. It is thus highly suitable for monitoring asymmetric syntheses including enzyme-catalyzed transformations.

Chromatography

Experimental

5 mg of the corresponding pharmaceutical is dissolved in 50% (v/v) aqueous acetonitrile (or dimethylformamide) containing 0.55% (v/v) triethyl amine (in this case, hydro chlorides are employed) to give a final volume of 10 mL. To 50 μ L of this stock solution 50 μ L of 066% (w/v) BGITC in acetonitrile is added. The resulting solution is shaken on a laboratory shaker for 30 min, after which 10 µL of 0.26% (v/v) ethanolamine (or hydrazine) in acetonitrile is added and shaking is continued for another 10 min. Ethanolamine (hydrazine) reacts with any excess of BGITC and the resulting thiourea derivative is eluted well behind any of the amino acid derivatives. The mixture is then diluted to a final volume of 1 mL and a 10 μ L aliguot is injected into the HPLC. (RP-18, mobile phase MeOH : H₂O [67 mM phosphate buffer (pH 7) = 65:27:8 up to 70:25:5 and 80:15:5], depending on the case, flow rate 0.5 mL/min, compare Tables).

References

 Nimura, N., Kasahara, Y., Kinoshita, T., resolution of enantiomers of norepinephrine and epinephrine by reversedphase high-performance liquidchromatography *J. Chromatogr.* 213 (1981) 327–330.

- [2] Allgire, J.F., Juenge, E.C., Adamo, C.P., Sullivan, G.M., Kirchhoefer, R.D. High-performance liquid chromatographic determination of d-/l-epinephrine enantiomer ratio in lidocain-epinephrine local anesthetics; *J. Chromatogr.* 325 (1985) 249–254.
- [3] Miller, K.J., Gal, J., Ames, M.M. High-performance liquid chromatographic resolution of enantiomers of 1-Phenyl-2-1minopropanes (amphetamines) with four chiral reagents, *J. Chromatogr.* 307 (1984) 335–342.
- [4] Gal, J. Resolution of the enantiomers of ephedrine, norephedrine and pseudoephedrine by high-performance liquid chromatography, J. Chromatogr. 307 (1984) 220–233.
- [5] Grech-Belanger, O., Turgeon, J. High-performance liquid chromatographic assay for mexiletine enantiomers in human plasma, J. Chromatogr. 337 (1985) 172-177.
- [6] Sedman, A.J., Gal, J. resolution of the enantiomers of propranolol and other beta-adrenergic antagonists by highperformance liquid chromatography, J. Chromatogr. Biomedical Appl. 278(1983) 199–203.
- [7] Lobell, M., Schneider, M., 2,3,4,6-Tetra-O-benzoyl-ß-Dglucopyranosyl isothiocyanante: an efficient reagent for the determination of enantiomeric purities of amino acids, ß-adrenergic blockers and alkyloxiranes by high-performance liquid chromatography using standard reversed phase columns, J. Chromatogr. 633 (1993) 287–294.
- [8] Schneider, M., recent results, unpublished.

Cat. No.	Brand	Description	Abbr.	Package Size
90245	Fluka®	2,3,4-Tri-O-acetyl-α-D-arabinopyranosyl-isothiocyanate	AITC	100 mg, 500 mg
T5783	Sigma	2,3,4,6-Tetra-O-acetyl-ß-D-glucopyranosyl-isothiocyanate	GITC	100 mg, 1 g
335622	Aldrich	2,3,4,6-Tetra-O-benzoyl-B-D-glucopyranosyl-isothiocyanate	BGITC	500 mg
44891	Fluka	2,3,4,6-Tetra-O-pivaloyl-ß-D-glucopyranosyl-isothiocyanate	PGITC	100 mg
88102	Fluka	2,3,4,6-Tetra-O-pivaloyl-ß-D-galactopyranosyl-isothiocyanate	PGalITC	100 mg, 500 mg
04669	Fluka	2,3,4,6-Tetra-O- (2-naphthoyl)-ß-D-galactopyranosyl-isothiocyanate	NGalITC	25 mg, 100 mg

Product Table for Carbohydrate-based Isothiocyanates



Ascentis[®] Express Fused-Core[™] Technology Redefining the Limits of Your HPLC

Based on Fused-Core particle technology, Supelco's Ascentis Express columns provide a breakthrough in HPLC column performance.

- Double the efficiencies of conventional 3 µm particles
- Equal efficiencies of sub-2 µm columns at half of the backpressure
- Rugged design capable of high pressure operation

For additional information, visit us on the web: *sigma-aldrich.com/express* Ascentis is a registered trademark of Sigma-Aldrich Biotechnology LP. Fused-Core is a trademark of Advanced Materials Technologies, Inc. Now available ...

Ascentis Express Peptide ES C18 for peptides and small proteins.



Highest Purity Acids for Ultra-Trace Analysis

Highest quality for the best price

TraceSELECT®Ultra

Ultra-pure acids, bases, and salts for smelting and wet digestion in environmental, water and food analyses

Sample preparation for trace analysis requires reagents of the highest purity. Our *TraceSELECT Ultra* acids for ultra-trace analysis **at ppb and even ppt levels** are produced by sub-boiling distillation. Sub-boiling is recognized as the best way to obtain high purity acids with the lowest blank values for ultra-trace analysis. The technique is based on the evaporation of liquid by infrared heating at the surface. It avoids violent boiling and the formation of liquid aerosols that can be transported with the distillate.

To maintain their high purity, *TraceSELECT* Ultra products are supplied in PFA (fluoropolymer) bottles. Water, ammonium hydroxide solution, TMAOH and ortho-phosphoric acid are supplied in specially pre-leached HDPE bottles. Recent process improvements have allowed us to reduce our impurity specifications in order to guarantee the lowest levels of trace impurities in our *TraceSELECT* Ultra products.



For more information, see our website: sigma-aldrich.com/traceselect

Cat. No.	Brand	Product Name	Concentration	Pack sizes
07692	Fluka®	Acetic acid	≥ 99.0%	250 mL/1 L
16748	Fluka	Ammonium hydroxide solution	≥ 25%	250 mL/1 L
23828	Fluka	Hydrobromic acid	≥ 44%	250 mL/1 L
96208	Fluka	Hydrochloric acid	≥ 30%	250 mL/1 L
02658	Fluka	Hydrofluoric acid	≥ 49%	250 mL/1 L
16911	Fluka	Hydrogen peroxide solution	≥ 30%	250 mL/1 L/5 L
02650	Fluka	Nitric acid	~ 65%	250 mL/1 L
64957	Fluka	Phosphoric acid	≥ 85%	250 mL/1 L
77239	Fluka	Sulfuric acid	≥ 95%	250 mL/1 L
14213	Fluka	Tetramethylammonium hydroxide solution (TMAOH)	$\sim 25\%$ in H_2O	250 mL
14211	Fluka	Water		1 L

Chromatography

SIGMA-ALDRICH®

Chromatography





SUPELCO

Otv/Pk

20

20

15

12

Otv/Pk

16

16

Silica Type

irregular

irregular

irregular

irregula

Supel Flash Cartridges for Isco and Analogix Flash Systems

Silica Type

irregular

irregular

irregular

irregular

irregular

irregular

irregular

irregular

spherical

spherical

spherical

spherical

Supel Flash FM Cartridges for Biotage FlashMaster Flash Systems

Tube Volume

40–63 µm

40–63 µm

Cat. No. Bed Wt.

FCISI004 4 g 40-60 µm

FCISI012 12 g 40-60 µm

FCISI025 25 g 40-60 µm

FCISI040 40 g 40-60 µm

FCISI080 80 g 40-60 µm

FCISI240 240 g 40-60 μm

97790-U 100 g 20-45 μm

Cat. No. Bed Wt.

120 g 40 – 60 µm

330 g 40–60 µm

11 a 20 – 45 um

23 g 20–45 µm

50 g 20–45 µm

10 g 70 mL

20 g 70 mL

97717-U 50 a 150 ml 40-63 um

97718-U 70 g 150 mL 40-63 μm

FCISI120

FCISI330

97787-U

97788-U

97789-U

97715-U

97716-U

Supel[™] Flash Cartridges

Supel Flash cartridges utilize proprietary column packing technology and are packed with high-purity silica gel. The column packing process results in a high-density, uniform bed insuring the elimination of channeling and an overall efficiency improvement. Supel Flash cartridges are available in two styles (Supel Flash and Supel Flash FM) for compatibility with different flash chromatography systems. The standard Supel Flash cartridge has a female luer lock inlet fitting and a male slip luer outlet fitting making it compatible with Isco, Analogix and some Biotage® flash systems. The Supel Flash FM cartridges are SPE-tube style cartridges designed for use in Biotage Flashmaster[™] systems.

Features include:

- High efficiency leading to sharper and more symmetric peaks
- No fronting, doubling or overlapping of peaks
- High yield and sample loading capacity
- No compound decomposition
- Reliable, consistent purification performance
- Single-piece cartridge design with Luer-Lok fittings

Sigma-Aldrich has a wide selection of bulk silica products and low pressure liquid chromatography columns and accessories. For more information visit **sigma-aldrich.com/versaflash**

TRADEMARKS: Biotage, Flashmaster – Biotage AB; CombiFlash, Companion, Retrieve, OptiX – Teledyne Isco; IntelliFlash, SimpliFlash – Varian, Inc.; Supel – Sigma-Aldrich Biotechnology LT



If you need more information please tick the box for our eCommerce brochure or get directly in contact with us: sigma-aldrich.com/ebusiness

SIGMA-ALDRICH

sigma-aldrich.com/register

Analysis of Bioethanol – Detection of Chloride and Sulfate with Dication and Trication Solutions

R. Köhling, N. Reichlin, Sigma-Aldrich Production GmbH (CH), R&D – LC/MS Applications rudolf.koehling@sial.com



Biodiesel and bioethanol are increasingly being used as renewable energy sources to replace the use of fossil fuels for common combustion engines. Diesel engines can be converted for use with plant oils or other fats, but significant changes must be made to the vehicle. Modern biodiesel is an alternative to plant oils and used increasingly in diesel vehicles, e.g. trucks. There is also a renewable fuel source for gasoline engines which can run on ethanol-gasoline blends instead of pure gasoline. This sourcing of fuels from biomass is a common technique in countries with a large corn, sugar cane or general biomass emergence (USA, Brazil). The bioethanol can be obtained in large amounts through fermentation and distillation and is added to gasoline in varying amounts, e.g. 85% ethanol in E85 fuel [1]. Bioethanol and biodiesel contain several impurities despite passing through several cleaning steps. Especially concerning are dissolved salts which can damage modern engines. Thus, the determination of sulfate and chloride in ethanol-based fuels is an important quality criterion. A standard detection method for these analytes is ion chromatography configured with a conductivity detector [2]. It is a very sensitive analytical technique, but it lacks the capability to definitively identify the compounds in addition to retention time. An alternative to IC can be Di- (75128) and tricationic organic compounds (08675), which form positively charged adducts with chloride and sulfate anions, making them detectable for mass spectroscopy in the highly sensitive positive ESI mode [3]. The use of LC/MS also has the advantage of easy sample handling. For a large number of matrices, the samples can be injected without further treatments.

Independent of analytical techniques, the quality of quantitative results strongly depends on the precision and accuracy of the analytical reference standards. A robust calibration method still can lead to results with a high uncertainty, when the content of the reference material is not well defined. In this case, matrix samples are spiked with *TraceCERT®* anion standards for ion chromatography in order to ensure the highest accuracy for the calibration data.



Figure 1 Schematic drawing of the LC/MS set-up and the data processing of the MS data

Method

There are two favored methods for the application of multiply charged cation solutions with LC/MS: (1) a direct injection of a mixture of cation solution, sample and solvent (e.g. methanol/water) into the MS, or (2) the injection of the sample into the constant flow of a mobile phase, which is mixed with a cation solution post-column via a T-connector (**Figure 1**).

The second method is preferred for the analysis of biofuels since it provides the highest sensitivity and opens the possibility to change all compositions during the method development. Additionally, an analytical column can remove parts of the matrix to prevent suppression effects. However, one should be aware that ionic analytes could interact with metal parts throughout the entire HPLC system, columns or tubing. If this happens and peak broadening occurs, then PEEK tubing, injectors and columns can prevent this negative influence on the peak height and shape.

Starting with the HPLC, the pump delivers a typical mobile phase water and methanol (90/10, v/v) at a low flow rate of 0.2 ml/min. Additives like formic acid, TFA or acetic acid should not be used, as they will form adducts with the multi-cationic reagent and will lower its efficiency. The solvents should have the best available quality to minimize bias and noise. Reducing the flow rate also reduces suppression effects caused by sample matrices [4]. The di- and tricationic fluoride solution is added to the mobile phase with a flow rate of 80 μ l/h, but can be increased in case of low signal intensity of the adduct (ion pair of di-/trication and anion).

The calibration samples are prepared according to DIN 38402 part 51 and 32645 with equidistant concentrations. 5 calibration levels and 1 blank sample cover a concentration range from 2 to 10 μ g/ml and result in a typical limit of detection (LOD) of 0.5 μ g/ml and a lowest limit of quantification (LLQ) of 1.2 μ g/ml (**Figure 2**). Ethanol and water serve as matrix for the preparation of the calibration standards. Finally a sample of the lab water supply is analyzed without any sample pre-treatments.

Discussion

Solutions of di- and tricationic compounds can easily be used for the determination of chloride and sulfate in ethanol or drinking water. The lab water supply was directly injected into the HPLC system and resulted in a typical chloride concentration of 4 ppm. The solutions of the



Figure 2 Typical calibration curve for sulfate ions in water detected as trication-sulfate adduct. Similar results are obtained for chloride ions in ethanol detected as a dication-chloride adduct. Correlation coefficients of 0.996 to 0.999 are obtained for the linear fit of the calibration data. The calculation of the confidence and prediction bands is based on a confidence level of 0.95.

di-/trications are supplied in a condition capable of being used directly in mixing with the mobile phase post-column. This set-up is installed on most LC/MS systems and does not need additional items except a syringe pump.

Mass spectroscopy with electrospray ionization (ESI) has the advantage of a very sensitive and selective detection, which is responsible for low LODs and LLQs as can be seen in **Figure 2**. Although 2 ppm is the lowest calibration level, it is possible to get lower LLQs even at a confidence level of 0.99 because of the steep regression curve. One important advantage of MS is its capability to definitively identify analytes by their mass (exact mass), isotopic pattern and MS/MS spectra.

In addition to IC, this LC/MS method can also be a useful tool for the analysis of anions in complex matrices like biofuels, especially for those laboratories equipped with an LC/MS system.

References

[1] www.bio-kraftstoffe.info

- [2] Norm ASTM D4806 "Standard Specification for Denatured Fuel Ethanol for Blending with Gasolines for Use as Automotive Spark-Ignition Engine Fuel".
- [3] R. Köhling, N. Reichlin, "Highly Sensitive Detection of Organic and Inorganic Anions with Di- or Tricationic LC/MS Additives", Analytix, 2, 2009.
- [4] F. Gosetti, E. Mazzucco, D. Zampieri, M. Gennaro, "Comparison of Matrix Effects in HPLC-MS/MS and UPLC-MS/MS Analysis of Nine Basic Pharmaceuticals in Surface Waters", J. Am. Soc. Mass. Spectrom., 2008, 19, 713–718.

Cat. No.	Brand	Description	Package Size
8675	Fluka®	1,3-Bis[6-(3-benzyl-1-imidazolio)-hexyl]imidazolium trifluoride solution	100 mL
75128	Fluka	1,9-Nonanediyl-bis(3-methylimidazolium) difluoride solution	100 mL
79254	Fluka	1,3,5-Tris[(3-butyl-imidazolium)methyl]-2,4,6-trimethylbenzene trifluoride solution	100 mL
78897	Fluka	1,3,5-Tris[(tripropylphosphonium)methyl]benzene trifluoride solution	100 mL
56618	Fluka	1,5-Pentanediyl-bis(1-butylpyrrolidinium) difluoride solution in water/methanol 1:1	100 mL
76507	Fluka	1,5-Pentanediyl-bis(3-benzylimidazolium) difluorid solution in water/methanol 1:1	100 mL

Increasing Reliability and Value of Analytical Results

Sigma-Aldrich achieves accreditation according to ISO/IEC 17025 for Water Standards and Karl Fischer Reagents

HYDRANAL® Technical Service accredited according to ISO/IEC 17025

Sigma-Aldrich's HYDRANAL® service accreditat laboratory in Seelze, Germany, defined ar

proved its expertise for testing and

calibration and achieved accreditation as a testing laboratory according to ISO/IEC 17025:2005 "General requirements for the competence of testing and calibration laboratories". It was audited by the German Accreditation Council (Deutscher Akkreditierungsrat DAR) and received accreditation for titration and gravimetry for water standards and reagents for Karl Fischer titration.

From the beginning of our HYDRANAL product line over 30 years ago, the HYDRANAL technical service in Seelze has applied a high standard and continual improvement of testing methods for production and quality control of HYDRANAL water standards and HYDRANAL Karl Fischer reagents. These high-quality testing methods have now been proven to comply with the high standards of the auditors and have led to the accreditation of HYDRANAL Technical Service according to ISO/IEC 17025. This accreditation details, along with the logo of the accreditation body, will now appear on the Certificate of Analysis (see **Figure 1**) and on the label for each produced lot.

Accreditation requirements and procedure

Accreditation means that an authoritative body formally recognizes that an organization or individual is competent to execute a specific service as described in the scope of

Helga Hoffmann, Technical Service HYDRANAL Manager helga.hoffmann@sial.com Andrea Felgner, Product Manager Analytical Reagents andrea.felgner@sial.com

accreditation. ISO/IEC 17025 is associated with a welldefined analytical technique and a stated measurement range comprising analyte, matrix and concentration range. This is called the scope of the accreditation. In this case, HYDRANAL Technical Service is accredited to perform titration and gravimetry for water standards and reagents for Karl Fischer titration according to the following methods/ norms: ISO 760:1978-12, ASTM E203-08 and HYDRANAL Technical Service's in-house methods 01-04.

To fulfill the requirements of ISO/IEC 17025, a laboratory must demonstrate its general competence and compliance in performing specific tests or calibrations, as well as establishing and documenting its quality management system. Beyond the basic management requirements of ISO 9001, ISO/IEC 17025 requires the following topics:

- Instrument qualification
- Validation of analytical methods
- Traceability statements
- Evaluation of measurement uncertainty
- Ongoing education of personnel
- Periodic participation in proficiency tests to demonstrate technical capability

Following the accreditation, evidence of continued improvement is required, including regular internal audits and the demonstration of efforts towards scientific and technological progress.

Certificate of Analysis sigma-aldrich.com
Analyzed for Sigma-Aldrich Laborchemikalien GmbH, Production HYDRANAL [®] , Wunstorfer Str. 40, D-30926 Seelze
HYDRANAL®-Standard Sodium tartrate dihydrate, reference no. 34696, is intended for the titre
determination of volumetric Karl Fischer reagents. The certified water content of this Lot
(SZE8347B, Exp.date Nov. 2013) is
15.63% (expanded uncertainty = 0,02%, k=2.78).
The water content is analyzed by loss on drying at 150 °C on 5 samples.
HYDRANAL [®] -Standard Sodium tartrate dihydrate can be used for standardization of KF reagents according to ISO 760 or ASTM E 203. This standard is traceable to SI-Unit (kg) and tested against NIST SRM 2890.
Sigma-Aldrich Laborchemikalien GmbH Helga Hoffmann Manager Technical Service
HYDRANAL® DGAPL6870.09 HELGA IF A PLANT page 1 of 1

Figure 1 Example of HYDRANAL Certificate of Analysis with accreditation details

Wunstorfer Str.40, D-30926 Seelze



Analytical results produced in laboratories fulfilling ISO/IEC 17025 can be labeled with a special quality logo that confirms the reliability of the results. Every country has at least one official body that is responsible for accreditation, and each body signifies its accreditation with its own unique logo. Since HYDRANAL® Technical Service was accredited by DAR, the DAR symbol and the individual registration number of the laboratory must appear on the certificate of Karl Fischer reagents and water standards. Beginning in 2010, the authoritative body for accreditations in Germany is the Deutsche Akkreditierungsstelle (DAkkS), which will also administrate accreditations issued by former authorities.

Helga Hoffmann **Technical Service HYDRANAL** Wunstorfer Straße 40 D-30926 Seelze, Germany Tel. +49 (0) 5137 8238-353 Fax +49 (0) 5137 8238-698 E-Mail: helga.hoffmann@sial.com



Figure 2 Official accreditation document for HYDRANAL Technical Service

Cat. No.	Product	Description	Pack Size
34849	HYDRANAL-Water Standard 10.0	Standard for volumetric KF titration	80 mL
		1 g (= 1 mL at 20 °C) contains 10.0 mg = 1.00% water	
		Contains 10 glass ampoules of 8 mL. Tested against NIST SRM 2890	
34828	HYDRANAL-Water Standard 1.00	Standard for coulometric KF titration	40 mL
		1 g (= 1 mL at 20 °C) contains 1.00 mg = 0.10% water	
		Contains 10 glass ampoules of 4 mL. Tested against NIST SRM 2890	
34847	HYDRANAL-Water Standard 0.10	Standard for coulometric KF titration	40 mL
		1 g contains 0.10 mg = 0.01% water	
		Contains 10 glass ampoules of 4 mL. Tested against NIST SRM 2890	
34694	HYDRANAL-Water Standard Oil	Mineral oil-based standard with water content in low ppm-range (CoA)	80 mL
		Contains 10 glass ampoules of 8 mL	
34693	HYDRANAL-Water Standard KF-Oven	Solid standard for control of KF Ovens, Water content ~5%	10 g
	140–160 °C		
34748	HYDRANAL-Water Standard KF-Oven	Solid standard for control of KF Ovens, Water content: 5.55 \pm 0.05%	10 g
	220−230 °C		

Table 1 HYDRANAL Water Standards (ISO/IEC 17025 accreditation is valid for all Water standards and Karl Fischer reagents)

Monthly Savings Program **SAVE 35%**

16

FIXAMAL

Fluka

Fluka

FIXAMAL Gruns

Volumetric Titration Reagents by Sigma-Aldrich -The variety you need for laboratory success

The following FIXANAL® concentrates for volumetric solutions are available with a HUGE SAVING OF 35% OFF until the end of July:

Carris & TUNU

Cat. No.	Brand	Description	
38272	Fluka®	Hydrochloric acid std., concentrate, pkg of 0.01 mol (0.3646 g HCl)	
38280	Fluka	Hydrochloric acid std., concentrate, pkg of 0.1 mol (3.646 g HCl)	
38287	Fluka	Hydrochloric acid std., concentrate, pkg of 0.2 mol (7.292 g HCl)	
38285	Fluka	Hydrochloric acid std., concentrate, pkg of 0.5 mol (18.231 g HCl)	
38282	Fluka	Hydrochloric acid std., concentrate, pkg of 1.0 mol (36.461 g HCl)	
38281	Fluka	Hydrochloric acid std., concentrate, pkg of 2.0 mol (72.922 g HCl)	
38283	Fluka	Hydrochloric acid std., concentrate, pkg of 10 mol (364.61 g HCl)	
38070	Fluka	Potassium hydroxide std., concentrate, pkg of 0.1 mol (5.611 g KOH)	
38073	Fluka	Potassium hydroxide std., concentrate, pkg of 1.0 mol (56.109 g KOH)	
38227	Fluka	Sodium hydroxide std., concentrate, pkg of 0.01 mol (0.400 g NaOH)	
38226	Fluka	Sodium hydroxide std., concentrate, pkg of 0.025 mol (1.000 g NaOH)	
38210	Fluka	Sodium hydroxide std., concentrate, pkg of 0.1 mol (4.000 g NaOH)	
32041	Fluka	Sodium hydroxide std., concentrate, pkg of 0.2 mol (7.999 g NaOH)	
32040	Fluka	Sodium hydroxide std., concentrate, pkg of 0.25 mol (9.999 g NaOH)	
38217	Fluka	Sodium hydroxide std., concentrate, pkg of 0.5 mol (19.999 g NaOH)	
38215	Fluka	Sodium hydroxide std., concentrate, pkg of 1.0 mol (39.997 g NaOH)	
38212	Fluka	Sodium hydroxide std., concentrate, pkg of 2.0 mol (79.994 g NaOH)	
38214	Fluka	Sodium hydroxide std., concentrate, pkg of 10 mol (399.97 g NaOH)	
38243	Fluka	Sodium thiosulfate std., concentrate, pkg of 0.01 mol (2.482 g $Na_2S_2O_3 \cdot 5H_2O$)	
38200	Fluka	Sodium thiosulfate std., concentrate, pkg of 0.1 mol (24.818 g $Na_2S_2O_3$ ·5H $_2O$)	
38308	Fluka	Sulfuric acid standard std., concentrate, pkg of 0.005 mol (0.4904 g $\rm H_2SO_4$)	
32043	Fluka	Sulfuric acid standard std., concentrate, pkg of 0.05 mol (4.904 g $\mathrm{H_2SO_4}$)	
38295	Fluka	Sulfuric acid standard std., concentrate, pkg of 0.25 mol (24.519 g $\rm H_2SO_4$)	
38294	Fluka	Sulfuric acid standard std., concentrate, pkg of 0.5 mol (49.039 g $\rm H_2SO_4$)	
38291	Fluka	Sulfuric acid standard std., concentrate, pkg of 1.0 mol (98.078 g $\rm H_2SO_4$)	
32044	Fluka	Sulfuric acid standard std., concentrate, pkg of 5 mol (490.39 g $\rm H_2SO_4$)	

To take advantage of this offer, please use promotion code 982. Offer valid until July 31, 2010.

SIGMA-ALDRICH®

Matthias Nold, Product Manager Analytical Standards matthias.nold@sial.com



Alkylresorcinols are naturally occurring phenolic lipids that are present in substantial amounts in the outer layer of cereal grains. They are especially abundant in the bran layer of wheat and rye. Under the Fluka® brand, Sigma-Aldrich now launches a series of analytical standards of alkylresorcinols with different chain lengths, as well as the metabolite 3,5-Dihydroxycinnamic acid (3-(3,5-dihydroxyphenyl)-1-propanoic acid (DHPPA)).

Alkylresorcinols are known to have antimicrobial activity and are believed to have an antitumor function as well. Since they are not present in the endosperm, they are virtually absent in white bread. Therefore, alkylresorcinols are often used as biomarkers for whole-grain intake in research investigating the health benefits of whole grain consumption.



Alkylresorcinols (General Formula) 3,5-Dihydroxyhydrocinnamic

acid

Cat. No.	Brand	English desc	Pack Size	n
56453	Fluka	5-Tridecylresorcinol	10 mg	11
91822	Fluka	5-Pentadecylresorcinol	10 mg	13
97001	Fluka	5-Heptadecylresorcinol	10 mg	15
57981	Fluka	5-Nonadecylresorcinol	10 mg	17
49519	Fluka	5-(Nonadecyl-1,1,2,2-d4)resorcinol	10 mg	17
53503	Fluka	5-Eicosylresorcinol	10 mg	18
50851	Fluka	5-Heneicosylresorcinol	10 mg	19
03422	Fluka	5-Tricosylresorcinol	10 mg	21
56452	Fluka	3,5-Dihydroxyhydrocinnamic acid	10 mg	-

New Standards for Regulated Food Dyes

Matthias Nold, Product Manager Analytical Standards matthias.nold@sial.com



Sigma-Aldrich offers a wide range of reference standards for the accurate detection of regulated dyes. These dyes, although banned worldwide for use as food additives, are still being used illicitly as additives in food products. We recently expanded our portfolio to include reference standards for the detection of additional food dyes. Please visit our webpage at **sigma-aldrich.com/fooddyes** for a complete product listing.

Cat. No.	Description	Brand	Package Size
11298	Carminic Acid	Fluka	25 mg
89774	Citrus Red 2	Fluka	25 mg
68775	Sunset Yellow	Fluka	25 mg
59659	Toluidine Red	Fluka	25 mg



Sigma-Aldrich Worldwide Offices

Argentina

Free Tel: 0810 888 7446 Tel: (+54) 11 4556 1472 Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097 Free Fax: 1800 800 096 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10 Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747 Free Fax: 0800 14745 Tel: (+32) 3 899 13 01 Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425 Tel: (+55) 11 3732 3100 Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400 Free Fax: 1800 265 3858 Tel: (+1) 905 829 9500 Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395 Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336 Tel: (+86) 21 6141 5566 Fax: (+86) 21 6141 5567

Czech Republic Tel: (+420) 246 003 200 Fax: (+420) 246 003 291

Denmark Tel: (+45) 43 56 59 00 Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250 Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408 Free Fax: 0800 031 052 Tel: (+33) 474 82 28 88 Fax: (+33) 474 95 68 08

Germany Free Tel: 0800 51 55 000 Free Fax: 0800 64 90 000 Tel: (+49) 89 6513 0

Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355 Ingyenes fax szám: 06 80 344 344 Tel: (+36) 1 235 9063 Fax: (+36) 1 269 6470

India

Telephone Bangalore: (+91) 80 6621 9400 New Delhi: (+91) 11 4358 8000 Mumbai: (+91) 22 2570 2364 Hyderabad: (+91) 40 4015 5488 Kolkata: (+91) 33 4013 8003

Fax

Bangalore: (+91) 80 6621 9650 New Delhi: (+91) 11 4358 8001 Mumbai: (+91) 22 4087 2364 Hyderabad: (+91) 40 4015 5488 Kolkata: (+91) 33 4013 8000

Ireland

Free Tel: 1800 200 888 Free Fax: 1800 600 222 Tel: (+353) 402 20370 Fax: (+ 353) 402 20375

Israel

Free Tel: 1 800 70 2222 Tel: (+972) 8 948 4100 Fax: (+972) 8 948 4200

Italy

Tel: (+39) 02 3341 7310 Fax: (+39) 02 3801 0737

Japan Tel: (+81) 3 5796 7300 Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111 Free Fax: (+82) 80 023 8111 Tel: (+82) 31 329 9000 Fax: (+82) 31 329 9090

Malaysia Tel: (+60) 3 5635 3321 Fax: (+60) 3 5635 4116

Mexico Free Tel: 01 800 007 5300 Free Fax: 01 800 712 9920 Tel: (+52) 722 276 1600 Fax: (+52) 722 276 1601

The Netherlands Free Tel: 0800 022 9088 Free Fax: 0800 022 9089 Tel: (+31) 78 620 5411 Fax: (+31) 78 620 5421

New Zealand Free Tel: 0800 936 666 Free Fax: 0800 937 777 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Norway Tel: (+47) 23 17 60 00 Fax: (+47) 23 17 60 10

Poland Tel: (+48) 61 829 01 00 Fax: (+48) 61 829 01 20

Portugal Free Tel: 800 202 180 Free Fax: 800 202 178 Tel: (+351) 21 924 2555 Fax: (+351) 21 924 2610

Russia Tel: (+7) 495 621 5828 Fax: (+7) 495 621 5923

Singapore Tel: (+65) 6779 1200 Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562 Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75 Free Fax: 0800 1100 79 Tel: (+27) 11 979 1188 Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376 Free Fax: 900 102 028 Tel: (+34) 91 661 99 77 Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200 Fax: (+46) 8 742 4243

Switzerland Free Tel: 0800 80 00 80 Free Fax: 0800 80 00 81 Tel: (+41) 81 755 2828 Fax: (+41) 81 755 2815

United Kingdom Free Tel: 0800 717 181 Free Fax: 0800 378 785

Tel: (+44) 1747 833 000 Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010 Toll-Free Fax: 800 325 5052 Tel: (+1) 314 771 5765 Fax: (+1) 314 771 5757

Vietnam Tel: (+84) 3516 2810 Fax: (+84) 6258 4238

Internet sigma-aldrich.com



Mixed Sources Product group from well-managed forests, controlled sources and recycled wood or fiber FSC www.fsc.org Cert no. GFA-COC-001446 © 1996 Forest Stewardship Council

Accelerating Customers' Success through Innovation and Leadership in Life Science, High Technology and Service

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052 Technical Service (800) 325-5832 • sigma-aldrich.com/techservice Development/Custom Manufacturing Inquiries SAFC° (800) 244-1173 Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters 3050 Spruce St. St. Louis, MO 63103 (314) 771-5765 sigma-aldrich.com

Date: 06/2010; SAMS Code: MIO

©2010 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, LP. Sigr brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip. Fluka, TraceSELECT, TraceSELECT ULTRA, PERDROGEN, CHROMASOLV, FIXANAL, HYDRANAL, Riedel-de Haën, IDRANAL, SPECTRANAL, and VOLPAC are trademarks of Honeywell Specialty Chemicals Seelze GmbH Sigma-Aldrich Chemie GmbH is a subsidiary of Merck

SIGMA-ALDRICH[®]