



Product Information

SIGMA QUALITY CONTROL TEST PROCEDURE

Enzymatic Assay of LYSOSTAPHIN

PRINCIPLE:

Staphylococcal Cells (intact) $\xrightarrow{\text{Lysostaphin}}$ Staphylococcal Cells (lysed)

CONDITIONS: T = 37°C, pH 7.5, $A_{620\text{nm}}$, Light path = 1 cm

METHOD: Turbidimetric Rate Determination

REAGENTS:

- A. 50 mM Tris HCl Buffer with 145 mM Sodium Chloride, pH 7.5 at 37°C
(Prepare 1 liter in deionized water using Trizma Base, Sigma Prod. No. T-1503 and Sodium Chloride, Sigma Prod. No. S-9625. Adjust to pH 7.5 at 37°C with 1 M HCl.)
- B. Cell Suspension¹
(Prepare using the following steps:
 1. Prepare 100 ml of a 30% (w/v) Tryptic Soy Broth in Reagent A using Tryptic Soy Broth, Sigma Prod. No. T-8261. In a 250 ml flask, inoculate 50 ml of the Tryptic Soy Broth with 1 ml of a 1 mg/ml culture of *Staphylococcus aureus*,² American Type Culture Collection #29737. Incubate, with gentle shaking, at 37°C for 18 - 24 hours.
 2. Centrifuge the culture to pellet the cells. Remove the supernatant liquid and suspend the cells in cold Reagent A. Repeat the centrifugation and resuspend the cells in a minimum volume of cold Reagent A.
 3. Adjust the cell density to an $A_{620\text{nm}}$ of 1.3 - 1.5 by diluting with cold Reagent A. Keep the diluted suspension on ice until ready to assay samples.

THE FINAL DILUTED SUSPENSION (STEP 3) MUST BE PREPARED FRESH. USE THE STEP 3 CELL SUSPENSION WITHIN 15 - 20 MINUTES.)

- C. Lysostaphin Enzyme Solution
(Immediately before use, prepare a solution containing 0.05 - 0.25 unit/ml of Lysostaphin in cold Reagent A.)

Enzymatic Assay of LYSOSTAPHIN

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable vials:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	-----	5.00
Reagent C (Enzyme Solution)	5.00	-----

Equilibrate to 37°C for approximately 5 minutes. Then add:

Reagent B (Cell Suspension)	1.00	1.00
-----------------------------	------	------

Mix by swirling and immediately determine at zero time, (T_0), the ΔA_{620nm}^3 , for both the Test and the Blank using Reagent A as a reference.

Incubate the cuvette at 37°C for exactly 10 minutes. Mix by swirling and determine at ten minutes, (T_{10}), the A_{620nm} for both the Test and the Blank.

CALCULATIONS:

$$\Delta A_{620nm} = A_{620nm} \text{ Test} - A_{620nm} \text{ Blank}$$

$$\% \text{ Lysis} = \frac{\Delta A_{620nm} T_0 - \Delta A_{620nm} T_{10}}{\Delta A_{620nm} T_0} \times 100\%$$

$$\text{Units/ml} = \frac{\Delta A_{620nm} T_0 - \Delta A_{620nm} T_{10}}{(0.125)(5)} \text{ (df)}$$

df = Dilution factor

5 = Volume (in milliliters) of enzyme assay

0.125 = The change in turbidity produced by 1 unit of of Lysostaphin (as per the Unit Definition)

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

Enzymatic Assay of LYSOSTAPHIN

UNIT DEFINITION:

One unit will reduce the turbidity ($A_{620\text{nm}}$) of a suspension of *S. aureus* cells from 0.250 to 0.125 in 10 minutes at pH 7.5 at 37°C in a 6.0 ml reaction mixture.

FINAL ASSAY CONCENTRATION:

In a 6.00 ml reaction mix, the final concentrations are 50 mM Tris, 145 mM sodium chloride and 0.25 - 1.25 units lysostaphin.

REFERENCE:

Iversen, O.-J. and Grov, A. (1973) *Eur. J. Biochem.* **38**, 293-300

NOTES:

1. Reproducibility of this assay is very dependent on the freshness and quality of the cell suspension.
2. This procedure incorporates the use of a potentially infectious microorganism. Containment and/or disinfection of all such products, solutions and test reagents is imperative. Heat or chemical sterilization of all contaminated reagents and any spills or releases, with subsequent disposal as Biohazardous Waste, must be employed.
3. The $\Delta A_{620\text{nm}}$ at zero time (T_0) should fall between 0.2 - 0.3 absorbance unit.
4. This assay is based on the cited reference.
5. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.