

Technical Data Sheet

CE Standard I Nutrient Agar

Ordering number: 1.07881.0500 / 1.07881.5000

These culture media are suitable for the cultivation of fastidious bacteria; after addition of blood, ascites fluid or serum they can also be used to cultivate *streptococci*, *pneumococci* and *erysipelas* organisms etc. They are employed for the enumeration, isolation and enrichment of bacteria and also as high-grade bases for preparing special culture media.

IVD in vitro diagnosticum - For professional use only

Typical Composition

Peptones	15 g/l
Yeast Extract	3 g/l
NaCl	6 g/l
D(+)-Glucose	1 g/l
Agar-Agar	12 g/l

Preparation

Suspend 37 g/l. Autoclave 15 min at 121 °C.

The appearance of the prepared media is clear and yellowish-brown.

The pH value at 25 °C is in the range of 7.3-7.7.

Experimental Procedure and Evaluation

Depend on the purpose for which the media are used.

Incubation: 24 h at 35 °C aerobically.

Storage

Usable up to the expiry date when stored dry and tightly closed at +15 to +25 °C. Protect from light.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed at +15 to +25 °C.

Specimen

e.g. Blood.

Clinical specimen collection, handling and processing. See general instructions of use.

Quality Control

Control Strains	ATCC #	Incubation	Expected Results
<i>Escherichia coli</i>	25922	24 h at 35 °C	Recovery ≥ 70 %
<i>Staphylococcus aureus</i>	25923	24 h at 35 °C	Recovery ≥ 70 %
<i>Streptococcus pyogenes</i>	12344	24 h at 35 °C	Recovery ≥ 70 %
<i>Streptococcus pneumoniae</i>	6301	24 h at 35 °C	Recovery ≥ 70 %
<i>Listeria monocytogenes</i>	19118	24 h at 35 °C	Recovery ≥ 70 %
<i>Shigella flexneri</i>	12022	24 h at 35 °C	Recovery ≥ 70 %
<i>Enterococcus faecalis</i>	11700	24 h at 35 °C	Recovery ≥ 70 %

Please refer to the actual batch related Certificate of Analysis.



Escherichia coli ATCC 25922

Ordering Information

Product	Cat. No.	Pack size
Standard I Nutrient Agar	1.07881.0500	500 g
Standard I Nutrient Agar	1.07881.5000	5 kg
Standard I Nutrient Agar	1.46256.0020	20 plates
Standard I Nutrient Agar	1.46256.0100	100 plates

Merck KGaA, 64271 Darmstadt, Germany
Fax: +49 (0) 61 51 / 72-60 80
mibio@merckgroup.com
www.merckmillipore.com/biomonitoring

Find contact information for your country at:
www.merckmillipore.com/offices
For Technical Service, please visit:
www.merckmillipore.com/techservice



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Merck Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany. ATCC is a registered trademark of ATCC, Manassas, VA, USA. Lit. No. TN1473EN00