

Detection of Shiga Toxin *Escherichia coli* (STEC) in Large Test Portions of Raw Ground Meat and Vegetables by Immunomagnetic Separation and qPCR

Marion Bouvier¹, Mélissa Canizares¹, Lionel Meyer², and Delphine Thevenot-Sergentet¹
 (1) Université de Lyon, VetAgro Sup, LNR *E. coli* (incl. STEC), Marcy l’Etoile, France,
 (2) Millipore SAS, Molsheim, France.

Introduction


The aim of this study was to evaluate the performance of the Assurance® GDS method combining immunomagnetic separation (IMS) and qPCR for the analysis of large test portions (up to 375 g) of vegetables and meat. Performance was compared with the reference method ISO/TS-13136:2012 and the validated 25 g sample size.

The alternative method includes short, 8 h enrichment in proprietary broth (1/5 dilution) at 41.5 °C followed by screening of *eae*, *stx1* & *stx2*, O157:H7 markers, and other serogroup identification. Following, cultural confirmation was performed on two selective agars by direct plating and using IMS beads.

Methods

A total of 120 samples of meat and vegetables were analyzed using both methods, including 20 uninoculated and 100 samples spiked with stressed cells at positive and fractional levels of different STEC strains.

For the reference method and sample size, enrichment was performed in BPW (1/10 dilution) at 41.5 °C for 18 h.



Pathogen Screening Kit: Assurance® GDS MPX Top 7 STEC (Cat. No. 71015BC)

Results

- Results showed that 70% and 62% of the all samples spiked were detected and confirmed by the alternative method and the reference method, respectively (**Table 1**).
- The RLOD₅₀ ratios (= LOD₅₀ alternative/LOD₅₀ reference) ranged between 0.35 and 0.96 for the four food items tested.
- The performance of the methods and sample sizes were considered equivalent (RLOD₅₀ < 2.5 according ISO 16140-2:2020).
- 100% of PCR hits were confirmed by strain isolation on agar plates.

Conclusion

The alternative method using IMS + PCR technology allowed rapid STEC detection (8 h) of larger sample sizes (375 g) from beef and vegetables with equivalent performance to the reference method.

The IMS step also contributed to improve the cultural confirmation of STEC positive samples.

| Matrix | Serotype | Method | Test Sample | LOD ₅₀ * | ROLD |
|--------------------------------|----------|---------------------------------|-------------|---------------------|-------|
| Lettuce | 0111:H8 | LNR method adapted from TS13136 | 25 g | 0.418 | 0.742 |
| | | Assurance® GDS Method | 375 g | 0.31 | |
| Frozen Minced Meat | 0145:H28 | LNR method adapted from TS13136 | 25 g | .0544 | .035 |
| | | Assurance® GDS Method | 375 g | 0.12 | |
| Spinach | 0157:H7 | LNR method adapted from TS13136 | 25 g | 0.515 | 0.96 |
| | | Assurance® GDS Method | 375 g | 0.522 | |
| Fresh Ground Meat with 20% fat | 0111:H8 | LNR method adapted from TS13136 | 25 g | 0.343 | 0.726 |
| | | Assurance® GDS Method | 375 g | 0.374 | |

Table 1: LOD results.

* LOD and RLOD were calculated using the Excel file on the ISO website (<http://standards.iso.org/iso/16140>)

Key Features

Sample prep

375 g = 25 g

LOD₅₀ Assurance® GDS method
[0.12 to 0.522]
CFU/375 g

Enrichment

41.5 °C

Sufficient growth of STEC at this temperature

Screening & Confirmation

% of recovery

IMS method

100%

Direct spread of the enrichment

94%

Use 2 different selectivity agar plates

Assurance® GDS method shows similar data to the reference method

Request more information
SigmaAldrich.com/assurance-gds-info

MilliporeSigma
400 Summit Drive
Burlington, MA 01803

We have built a unique collection of life science brands with unrivalled experience in supporting your scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

