

## Product Information

### FSL-Tyr1

Catalog Number **F7558**

Storage Temperature  $-20\text{ }^{\circ}\text{C}$

Synonym: FSL-Tyr1(forY)-SC2-L1

#### Product Description

Molecular formula:  $\text{C}_{91}\text{H}_{140}\text{N}_{16}\text{O}_{33}\text{PS}$

Molecular weight: 2017.15

FSL-Tyr1 is a KODE™ technology construct designed to label hydrophobic surfaces including living cells and virions with  $^{125}\text{I}$ . All KODE FSL constructs consist of three essential designable features:

- functional component (F)
- spacer (S)
- diacyl lipid (L)

FSL-Tyr1 is comprised predominantly of an *N*-formyl-tyrosine monomer, which can be radiolabeled with  $^{125}\text{I}$ , representing F, conjugated via a carboxymethylglycine based linker (SC2) to an activated adipate derivative of dioleoylphosphatidylethanolamine (L). All FSL constructs disperse in biocompatible media, and spontaneously and stably incorporate into cell/virion membranes. Cells and virions modified with KODE constructs are known as kodeocytes<sup>1</sup> or kodevirions,<sup>2</sup> respectively, and usually maintain their normal vitality and functionality.

Following  $^{125}\text{I}$  radiolabeling, FSL-Tyr1 has been specifically designed to create  $^{125}\text{I}$  labeled kodeocytes or kodevirions.<sup>2</sup>

#### Reagents Required but Not Provided for $^{125}\text{I}$ Radiolabeling Procedure

- Iodination tubes (Cat. 28601, Thermo Fisher Scientific Inc.)
- $\text{Na}^{125}\text{I}$  (Cat. NEZ-033A, Perkin-Elmer)
- Microcon-10 filter device (Millipore Corporation)
- Tris-HCl-NaCl buffer (25 mM Tris-HCl, pH 7.5, with 0.4 M NaCl)

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

Stock Solution is prepared by reconstituting the product at a concentration of 2 mg/mL in saline or PBS. Buffered solutions are preferred for long-term storage. The product should not be reconstituted in water, unless used immediately as the product is unstable when stored in water.

The 2 mg/ml Stock Solution can be frozen in aliquots for later usage. Thawed product should be briefly sonicated before use.

#### Storage/Stability

Store unopened product at  $-20\text{ }^{\circ}\text{C}$ . Store the Stock Solution in aliquots at  $-20\text{ }^{\circ}\text{C}$ . Avoid repeated freezing and thawing of solutions. Solutions in PBS, pH 7, can be stored at  $2\text{--}8\text{ }^{\circ}\text{C}$  for up to 2 weeks.

#### Procedure

Iodination of FSL-tyrosine to FSL- $^{125}\text{I}$

1. Wet interior surface of iodination tube with 1 mL of Tris-HCl-NaCl buffer and decant.
2. Add 100  $\mu\text{L}$  of Tris-HCl-NaCl buffer to the tube, followed by 1 mCi of  $\text{Na}^{125}\text{I}$ . Incubate at room temperature for 6 minutes with swirling every 30 seconds.
3. Transfer the activated  $\text{Na}^{125}\text{I}$  solution to a new centrifuge tube containing 150  $\mu\text{L}$  of FSL-Tyr1 Stock Solution (0.3 mg). Incubate at room temperature for 8 minutes, with swirling every 30 seconds.
4. Transfer reaction solution to a Microcon-10 filter device. Rinse the reaction tube with 250  $\mu\text{L}$  of Tris-HCl-NaCl buffer and add to reaction solution. Centrifuge the filter device for 30 minutes at  $15,000 \times g$ . Collect flow-through solution.
5. Rinse the reaction tube with an additional 300  $\mu\text{L}$  of Tris-HCl-NaCl buffer and add to Microcon-10 filter device. Centrifuge for 30 minutes at  $15,000 \times g$ . Collect second flow-through solution.

The flow-through solutions contain the radiolabeled FSL-Tyr1 (FSL- $^{125}\text{I}$ ). Quantify all fractions using a dose calibrator prior to pooling, then aliquot solution and store at  $-20\text{ }^{\circ}\text{C}$ .

FSL radioiodination can be confirmed by SDS-PAGE on a 15% acrylamide gel in non-reducing sample buffer.

Cell/Virion labeling - Add 1 volume of FSL-<sup>125</sup>I (1–100 µg/mL diluted in PBS) to 1 volume of cells/virions. Incubate for 1 hour at 37 °C. Wash with PBS or other appropriate buffer (may be optional). Store kodecytes/kodevirions in serum free medium.  
Note: Rate of FSL insertion is primarily determined by FSL concentration, incubation time, and temperature.

## References

1. Henry, S.M., Modification of red blood cells for laboratory quality control use. *Curr. Opin. Hematol.*, **16**, 467-472 (2009).
2. Hadac, E.M. et al., Fluorescein and radiolabeled Function-Spacer-Lipid constructs allow for simple *in vitro* and *in vivo* bioimaging of enveloped virions. *J. Virol. Methods*, **176**, 78-84 (2011).

## KODE Construct End User License (Non-Commercial Research Purposes)

KODE Constructs are supplied under license from KODE Biotech Limited. As the entity, which has purchased a KODE Construct, you have the non-transferable right to use the purchased KODE Construct for research conducted by you, for your own internal purposes, in accordance with the provisions of this document. You may not use the purchased KODE Construct, or any composition comprising the purchased KODE Construct (including modified cells or "kodecytes"), for any commercial purpose. 'Commercial purposes' mean any activity for which you obtain a financial or in-kind benefit from another entity, including manufacture of products for sale or other disposition, and the conduct of research or other services for or on behalf of a third party or where a third party has any option or right of refusal over research outcomes. Transfer to a third party is permitted where the third party is a genuine scientific collaborator and has agreed in writing to be bound by the provisions of this document prior to the transfer. You are responsible for ensuring that your personnel and scientific collaborators comply with the restrictions provided in this document. Your purchase of the KODE Construct does not give you any license of, or interest in, any intellectual property rights of KODE Biotech Limited other than as may be necessary for your internal research as permitted by this document. The KODE Construct is supplied "as is" and neither KODE Biotech Limited nor the supplier gives any warranty that the purchased KODE Construct is fit for any particular purpose or that it has any particular qualities or characteristics.

KODE is a trademark of KODE Biotech Materials Limited.

VNC,MAM 01/12-1