

Technical Data Sheet

CE Standard I Nutrient Broth

Ordering number: 1.07882.0500

These culture media are suitable for the cultivation of fastidious bacteria; after addition of blood, ascites fluid or serum they can also be used to cultivate *streptococci*, *pneumococci* and *erysipelas* organisms etc. They are employed for the enumeration, isolation and enrichment of bacteria and also as high-grade bases for preparing special culture media.

Standard I Nutrient Broth is a nonselective medium, useful in routine cultivation of microorganisms.

VD in vitro diagnosticum - For professional use only

Typical Composition

Peptones	15 g/l
Yeast Extract	3 g/l
NaCl	6 g/l
D(+)-Glucose	1 g/l

Preparation

Suspend 25 g/l. Autoclave 15 min at 121 °C.

The appearance of the prepared media is clear and yellowish-brown.

The pH value at 25 °C is in the range of 7.3-7.7.

Experimental Procedure and Evaluation

Depend on the purpose for which the media are used.

Incubation: 24 h at 35 °C aerobically.

Storage

Usable up to the expiry date when stored dry and tightly closed at +15 to +25 °C. Protect from light.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed at +15 to +25 °C.

Specimen

e.g. Blood.

Clinical specimen collection, handling and processing. See general instructions of use.

Quality Control

Control Strains	ATCC #	Inoculum CFU	Incubation	Expected Results
<i>Escherichia coli</i>	25922	10-100	24 h at 35 °C	Growth good to very good
<i>Staphylococcus aureus</i>	25923	10-100	24 h at 35 °C	Growth good to very good
<i>Staphylococcus aureus</i>	6538	10-100	24 h at 35 °C	Growth good to very good
<i>Staphylococcus epidermidis</i>	12228	10-100	24 h at 35 °C	Growth good to very good
<i>Streptococcus pyogenes</i>	12344	10-100	24 h at 35 °C	Growth good to very good
<i>Streptococcus pneumoniae</i>	6301	10-100	24 h at 35 °C	Growth good to very good
<i>Listeria monocytogenes</i>	19118	10-100	24 h at 35 °C	Growth good to very good
<i>Shigella flexneri</i>	12022	10-100	24 h at 35 °C	Growth good to very good
<i>Enterococcus faecalis</i>	11700	10-100	24 h at 35 °C	Growth fair to good

Please refer to the actual batch related Certificate of Analysis.

Ordering Information

Product	Cat. No.	Pack size
Standard I Nutrient Broth	1.07882.0500	500 g

Merck KGaA, 64271 Darmstadt, Germany
Fax: +49 (0) 61 51 / 72-60 80
mibio@merckgroup.com
www.merckmillipore.com/biomonitoring

Find contact information for your country at:
www.merckmillipore.com/offices
For Technical Service, please visit:
www.merckmillipore.com/techservice



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Merck Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany. ATCC is a registered trademark of ATCC, Manassas, VA, USA. Lit. No. TN1474EN00