

# Technical Bulletin

# Microbe Selection for the SER-TAIN™ **Process Validation**

One of the critical aspects of designing an extraneous agent inactivation validation is the selection of specific microbes to evaluate. Those used in the SER-TAIN™ Gamma Irradiation Serum Validation were carefully selected to represent a wide range of organisms that are potential contaminants of serum products arising from either the starting material itself or from the processing environment.

When considering potential endogenous microbes in bovine serum, viruses are a main concern. Several viruses that can affect cattle are listed in Table 1, all of which are represented in the SER-TAIN™ validation.

BVD, PI<sub>3</sub>, IBR, BTV and BLV are all bovine viral species of significant industrial importance. Although Porcine Parvovirus (PPV) does not cause disease in cattle under natural conditions, it has been isolated from bovine herds2. MvM was intended as a potential environmental contaminate.

Viruses of the same type and similar physical characteristics show comparable vulnerability to inactivation by gamma

radiation. This allows model viruses to be used as a substitute for viruses that are difficult or hazardous to grow and titer in culture. The use of viral models in validation studies is a common practice and is recommended by the Food and Drug Administration (FDA) in evaluating the safety of human therapeutics<sup>3</sup>. For example, BLV is utilized in studies as a model for Human Immunodeficiency Virus (HIV).

Model viruses were also used in the SER-TAIN™ validation. Consistently propagating Bovine Leukemia Virus (BLV) in vitro is nearly impossible. By using Feline Leukemia Virus (FeLV), a virus that is physically similar and considerably less fastidious, the amount of protection gamma radiation provides against contamination by BLV can be accurately assessed. Based on the information gathered about Pl<sub>3</sub> (of the family Paramyxoviridea), assessments can be made about the inactivation of other paramyxoviruses, such as Bovine Respiratory Syncytial Virus (BRSV). Similarly, the use of PPV as a model provides confirmation of protection against Bovine Parvovirus (BPV).

Table 1 Viruses Used in the SER-TAIN™ Gamma Irradiation Serum Validation				
Virus	Family	Characteristics	Validation Virus	
Bovine Viral Diarrhea (BVD)	Flaviviridae	ss RNA enveloped	BVD	
Parainfluenza Type 3 (PI <sub>3</sub> )	Paramyxoviridae	ss RNA enveloped	PI₃	
Infectious Bovine Rhinotracheitis (IBR)	Herpesviridae	ds DNA enveloped	IBR	
Bluetongue (BTV)	Reoviridae	ds RNA non-enveloped	BTV	
Bovine Leukemia (BLV)	Retroviridae	ss RNA enveloped	FeLV	
Porcine Parvovirus (PPV)	Parvoviridae	ss DNA PPV non-enveloped		
Minute Virus of Mice (MvM)	Parvoviridae	ss DNA non-enveloped	M∨M	

#### **United States**

SAFC Biosciences, Inc. 13804 W. 107th St Lenexa, Kansas 66215 USA

Toll free-USA E-mai

£1 913-469-5580 1 800-255-6032 +1 913-469-5584

### Europe

SAFC Biosciences Ltd. Smeaton Road, West Portway Andover, Hampshire SP10 3LF UNITED KINGDOM

+44 (0)1264-333311 Phone +44 (0)1264-332412 F-mail info-eu@sial.com

#### Asia Pacific

SAFC Biosciences Ptv. Ltd. 18-20 Export Driv Brooklyn, Victoria 3025 AUSTRALIA

Phone Toll free-AUS E-mail

+61 (0)3-9362-4500 1 800-200-404 +61 (0)3-9315-1656 info-ap@sial.com

www.safcbiosciences.com

Non-viral contaminants from the environment, such as bacteria, fungi and mycoplasma, are also a major concern for users of serum products. Table 2 shows those organisms used in the SER-TAIN™ Validation.

## Table 2 Organisms Used in the SER-TAIN™ Gamma Irradiation Serum Validation

Organism	Type	Characteristics	
Escherishia coli	Bacterium	Gram Negative	
Bacillus pumilus	Bacterium	Gram Positive	
Candida albicans	Fungus	Dimorphic	
Acholeplasma laidlawii	Mycoplasma	Bovine specific	
JH Strauss, Strain MS2	Bacteriophage	Phage 15597-B1	

The bacteria and fungi species used were based on those recommended for sterility testing in the United States Pharmacopeia⁴. Although there are many types of mycoplasma, *A. laidlawii* was chosen because it is the variety that specifically affects cattle. Additionally, an *E. coli* specific bacteriophage (JH Strauss, Strain MS2) was included as a potential contaminant from the environment. The use of these organisms provides a look at the protection SER-TAIN™ processing provides against a broad range of potential environmental agents.

The unique SER-TAIN™ radiation process provides assurance against the presence of microbial contaminants in serum products. The validation of this process has been solidified by including a wide spectrum of organisms that have the greatest potential to be present in serum products and are of major industry concern. Our experience with the SER-TAIN™ process allows the validation to be expanded as new extraneous agents emerge that have the potential for transmission through bovine serum.

For more information about this subject or other SAFC Biosciences' products and services, please call our Technical Services department.

### References

- 1. Fenner, F.J., Gibbs, E.J., Murphy, F. A., et al. Veterinary Virology, Academic Press, Inc., 1993, pp. 622-623.
- 2. Ibid., pg. 308.
- 3. "Guide To Inspections of Viral Clearance Processes For Plasma Derivatives", Federal Drug Administration Internet homepage, www.fda.gov/ora/inspect\_ref/igs/viralcl.html.
- 4. Section 71 Sterility Tests, United States Pharmacopeia, 1995 Edition, pp. 1686-1687.

#### Warranty, Limitation of Remedies

SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, as copy of which is available upon request. All OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request

SER-TAIN™ is a trademark of SAFC Biosciences, Inc.

© 2006 SAFC Biosciences, Inc.

Issued April 2006 T035 0103 1005 1205