

Application Note

HepaRG™ cells and their application to diverse endpoints

Nicola J. Hewitt. SWS, Erzhausen, Germany.

Abstract

Reproducible measurements of drug uptake, metabolism and toxicity require robustly functioning hepatic cells with stable, high expression of transporters and metabolic enzymes. Fulfilling this need are HepaRG™ cells. These widely published, terminally differentiated hepatocellular carcinoma cells have facilitated numerous studies of uptake, metabolism, and disposition of drug candidates. In this review, we describe the current uses of HepaRG™ cells and their position on the market as a hepatic model used by the pharmaceutical, chemical and cosmetics industries.



Figure 1.

HepaRG™ cells are an easy, convenient, functionally relevant model for diverse applications traditionally requiring hepatocyte model systems.

Why are hepatocyte models needed?

Although it is crucial to assess the effects of preclinical drug candidates on the liver, efficient completion of this step is hindered by the limited, sporadic supply of primary human hepatocytes (PHHs), the traditionally used hepatic model system. PHHs have little, if any, proliferative capacity, a short lifespan and lose their differentiated functions when cultured even for a few days. There are also large donor-specific variations in initial and longer-term functions and enzyme activities

(especially CYPs). Therefore, a long-lived hepatic model that features good proliferation, stable phenotype, low variability and has the advantage of an abundant supply would be an extremely useful tool as an alternative to PHHs. There have been numerous reports of potential replacements for PHHs but until now no single cell model has exhibited a phenotype or functions close enough to PHHs to be sufficiently predictive¹. The origin and characterization of HepaRG™ cells is reviewed here, along with their application to metabolism and toxicology assays (summarized in Table 1).

Table 1.

Characteristics and applications of HepaRG™ cells and other cell lines (characteristics of HepaRG™ cells are referring to those in a differentiated state, unless otherwise specified).

Characteristic	HepaRG™ cells	Other cell lines
Hepatic-specific markers	<ul style="list-style-type: none"> • Good carbohydrate metabolism, produce albumin, and eliminate galactose and sorbitol at comparable rates to PHHs. • Express aldolase B (20% of PHH), cytokeratin 8 and 18 and hepatocyte-specific antigen, CD26, and E-cadherin (markers of apical and lateral polarized phenotype), ZO-1 (tight junctions marker), CD49a, and are p53 competent. • Poor urea production. • As with PHHs, HepaRG™ cells are negative for α-fetoprotein. 	<ul style="list-style-type: none"> • Lack a variable and substantial set of liver-specific functions, making them unsuitable as representative of <i>in vivo</i> PHHs. • No aldolase B mRNA in HepG2 cells and albumin and haptoglobin mRNA levels are much lower than those in HepaRG™ cells. • Poor urea production. • α-Fetoprotein is highly expressed.
Enzyme characteristics (See Table 2)	<ul style="list-style-type: none"> • Differentiated HepaRG™ cells retain many drug metabolizing enzymes (DMEs) at levels comparable to those in PHH cultures. • CYP2D6 and CYP2C9 polymorphic. 	<ul style="list-style-type: none"> • HepG2 cells, although retaining some DMEs, have low (CYP1A2, CYP2B6) or absent (CYP2C9, CYP2E1 and CYP3A4) CYP expression apart from foetal forms (CYP1A1 and CYP3A7). • CYP activities are also low or absent.
DME regulation, transcription factors and nuclear receptors	<ul style="list-style-type: none"> • The transcription factor, AhR, and nuclear receptors, CAR, PXR and PPAR, responsible for the regulation of many DMEs and endogenous compounds are all expressed in high levels in differentiated HepaRG™ cells. 	<ul style="list-style-type: none"> • AhR, PXR and PPAR are all expressed in HepG2 cells but (with the exception of PPAR) are present at much lower levels than differentiated HepaRG™ cells or PHHs. • CAR has never been reported to be expressed at high levels in any other hepatoma cell line.
Drug transporters	<ul style="list-style-type: none"> • Express functional sinusoidal and canalicular transporter functions and regulation. 	<ul style="list-style-type: none"> • Uptake transporters expression absent except for MDR1, MRP1, and breast cancer resistance protein, which were expressed similar levels to PHHs. • No detectable expression of OATP-C, NTCP, OCT1 or BSEP.
Induction of CYPs	<ul style="list-style-type: none"> • CYP activities are readily induced. • Can be used to classify compounds as inducers or non-inducers of CYP3A4 as well as to predict the extent of induction <i>in vivo</i> in humans. 	<ul style="list-style-type: none"> • CYP1A2 and CYP3A4 induced in Fa2N4 cells but not CYP2B6 because they lack CAR. EC₅₀ values for CYP3A4 mRNA induction in Fa2N4 cells by a number of compounds were ~3-fold lower than in PHHs but the EC₅₀ for rifampicin was 10-fold higher. • HepG2 cells are not suitable for CYP3A induction assays due to a lack of many DMEs, many liver-specific functions and CAR.
Cytotoxicity	<ul style="list-style-type: none"> • Sensitive to hepatotoxicants whose toxicity is metabolism- and/or transporter-dependent. • Suitable for higher throughput screening. • The only known hepatic model to predict steatosis. • Best model for phospholipidosis. 	<ul style="list-style-type: none"> • HepG2 cells fail to detect a number of known hepatotoxicants due to a lack of DMEs and transporters. • HepG2 cells are poor models for steatosis and phospholipidosis.
Genotoxicity and carcinogenicity	<ul style="list-style-type: none"> • Can be adapted to micronucleus and Comet assays to detect direct and bioactivated genotoxicants. • Gene expression analysis correlates well with known <i>in vivo</i> effects of carcinogens. 	<ul style="list-style-type: none"> • HepG2 cells are more relevant than rodent cell lines but they still lack bioactivating enzymes, resulting in the potential for false negative results. • Different HepG2 clones do not give reproducible results for dietary genotoxins.
Inflammation	<ul style="list-style-type: none"> • Responsive to lipopolysaccharide. • Express important inflammatory mediators. • Effects of IL-1β on transporter function similar to PHHs. Represent a good model for sepsis and cholestasis. 	<ul style="list-style-type: none"> • HepG2 cells have been used to model nitric oxide and IL-8 release in response to inflammatory stimuli; however, the responses may be different from those <i>in vivo</i> due to low DME activities.
Virus infection	<ul style="list-style-type: none"> • HepaRG™ cells retain functional and efficient signaling pathways involved in pro-inflammatory effects. • Uniquely susceptible to viral infection. 	<ul style="list-style-type: none"> • Hepatoma cells, such as HepG2, have a limited use in viral infection research since they have impaired antiviral responses. • Not susceptible to hepatitis B virus infection.
3D bioartificial livers	<ul style="list-style-type: none"> • Suitable for 3D culture and use in bioreactors. • Morphology and DME functions similar to that of PHHs in 3D. 	<ul style="list-style-type: none"> • Not suitable for bioreactors due to poor DME and transporter functions.

Characterization of HepaRG™ cells

Origin of HepaRG™ cells

The HepaRG™ cell line was established from a tumor of a female patient suffering from chronic hepatitis C infection and hepatocarcinoma. Although they are likely to have originated from ductular structures (rather than mature hepatocytes or bile ducts) associated with long-term HCV infection, HepaRG™ cells do not contain any part of the HCV genome or express any HCV protein². These cells are different from HepG2 cells in that the latter are derived from hepatoblastoma cells, whereas HepaRG™ cells are derived from a hepatocarcinoma. When passaged at low density, they are able to recover and differentiate into both hepatocytes and biliary epithelial cells and are thus considered to be progenitor cells². HepaRG™ cells do not grow in serum-free medium and show moderate anchorage-independent growth in soft agar. When injected into mice, p53-competent HepaRG™ cells did not form tumours and repopulated damaged livers in mice, and demonstrated differentiated hepatocyte functions 6 weeks after implantation^{3,4}.

Hepatocyte-like differentiation of HepaRG™ cells

In vitro, maximum cell differentiation is reached after 2 weeks of exposure to DMSO and 40 to 50% of the confluent cell population are hepatocyte-like in nature, with morphology close to that of PHHs. Differentiated HepaRG™ cells stop proliferating and retain their hepatocyte-like features. Genes up-regulated during differentiation are those related to cell cycle inhibition, increased susceptibility to apoptosis, innate immunity and liver-enriched transcripts involved in lipid homeostasis and drug metabolism⁵. The expression of different hepatic nuclear factors (HNFs) involved in hepatic-specific gene expression changes as the cells grow and differentiate³. The cells surrounding the hepatocyte-like cells are biliary epithelial cells.


Although HepaRG™ cells do not produce urea (due to poor or disturbed nitrogen elimination via the urea cycle), they are able to regulate carbohydrate metabolism (glycogenolysis and/or gluconeogenesis), produce lactate (a product of anaerobic metabolism) and albumin, and eliminate galactose and sorbitol at comparable rates to PHHs. These features, along with the limited chromosomal rearrangements in these cells, make HepaRG™ cells much more differentiated than any other human hepatocyte cell line reported thus far⁶.

Gene expression profiling has revealed that HepaRG™ cells are remarkably close to certain hepatocyte populations which have a low expression of several genes involved in xenobiotic (CYP3A4 and CYP2C9), lipid, and carbohydrate metabolism⁷. It is therefore suggested that HepaRG™ cells behave as a PHH population, and can be classified as an "average human hepatocyte population." However, some drug metabolizing enzymes (DMEs) may be expressed at different levels than in PHHs and are sensitive to DMSO (see Figure 2), making direct comparisons with some drugs not possible.

Expression of drug metabolizing enzymes (DMEs)

There have been a number of comprehensive analyses of Phase 1 and 2 DME expression in HepaRG™ cells. A summary of the DMEs, regulation pathways and transporter characteristics is shown in Table 2. A genome-wide gene expression profiles analysis showed that for most genes encoding phase 1 and 2 DMEs and drug transporters, the differences between HepaRG™ cells and PHHs were much smaller than between HepG2 cells and PHHs⁸. The expression of DMEs is greater when they are cultured for longer times and high densities⁹. Moreover, inclusion of 2% DMSO in the culture medium caused the expression of some DMEs (CYP2B6, CYP2C9, CYP2E1, CYP3A4, UGT1A1, GSTA4, GSTM1, and GSTA1/A2) to be increased further, such that the levels were comparable with those measured in 3 to 5 day cultures of PHHs (which themselves are only 10 to 60% of the expression levels in freshly isolated hepatocytes, depending on the CYP). The expression of the CYP enzymes is maintained in the differentiated HepaRG™ cells for at least one month when cultured in the presence of DMSO^{9,10,11}. Removal of 2% DMSO reverses the effect; whereby the expression and activities of major CYPs, especially CYP3A4, decreases significantly¹².

Figure 2. HepaRG™ cells express phase 1 (A.) and phase 2 (B.) drug metabolism enzymes at higher levels than do HepG2 cells. Gene expression studies were conducted to measure the transcript levels of DMEs in HepaRG™ and HepG2 cells, shown here as a percentage of the transcript levels in freshly isolated PHHs. (Data from Aninat C. et al, 2006.⁹)



CYP2D6 expression in highly differentiated HepaRG™ cells is less than 10% of PHHs, suggesting this donor is polymorphic for this CYP¹³. These cells are also polymorphic for CYP2C9 which is reflected in the relatively low CYP2C9-mediated metabolism of diclofenac and tolbutamide^{9,12}. CYP2E1 and CYP1A2 activities are also relatively low in the differentiated HepaRG™ cells cultured with or without DMSO when compared with freshly isolated PHHs^{9,12}. The levels of these CYPs can be increased in HepaRG™ cells without altering the levels of other highly expressed CYPs (such as CYP3A4) using a transient transfection method

involving electroporation, as demonstrated for CYP2E1¹⁴. CYP4F3B and CYP4A11 are expressed and functional in HepaRG™ cells^{15,16}. Although these are not drug metabolizing CYPs, they are important to note because they catalyse ω -hydroxylation, a pathway accounting for 5–10% of total fatty acid oxidation¹⁷. The presence of these enzymes is especially relevant to the use of these cells as a model for steatosis (See section entitled "Drug-induced liver injury").

Table 2. DMEs, transporters and transcription factors in HepaRG™ cells and other cell lines

Characteristic	HepaRG™ cells		HepG2 cells	
	Basal expression	Inducibility	Basal expression	Inducibility
Phase 1 DMEs				
CYP1A1	+	✓	+	✓
CYP1A2	-/+	✓	+	✓
CYP2B6	+++	✓	+	Poor
CYP2C8	+	✓	ND	ND
CYP2C9	+(PM)	✓	-	ND
CYP2C19	+++	✓	ND	ND
CYP2E1	+	✓	-	x
CYP3A4	+++	✓	-	x
CYP2D6	+(PM)	NI	+	NI
CYP4F3B/CYP4A11	++	✓	ND	ND
FMO1/2/3	+++	x	++	x
Aldehyde dehydrogenases	+++	ND	ND	ND
Alcohol dehydrogenases	+++	ND	ND	ND
Phase 2 DMEs				
UGT1A1	+++	✓	+	ND
GSTM1	++	ND	+	ND
GSTA4	+++	ND	+	ND
GSTA1/A2	+++	✓	+	ND
SULT1B1/SULT1E1	+++	x	++	x
NAT1/2	+++	ND	ND	ND

Characteristic	HepaRG™ cells		HepG2 cells	
	Basal expression	Inducibility	Basal expression	Inducibility
DME regulation transcription factors				
AhR	+++	ND	++	ND
CAR	++	ND	+/-	ND
PXR	+++	ND	++	ND
PPARs	+++	ND	+++	ND
Drug transporters				
MRP1/MPR3	+++	ND	ND	ND
MRP2	++	✓	++	✓
MDR1	+++	✓	++	x
MDR3	++	ND	ND	ND
BSEP	+	ND	ND	ND
BCRP	+++	✓	+++	ND
OCT1	+++	x	-	x
OATP1B1/OATP1B3/OATP2B1	+	ND	-	ND
OAT2	+++	ND	++	ND
OCT1	+++	ND	-	x
NTCP	++	x	+	x
PepT1	++	ND	ND	ND

Taken from Kanebratt and Andersson (12,19); Hart et al. (8); Wilkening et al. (69) and Gerets et al. (70).
 - = absent, + = present at low levels, ++ = present at measureable levels, +++ = present at levels comparable to PHHs,
 ✓ = inducible; PM = polymorphic; NI = not inducible in PHHs; ND = Not determined, x = little or no induction measured.

DME activity

Because relevant DMEs are generally expressed and active in HepaRG™ cells, this model has been successfully used for studies of enzyme induction and inhibition. A panel of DME activities in HepaRG™ cells have also been determined (Table 2)^{9,11,18,19,20}. The levels of activities correlate with the expression level of each CYP. A major difference between HepaRG™ cells and PHHs is that CYP1A1 is continuously expressed in the former; but is attributed to the transformed state of HepaRG™ cells rather than to the presence of biliary cells¹⁰. CYP2C19 activity is much higher in differentiated HepaRG™ cells than PHHs and CYP2C8 and CYP2C9 activities were also present at levels equivalent to those in PHHs used for comparison in the study¹⁸. Likewise, bupropion 6-hydroxylation (CYP2B6) and testosterone 6β-hydroxylation (CYP3A4) are present in HepaRG™ cells at levels reported to be present in PHH cultures^{9,18}. In accordance with the low expression of CYP2A6 and CYP2E1, metabolism of coumarin and chlorzoxazone was equally low or lacking. The characterization of Phase 2 DMEs in HepaRG™ cells is relatively limited compared to that of CYPs. GSTs are expressed at levels comparable to PHHs (Table 2) and were stable for up to 1 month^{9,10}. GSTA1/2 and UGT1A1 expression was induced by omeprazole, phenobarbital and rifampicin in low- and high-density cultures.

Transcription factors

Any cell model proposed as an alternative to PHHs should express transcription factors involved in the regulation of DMEs and transporters. The nuclear receptor, CAR, is present in HepaRG™ cells and the addition of 2% DMSO causes CAR to be increased further; which is important because this nuclear receptor is responsible for "phenobarbital-like" induction of CYP2B6, an enzyme not expressed in Fa2N4 and HepG2 cells^{9,21}. FDA compliance requires that CAR-mediated CYP2B6 induction should be demonstrated by the hepatic test model²². Other important transcription factors include AhR, PXR and PPAR, and all are expressed in high levels in differentiated HepaRG™ cells. In contrast to HepaRG™ cells and PHHs, HepG2 cells lack or have only low expression of DMEs and express less than 50% of PHH levels of AhR, CAR and PXR.

Drug transporters

The presence of drug transporters in HepaRG™ cells was investigated at the mRNA level and compared with the levels in PHHs and HepG2 cells (Table 2)¹². The expression of many efflux transporters was equivalent to or higher than PHHs; whereas mRNA levels of bile salt export pump were at 5- to 100-fold lower than in PHHs. The expression levels of the uptake transporters were also 5-fold lower than in PHHs. In order to achieve significant expression of some liver-specific transporters (OCT1, OATP-C, NTCP and BSEP) the HepaRG™ cells had to be grown to confluence. Functional bile canalicular structures in differentiated HepaRG™ cells have been demonstrated using fluorescent substrates^{3,10}.

Applications

Initially, HepaRG™ cells required a number of weeks of culture to bring them to a differentiated state, which was considered somewhat of a disadvantage in terms of end user planning and resources; however, now HepaRG™ cells are available as cryopreserved differentiated cells (i.e. they are differentiated in 2% DMSO and then frozen as a suspension). The freeze/thaw process does not alter their functional activities and the inter-batch reproducibility is excellent; therefore, these cells are ready to use in diverse applications. Some of these are outlined below.

Metabolism studies

PHHs are the standard by which to measure metabolism, since they have a full complement of DMEs, cofactors and transporter proteins²³. Hepatoblastoma cell lines have low or lack important DMEs required and therefore are poor substitutes for PHHs. The high levels of relevant DMEs and transporters present in HepaRG™ cells makes them a much better model than any other cell line currently available. The relative mRNA content of CYPs in HepaRG™ cells is similar to that in PHHs, suggesting that HepaRG™ cells represent a model which is not overabundant in a single CYP isoform and therefore representative of the distribution of CYPs present in PHHs¹². The intrinsic clearance of a number of compounds in differentiated HepaRG™ cells has been compared with that in PHHs^{12,24,25}. In line with basal CYP3A4 activities, midazolam clearance (CYP3A4) was much higher in differentiated HepaRG™ cells than PHHs, but was equivalent in HepaRG™ cells which had been cultured without DMSO. The clearance of compounds metabolized by UGTs, CYP1A2, CYP2C19, CYP3A4, and UGT1A4 in differentiated HepaRG™ cells was comparable to their clearance in PHHs. While the prediction of *in vivo* intrinsic clearance was mostly within 2-fold, there was a greater under-prediction of more highly cleared compounds by HepaRG™ cells, although this observation is also true for PHHs and is thought to be an inherent feature of cell-based clearance models²⁵. The clearance of compounds metabolized mainly by CYP2D6 and CYP2A6 was lower in HepaRG™ cells than in PHHs, but correlated with the lower mRNA expression of the DMEs responsible for the metabolism of these compounds. The long-term stability of CYPs in differentiated HepaRG™ cells makes them a suitable model for metabolism of low clearance compounds which require at least one day of incubation.

Enzyme inhibition

HepaRG™ cells are unique among hepatocyte cell lines in that they can be used for CYP inhibition studies, because they have sufficient levels of DMEs. This is not possible for cell lines which already have low or no activities.

Early studies showed a very good correlation between the IC_{50} values of eight inhibitors of different CYPs in HepaRG™ cells and PHHs, suggesting that HepaRG™ cells represent a promising model for this important endpoint¹⁸.

Enzyme induction

HepaRG™ cells represent a relevant alternative model for induction since they express functional transcription factors involved in the induction of the main CYP enzymes, as well as good DME activities and transporter functions^{9,26}. Another advantage of a cell model that has reproducible and marked CYP induction responses is that the effects of different compounds tested at different times can be directly compared, and inducers with low induction responses can be detected, thus avoiding false negative results.

HepaRG™ cells respond to prototypical inducers of CYP1A1, CYP1A2, CYP2B6, CYP2C8, CYP2C9, CYP2C19, and CYP3A4 (Table 2)^{10,18,19}. As with PHHs, CYP2D6 activities were not inducible in HepaRG™ cells¹⁸. The CYP induction responses were also shown to be stable over one month¹⁰. Importantly, responsiveness to prototypical inducers of CYPs was similar in different passages of HepaRG™ cells; therefore, large amounts of these cells can be produced at later passages which retain the same fold-induction values as earlier passages. The inter-batch reproducibility in these induction studies was shown to be very good¹⁹. The known lower level of induction of CYP2C9 compared with that of CYP3A4 *in vivo* was also reflected in HepaRG™ cells (CYP2C9 was induced by ~2-fold and CYP3A4 by 24-fold by rifampicin). CYP3A4 was also inducible by rifampicin in HepaRG™ cells cultured in conditions resulting in sub-maximal basal CYP3A4 activities. There was no CYP3A4 induction response in fully differentiated HepaRG™ cells, most likely because the basal activities had reached the maximal level obtainable (a normal phenomenon when maximum transcription activity is reached)²⁷. Kanebratt and Andersson¹² recommended the removal of DMSO for one day to decrease CYP3A4 activities before treatment with test compounds and prototypical inducers. With this in mind, HepaRG™ cells which are cultured under conditions that result in lower basal activities will likely result in induction responses that are much higher than in PHHs⁹.

There are reports of an excellent correlation between the *in vitro* and *in vivo* induction of CYP3A4 selected compounds. These predictions of *in vivo* induction were determined using different calculation models, namely the AUC/F₂ prediction model, the Relative Induction Score model and the "Relative Factor" model^{12,28,29}. EC_{50} values for CYP3A4 induction in HepaRG™ cells treated with rifampicin also correlated well with those reported in

PHHs and PXR reporter gene assays³⁰. The time-dependent inhibition of CYP3A4 could also be monitored in HepaRG™ cells by measuring the ratio between the mRNA and activity of this CYP. The prediction model incorporating opposing effects of induction and inhibition was developed further using PHHs and HepaRG™ cells²⁸. Furthermore, these studies underlined the fact that HepaRG™ cells and PHHs could be used interchangeably. A method measuring CYP1A2, CYP2B6 and CYP3A4 induction, using a substrate cocktail incubation for activity and real-time PCR for mRNA detection, has been described in detail by Andersson³¹. The method for CYP induction using HepaRG™ cells is currently part of a validation phase at the European Centre for the Validation of Alternative Methods³².

As with basal activities, the measurement of Phase 2 enzyme induction in HepaRG™ cells and PHHs is still relatively limited³³. Phenobarbital, rifampicin and omeprazole increased both UGT1A1 and GSTA1/A2 transcripts in differentiated HepaRG™ cells, which was comparable to the level of induction of these DMEs measured in PHHs¹⁰. Likewise, the levels of induction of the transporters, MDR1, MRP2 and BCRP, were also much lower than CYP induction¹⁰ and comparable with the induction levels measured in PHHs³⁴.

Drug-induced liver injury

The stable metabolic activities in HepaRG™ cells make them ideal for long-term and repeated dose toxicity studies in which the toxic effects are metabolism-dependent and/or only evident after days or weeks^{12,35}. In addition, the presence of both biliary cells and hepatocytes provides information as to whether toxicity is specific to one cell type or whether both are affected. There have been a number of studies which have shown that HepaRG™ cells were able to predict drug-induced liver toxicity in humans due to their unique stable differentiated functions and ability to mirror *in vivo* effects. The toxicity of many compounds has been tested using rodent models but important species differences in metabolism and DME regulation result in conclusions which are not human-specific. For example, phenobarbital hepatotoxicity is much lower in humans than in rodents; and reasons for this are related to the gene expression profile over time, as shown in PHHs and HepaRG™ cells³⁶. Moreover, the expression profiles measured in HepaRG™ cells were comparable to PHHs, supporting the use of these cells for dose- and time-dependent studies. The use of HepaRG™ cells as a model for cytotoxicity assays has recently been investigated using a number of compounds. For example, two hepatotoxicants which require bioactivation are aflatoxin B₁ and acetaminophen and the IC₅₀s for these in HepaRG™ cells were similar to those measured in PHHs⁹.

Moreover, HepaRG™ cells were sensitive to the cytotoxic effect of these compounds at concentrations which were completely non-toxic to HepG2 cells^{13,20}. Aflatoxin B₁ cytotoxicity was time-dependent and specific to hepatocyte-like cells in the culture, such that depletion of ATP did not exceed 60% of the control levels (equivalent to the percentage of hepatocyte-like cells in the culture)²⁰. Thus HepaRG™ cells were able to differentiate between metabolism-dependent effects of compounds within the same culture. Studies on HepaRG™ cells treated with acetaminophen showed these cells demonstrated the key features of toxicity of this compound, namely marked glutathione depletion, adduct formation and mitochondrial dysfunction. Moreover, the release of LDH over time reflected the *in vivo* release of aminotransferase activity seen in humans as a result of acetaminophen overdose³⁷. Other hepatotoxicants, such as iron citrate, have been shown to be accumulated in HepaRG™ cells and PHHs, which only occurs in cells which have hepatocyte-like phenotypes (such as the presence of CYP2E1 and 3A4, lipid metabolism (e.g. fatty acid binding protein 1) and stress responses (e.g. superoxide dismutase)). Unlike the HepaRG™ hepatocyte-like cells, the HepaRG™ biliary cells do not accumulate iron³⁸. In another study, the relative metabolism-dependent hepatotoxic effects of three different thiopurines were demonstrated to be similar in HepaRG™ cells and PHHs³⁹. The sensitivity of both PHHs and HepaRG™ cells was lower than in rat hepatocytes, demonstrating a significant species difference in toxicity and the importance of using a human cell type. However, the human cell line, HepG2, was insensitive to the thiopurine, azathioprine, probably due to the lack of GST activity (mainly GSTM1) which catalyses the first metabolic step in the bioactivation of this compound⁴⁰. HepaRG™ cells have been used to determine the anti-tumor cytotoxicity of doxorubicin, delivered as free compound or loaded into nanoparticles⁴¹. The nanoparticles bypass the chemoresistant effects of MDR1-mediated efflux of the compound and thus making the cells more sensitive to toxicity. HepaRG™ cultures have also been shown to be adaptable to higher throughout cytotoxicity screening⁴².

Until now, there have been no *in vitro* hepatic models to predict steatosis (the accumulation of triglycerides) caused by drugs. However, this process was demonstrated in HepaRG™ cells treated with a number of polyunsaturated fatty acids and derivatives, suggesting that these cells may represent the first model for steatosis^{43,44}. HepaRG™ cells were shown to accumulate neutral lipids after treatment with the steatogenic drugs, tetracycline and amiodarone⁴³. Lipid accumulation increased after repeated dosing with relatively low concentrations of these drugs for two weeks.

Moreover, there was a direct correlation between lipid accumulation and levels of mRNA encoding for lipid synthesis markers in response to steatogenic drugs⁴³. The *in vitro* effects shown in HepaRG™ cells were predictive of the *in vivo* microvesicular steatosis reported in patients chronically treated with amiodarone and tetracycline.

Another mechanism of liver damage is phospholipidosis (accumulation of phospholipids and formation of lamellar bodies) and this liver injury has also been shown to occur in HepaRG™ cells after treatment with amiodarone^{11,43}. The HepaRG™ cells were able to discriminate between amiodarone which causes steatosis (after 24 h) and phospholipidosis (evident after 2 weeks) and tetracycline, which causes steatosis only. Similar to steatosis, the accumulation of phospholipids was concomitant with an up-regulation of lipid synthesis genes (and genes encoding proteins involved in the formation of vesicles) and was cumulative when the cells were dosed chronically over 2 weeks with non-cytotoxic concentrations of this drug. Similar studies using HepG2 cells showed that none of the genes involved in the formation of droplets were altered after amiodarone treatment, suggesting that HepG2 cells are poor models for this liver injury⁴⁵. It has been suggested that the gene markers involved in the formation of lamellar bodies and in the synthesis of phospholipids could represent potential biomarkers of drug-induced phospholipidosis in humans¹¹.

HepaRG™ cells have been used to investigate the effect of "glitazones" and "glitazars" on gene expression profiles using a whole genome transcriptomic approach⁴⁶. The first PPAR γ agonists with the intended therapeutic target of being antidiabetics were ciglitazone and troglitazone which were dropped during clinical trials and from the market, respectively, due to severe liver failure and death. Subsequent follow-up compounds, rosiglitazone and pioglitazone, were much less hepatotoxic. The dual PPAR α and PPAR γ agonists, muraglitazar and tesaglitazar, were not hepatotoxic but they were cardiotoxic and nephrotoxic⁴⁷. Incubations of troglitazone, rosiglitazone, muraglitazar and tesaglitazar with PHHs and HepaRG™ cells showed that troglitazone was the most toxic, causing loss of ATP, a decrease in NTCP uptake activity and inducing apoptosis. All four drugs altered genes regulating lipid, carbohydrate, xenobiotic and cholesterol metabolism, as well as inflammation and immunity in PHHs and HepaRG™ cells. These studies provided a comprehensive assessment of the effects of these drugs on both PHHs and HepaRG™ cells and demonstrated the relative contribution of different gene pathways involved in hepatotoxicity.

Genotoxicity and carcinogenicity

There are a number of studies which have focused on the possibility of using HepaRG™ cells for genotoxicity assays, since they are able to proliferate and they express significantly higher levels of DMEs than rodent cell lines or HepG2 cells. These cells have been developed and successfully used in the micronucleus and Comet assays^{20,48,49,50}. HepaRG™ cells therefore represent a promising cell model for the micronucleus test and have been recommended to be tested for use in this assay by ECVAM and IWGT workshops^{51,52}.

Phenotypic changes to HepaRG™ cells leading to malignant cell clones may also provide information on potential carcinogens⁵³. Whole genome expression analysis of HepaRG™ cells has been used to identify a set of genes linked to genotoxic and non-genotoxic compounds⁴⁹. Using this technology, HepaRG™ cells may be used as an alternative to PHHs to discriminate between genotoxic and non-genotoxic compounds and to provide information on compound and time dependent effects on cell cycle and apoptosis signaling pathways.

Transporter studies

The presence and function of important efflux and uptake transporters, together with the formation of tight junctions and correct location of canaliculi in cultured HepaRG™ cells, makes these cells suitable for biliary secretion studies². These studies were designed to reflect *in vivo*-like uptake, metabolism and secretion of both the parent and metabolites. Although the expression and activity of sinusoidal drug transporters such as OCT1, OATPs and NTCP are lower in HepaRG™ cells than in PHHs, the levels expressed are far higher than other cell lines, such as HepG2 cells, in which these transporters are low or lacking⁵⁴.

Transporters and their inhibition are thought to play a key role in cholestasis, in which the liver retains compounds that would normally be excreted into the bile. There is a lack of a model that can accurately predict the cholestatic effects of drugs in humans. Chlorpromazine, a known cholestatic drug, inhibited bile acid efflux in HepaRG™ cells, resulting in the accumulation of tritiated-taurocholic acid and oxidative stress (Antherieu et al., unpublished). Cholestasis can also be caused by inflammatory mediators, via down-regulation of transporters, as described below. Therefore, drugs causing cholestasis indirectly via inflammation can potentially be detected using HepaRG™ cells.

Inflammation

The liver plays an important role in sepsis and subsequent multiple organ failure. Sepsis is accompanied by a release of lipopolysaccharide (LPS), which compromises hepatocellular functions by upregulating

Kupffer cell tumor necrosis factor (TNF α), interleukin-1 β , (IL-1 β), and IL-6 production via the Toll-like receptor-4 and its co-receptor CD14 (which are both expressed in HepaRG™ cells⁵⁵). HepaRG™ cells have been shown to be responsive to LPS and exhibit classical responses of PHHs to inflammatory stimulation⁵⁵. Moreover, HepaRG™ cells are able to differentiate between the β -receptor-specific modulations of inflammatory responses and provide important information on the most suitable catecholamine to use for treatment of sepsis.

Human-specific effects of IL-1 β on transporter expression, function and liver cholestasis have been studied in HepaRG™ cells and compared with those in PHHs⁵⁶. The changes in uptake and efflux transporter profiles in PHHs were largely mirrored by those in HepaRG™ cells. Notably, NTCP and BSEP mRNA and protein expression and uptake of taurocholate were markedly decreased in both cell types by IL-1 β , suggesting that this inflammatory mediator participates in both cholestasis and perturbation of hepatic detoxification pathways caused in humans.

Viral infections

There is a huge demand for better viral infection models. Specifically, more information on the processes involved in virus receptor binding, uptake, and membrane fusion is needed so that treatments to prevent virus entry into cells can be developed. Hepatoma cells, such as HepG2, have a limited use in viral infection research since they have impaired antiviral responses. Animal models are less relevant and, for ethical reasons, should be avoided. Until now, no hepatoma cell line has been shown to be susceptible to viral infection and support viral replication. Although Huh-7 cells are used for hepatitis C virus (HCV) infection and replication they are poorly differentiated and so do not reflect effects seen *in vivo*. By contrast, HepaRG™ cells retain functional and efficient signalling pathways involved in pro-inflammatory effects and are susceptible to viral infection^{5,57}. Markers of viral infections, such as IFN- β and a number of chemokines involved in inflammation and chemotaxis were strongly induced in HepaRG™ cells exposed to extracellular stimuli and the appearance of newly synthesized hepatitis B virus (HBV) transcripts was demonstrated^{6,58,59}. Modulation of the small surface protein of the HBV envelope has been shown to interfere with HBV infections of PHHs and HepaRG™ cells, thus giving vital information on the mechanisms of infection and the cell surface receptor^{60,61}. Studies investigating the infection of HepaRG™ cells with HCV have provided information on the mechanisms involved in the internalization and replication of this virus⁶². Unlike HBV, the permissiveness of HepaRG™ cells to HCV is much lower and not as long-lasting, because HepaRG™ cells

have a strong type I interferon-mediated anti-viral response⁴. These findings regarding the permissiveness of HepaRG™ cells to HBV and HCV infection support the use of HepaRG™ cells as an important model for viral hepatitis research.

Drug target identification

HepaRG™ cells have been used in investigations into lipid metabolism and its regulation by PPAR α in order to identify novel therapeutic targets for dyslipidemia⁶³. In these studies, expression of the apoA-IV gene was significantly increased by the compound KRP-101 via PPAR α upregulation in HepaRG™ cells. The therapeutic relevance was confirmed *in vivo*, in which dogs treated with KRP-101 had decreased serum triglyceride levels and increased serum apoA-IV levels.

3D cultures and bioreactors

HepaRG™ cells have been cultured in three-dimensional (3D) bioreactors in volumes ranging from 800 mL to 2 mL^{64,65}. These cultures were subsequently used for drug metabolism studies. The morphology of the cells inside the bioreactor resembled the *in vivo* arrangement of the cells and the expression of transporter proteins reflected appropriate cell polarization. CYP2C9, CYP1A1/2, and CYP3A4-mediated activities were similar 2 and 17 days after the cells were differentiated with DMSO, while CYP2B6 activities increased over this time. The bioreactor model was also able to demonstrate the CYP3A4 selective inhibition by ketoconazole (a decrease of 69%) and a 6-fold induction of CYP3A4 and CYP2B6 activities. Both these inhibitory and induction responses were similar to those measured *in vivo*; moreover, all the results from this study were carried out in one bioreactor, indicating that sequential experiments can be successfully conducted using the same dynamic cell model.

HepaRG™ cells have been tested for their application to bioartificial liver (BAL) models⁶⁶. A main function of a BAL model is to remove ammonia via the urea cycle and conjugation with amino acids, especially glutamine by glutamine synthetase. HepG2 cells have a limited use in BAL models because they have such low DME activities and lack a functional urea cycle⁶⁷. The low urea production measured in 2D monolayer cultures of HepaRG™ cells was improved by preconditioning the BAL with carbamoylglutamate, an analogue of N-acetylglutamate. This compound activates carbamoylphosphate synthetase, which is the rate-determining enzyme of the urea cycle under normal physiological conditions⁶⁸. HepaRG™ cells are a promising cell type for use in BAL models; however, in order to generate a fully effective BAL, the urea cycle in these cells would have to be improved further⁶⁸.

References

1. Coecke, S. et al. 2006. Metabolism: a bottleneck in *in vitro* toxicological test development. The report and recommendations of ECVAM workshop 54. Alternatives to laboratory animals: ATLA 34, 49-84.
2. Parent R. et al. Origin and characterization of a human bipotent liver progenitor cell line. Gastroenterology 126, 1147-56 (2004).
3. Cerec V. et al. Transdifferentiation of hepatocyte-like cells from the human hepatoma HepaRG™ cell line through bipotent progenitor. Hepatology 45, 957-67 (2007).
4. Marion MJ et al. The HepaRG™ cell line: biological properties and relevance as a tool for cell biology, drug metabolism, and virology studies. Methods Mol Biol 640, 261-272 (2010).
5. Parent R, Beretta L. Translational control plays a prominent role in the hepatocytic differentiation of HepaRG™ liver progenitor cells. Genome Biol 9, R19 (2008).
6. Gripon P et al. Infection of a human hepatoma cell line by hepatitis B virus. Proc Natl Acad Sci U S A 99(24), 15655-15660 (2002).
7. Rogue A et al. Interindividual variability in gene expression profiles in human hepatocytes and comparison with HepaRG™ cells. Drug Metab Dispos 40(1), 151-8 (2012).
8. Hart SN et al. A comparison of whole genome expression profiles of HepaRG™ cells and HepG2 cells to primary human hepatocytes and human liver tissues. Drug Metab Dispos (2010).
9. Aninat C et al. Expression of cytochromes P450, conjugating enzymes and nuclear receptors in human hepatoma HepaRG™ cells. Drug Metab Dispos 34 75-83 (2006).
10. Anthérieu S et al. Stable Expression, Activity, and Inducibility of Cytochromes P450 in Differentiated HepaRG™ Cells. Drug Metab Dispos 38, 516-525, (2010).
11. Anthérieu S et al. Optimization of the HepaRG™ cell model for drug metabolism and toxicity studies. Toxicol In Vitro May 27 (2012) [Epub ahead of print]
12. Kanebratt KP, Andersson TB. Evaluation of HepaRG™ cells as an *in vitro* model for human drug metabolism studies. Drug Metab Dispos 36, 1444-1452 (2008).
13. Guillouzo A et al. The human hepatoma HepaRG™ cells: a highly differentiated model for studies of liver metabolism and toxicity of xenobiotics. Chem Biol Interact 168, 66-73 (2007).
14. Laurent V et al. Highly efficient gene transfer into hepatocyte-like HepaRG™ cells: new means for drug metabolism and toxicity studies. Biotechnol J 5(3), 314-320 (2010).
15. Antoun J et al. CYP4F3B is induced by PGA1 in human liver cells: a regulation of the 20-HETE synthesis. J Lipid Res 49, 2135-2141 (2008).
16. Antoun J et al. CYP4A11 is repressed by retinoic acid in human liver cells. FEBS Lett 580, 3361-3367 (2006).
17. Bjorkhem, I. On the quantitative importance of omega-oxidation of fatty acids. J Lipid Res 19, 585-590 (1978).
18. Turpeinen M et al. Functional Expression, Inhibition, and Induction of CYP Enzymes in HepaRG™ cells. Toxicol In Vitro 23(4), 748-53 (2009).
19. Kanebratt KP, Andersson TB. HepaRG™ cells as an *in vitro* model for evaluation of cytochrome P450 induction in humans. Drug Metab Dispos 36, 137-145 (2008).
20. Jossé R. et al. Long-term functional stability of human HepaRG™ hepatocytes and use for chronic toxicity and genotoxicity studies. Drug Metab Dispos 36, 111 -1118 (2008).
21. Hariparsad N et al. Comparison of immortalized Fa2N-4 cells and human hepatocytes as *in vitro* models for cytochrome P450 induction. Drug Metab Dispos 36, 1046-1055 (2008).
22. Huang S-M, Stifano T. (2006). Guidance for industry: Drug interaction studies – study design, data analysis, and implications for dosing and labelling. Draft Guidance.
23. Hewitt NJ et al. Primary hepatocytes: current understanding of the regulation of metabolic enzymes and transporter proteins, and pharmaceutical practice for the use of hepatocytes in metabolism, enzyme induction, transporter, clearance, and hepatotoxicity studies. Drug Metab Rev 39(1), 159-234 (2007).
24. Lübberstedt M et al. HepaRG™ human hepatic cell line utility as a surrogate for primary human hepatocytes in drug metabolism assessment *in vitro*. J Pharmacol Toxicol Methods 63(1), 59-68 (2011).
25. Zanelli U et al. Comparison of cryopreserved HepaRG™ cells with cryopreserved human hepatocytes for prediction of clearance for 26 drugs. Drug Metab Dispos 40(1), 104-10 (2012).
26. Le Vee M et al. Functional expression of sinusoidal and canalicular hepatic drug transporters in the differentiated human hepatoma HepaRG™ cell line. Eur J Pharm Sci 28, 109-17 (2006).
27. LeCluyse E et al. Expression and regulation of cytochrome P450 enzymes in primary cultures of human hepatocytes. J Biochem Mol Toxicol 14, 177-188 (2000).
28. Grime K et al. The use of HepaRG™ and human hepatocyte data in predicting CYP induction drug-drug interactions via static equation and dynamic mechanistic modelling approaches. Curr Drug Metab 11(10), 870-85 (2010).
29. Templeton IE, Houston JB, Galetin A. Predictive utility of *in vitro* rifampin induction data generated in fresh and cryopreserved human hepatocytes, Fa2N-4, and HepaRG™ cells. Drug Metab Dispos 39(10), 1921-1299 (2011).
30. Kaneko A et al. *In vitro* model for the prediction of clinical CYP3A4 induction using HepaRG™ cells. Xenobiotica 39(11), 803-10 (2009).
31. McGinness DF et al. Evaluation of multiple *in vitro* systems for assessment of CYP3A4 Induction in Drug Discovery : Human hepatocytes, PXR reporter gene, Fa2N-4, and HepaRG™ cells. Drug Metab Dispos 37, 1259-1268 (2009).
32. Andersson TB. The application of HepRG™ cells in evaluation of cytochrome P450 induction properties of drug compounds. Methods Mol Biol 640, 375-387 (2010).
33. ECVAM Technical report on the Status of Alternative Methods for Cosmetics Testing, 2008-2009; <http://ecvam.jrc.ec.europa.eu/>
34. Nishimura M et al. Comparison of inducibility of sulfotransferase and UDP-glucuronosyltransferase mRNAs by prototypical microsomal enzyme inducers in primary cultures of human and cynomolgus monkey hepatocytes. Drug Metab Pharmacokinet 23(1),45-53 (2008).
35. Jigorel E et al. (2006) Differential regulation of sinusoidal and canalicular hepatic drug transporter expression by xenobiotics activating drug-sensing receptors in primary human hepatocytes. Drug Metab Dispos 34:1756-1763.

36. Guguen-Guillouzo C et al. Maintenance and reversibility of active albumin secretion by adult rat hepatocytes co-cultured with another liver epithelial cell type. *Exp Cell Res* 143(1), 47-54 (1983).
37. Lambert et al, Guillouzo A. Reproducible chemical-induced changes in gene expression profiles in human hepatoma HepaRG™ cells under various experimental conditions. *Toxicol In Vitro* 23, 466-475 (2009).
38. McGill MR et al.. HepaRG™ cells: a human model to study mechanisms of acetaminophen hepatotoxicity. *Hepatology Mar*;53(3), 974-978 (2011).
39. Troadec MB et al. Hepatocyte iron loading capacity is associated with differentiation and repression of motility in the HepaRG™ cell line. *Genomics* 87, 93-103 (2006).
40. Petit E et al. Differential toxic effects of azathioprine, 6-mercaptopurine and 6-thioguanine on human hepatocytes. *Toxicol In Vitro* 22, 632-642 (2008).
41. Lee AU, Farrell GC. Mechanism of azathioprine-induced injury to hepatocytes: roles of glutathione depletion and mitochondrial injury. *J Hepatol* 35(6), 756-64 (2001).
42. Barraud L et al. Increase of doxorubicin sensitivity by doxorubicin-loading into nanoparticles for hepatocellular carcinoma cells *in vitro* and *in vivo*. *J Hepatol* 42, 736-73 (2005).
43. Pernelle K et al. Automated detection of hepatotoxic compounds in human hepatocytes using HepaRG™ cells and image-based analysis of mitochondrial dysfunction with JC-1 dye. *Toxicol Appl Pharmacol* 254(3), 256-266 (2011).
44. Anthérieu S et al. Induction of vesicular steatosis by amiodarone and tetracycline is associated with up-regulation of lipogenic genes in HepaRG™ cells. *Hepatology* 53(6), 1895-1890 (2005).
45. Madec S et al. CYP4F3B expression is associated with differentiation of HepaRG™ human hepatocytes and unaffected by fatty acid overload. *Drug Metab Dispos* 39(10), 1987-1896 (2011).
46. Sawada H et al. A toxicogenomic approach to drug induced phospholipidosis: analysis of its induction mechanism and establishment of a novel *in vitro* screening system. *Toxicol Sci* 83, 282-292 (2005).
47. Rogue A et al. Comparative gene expression profiles induced by PPARγ and PPARα/γ agonists in human hepatocytes. *PLoS One* 6(4), e18816 (2011).
48. Fievet C et al.. PPARalpha and PPARgamma dual agonists for the treatment of type 2 diabetes and the metabolic syndrome. *Curr Opin Pharmacol* 6, 606-614 (2006).
49. Bazin E et al. Genotoxicity of a freshwater cyanotoxin, cylindrospermopsin, in two human cell lines: Caco-2 and HepaRG. *Environ Mol Mutagen* 51(3), 251-259 (2010).
50. Jossé R et al. An adaptation of the human HepaRG™ cells to the *in vitro* micronucleus assay. *Mutagenesis* 27(3), 295-304 (2012).
51. Le Hegarat L et al. Assessment of the genotoxic potential of indirect chemical mutagens in HepaRG™ cells by the comet and the cytokinesis-block micronucleus assays. *Mutagenesis* 25(6), 555-560 (2010).
52. Kirkland D et al. How to reduce false positive results when undertaking *in vitro* genotoxicity testing and thus avoid unnecessary follow-up animal tests: Report of an ECVAM Workshop. *Mutat Res* 628, 31-55 (2007).
53. Pfuhrer S, et al. *In vitro* genotoxicity test approaches with better predictivity: summary of an IWGT workshop. *Mutat Res* 723(2), 101-107 (2011).
54. Dumont J et al. Preferential induction of the AhR gene battery in HepaRG™ cells after a single or repeated exposure to heterocyclic aromatic amines. *Toxicol Appl Pharmacol* 249(1), 91-100 (2010).
55. Cui Y et al. Detection of the human organic anion transporters SLC21A6 (OATP2) and SLC21A8 (OATP8) in liver and hepatocellular carcinoma. *Lab Invest* 83(4), 527-538 (2003).
56. Aninat C et al. Catecholamines induce an inflammatory response in human hepatocytes. *Crit Care Med* 36(3), 848-854 (2008).
57. Le Vee M et al. Down-regulation of organic anion transporter expression in human hepatocytes exposed to the proinflammatory cytokine interleukin 1beta. *Drug Metab Dispos* 36, 217-222 (2008).
58. Hantz, O et al. (2009) Persistence of the hepatitis B virus covalently closed circular DNA in HepaRG™ human hepatocyte-like cells. *J. Gen. Virol.* 90, 127-135.
59. Maire M et al, Petit MA. Characterization of the double-stranded RNA responses in human liver progenitor cells. *Biochem Biophys Res Commun* 368, 556-562 (2008).
60. Gripon P et al. Efficient inhibition of hepatitis B virus infection by acylated peptides derived from the large viral surface protein. *J Virol* 79, 1613-1622 (2005).
61. Jaoudé GA, Sureau C. Role of the antigenic loop of the hepatitis B virus envelope proteins in infectivity of hepatitis delta virus. *J Virol* 79 10460-10466 (2005).
62. Ndongo-Thiam N et al.. Long-term propagation of serum hepatitis C virus (HCV) with production of enveloped HCV particles in human HepaRG™ hepatocytes. *Hepatology* 54(2), 406-417 (2011).
63. Nagasawa M et al. Highly sensitive upregulation of apolipoprotein A-IV by peroxisome proliferator-activated receptor alpha (PPARalpha) agonist in human hepatoma cells. *Biochem Pharmacol* 74, 1738-1746 (2007).
64. Darnell M et al.. Cytochrome P450-dependent metabolism in HepaRG™ cells cultured in a dynamic three-dimensional bioreactor. *Drug Metab Dispos* 39(7), 1131-8 (2011).
65. Gerlach JC et al. Bioreactor for a larger scale hepatocyte *in vitro* perfusion. *Transplantation* 58, 984-988 (1994).
66. Hoekstra R et al. The HepaRG™ cell line is suitable for bioartificial liver application. *Int J Biochem Cell Biol* 43(10), 1483-1489 (2011).
67. Mavri-Damelin D et al. Ornithine transcarbamylase and arginase I deficiency are responsible for diminished urea cycle function in the human hepatoblastoma cell line HepG2. *Int J Biochem Cell Biol* 39(3), 555-64. (2007).
68. Kim JS et al. Mitochondrial permeability transition: a common pathway to necrosis and apoptosis. *Biochem. Biophys. Res. Commun* 304 (3), 463-470 (2003).
69. Wilkening S et al. Comparison of primary human hepatocytes and hepatoma cell line Hepg2 with regard to their biotransformation properties. *Drug Metab Dispos* 31(8), 1035-42 (2003).
70. Gerets HH et al. Characterization of primary human hepatocytes, HepG2 cells, and HepaRG™ cells at the mRNA level and CYP activity in response to inducers and their predictivity for the detection of human hepatotoxins. *Cell Biol Toxicol* 28(2), 69-87 (2012).

Related Products

Product	Size	Catalogue No.
Cryopreserved HepaRG™ Cells	1 mL in vial with $\geq 8 \times 10^6$ viable cells	MMHPR116
HepaRG™ Thawing/Plating Medium Supplement	12.5 mL	MMADD671
HepaRG™ Tox Medium Supplement	12.5 mL	MMADD631
HepaRG™ Induction Medium Supplement	2.6 mL	MMADD641
HepaRG™ Serum Free Induction Medium Supplement	0.6 mL	MMADD651
HepaRG™ Culture Medium Supplement	14.0 mL	MMADD621

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe,
please call: +44 (0) 115 943 0840

Or visit: www.merckmillipore.com/offices

For Technical Service visit:

www.merckmillipore.com/techservice



www.merckmillipore.com/heparg

Merck Millipore and the M logo are trademarks of Merck KGaA, Darmstadt, Germany.
HepaRG is a trademark of Biopredic International. Lit No. AN4455ENEU BS GEN-12-07120
10/2012 © 2012 Merck KGaA, Darmstadt, Germany. All rights reserved.

Get Connected!

Join Merck Millipore Bioscience on your favorite social media outlet for the latest updates, news, products, innovations, and contests!



[facebook.com/MerckMilliporeBioscience](https://www.facebook.com/MerckMilliporeBioscience)



twitter.com/Merck4Bio