



New Workflow for Lentivirus Purification, Concentration, and Immunodetection

Janet Smith¹, Sara Gutierrez¹, Charles Neville¹, Mikhail Kozlov, PhD², and Kathleen Ongena, PhD¹

Millipore Corporation, Billerica, MA USA

Presented at the 2008 Gene Therapy Conference

May, 2008-Boston, MA USA



www.millipore.com

Millipore is a registered trademark of Millipore Corporation. Amicon, Steriflip, Ultrafree, and Immobilon are registered trademarks of Millipore Corporation. Fast-Trap and SNAP id are trademarks of Millipore Corporation. MagicMark and Mark 12 are trademarks of Invitrogen Corporation. NuPage is a registered trademark of Invitrogen Corporation. Quick Titer is a trademark of Cell Biolabs. Sypro is a trademark of Molecular Probes, Inc. Transblot is a trademark of BioRad Laboratories Inc. Tween is a trademark of ICI Americas Inc. Lit. No. PS2076ENUS Printed in U.S.A. 2/08
©Millipore Corporation, Billerica, MA 01821 USA All rights reserved.

ADVANCING LIFE SCIENCE TOGETHER™
Research. Development. Production.

upstate | CHEMICON | *Linco*
THE EXPERTISE OF UPSTATE®, CHEMICON® AND LINCO®
IS NOW A PART OF MILLIPORE

New Workflow for Lentivirus Purification, Concentration, and Immunodetection

Janet Smith¹, Sara Gutierrez¹, Charles Neville¹, Mikhail Kozlov, PhD² and Kathleen Ongena, PhD¹

¹Millipore Corporation, BioScience Division, Danvers, MA, United States, 01923 and ²Millipore Corporation, Bedford, MA, United States, 01730



Abstract

Lentivirus is a negatively charged, enveloped, single stranded RNA virus from the *Retroviridae* family that is often used as a vector to transport genetic material into cells. These viral vectors can be used for genetic modification, RNAi, gene therapy, and vaccine production. Before viral preparations/propagations can be used for any of these applications, researchers need to purify their virus sample.

Traditionally, time-consuming density gradient centrifugation separation and/or chromatographic techniques have been used. The lab-scale chromatographic devices are commonly syringe or column-based. To process the virus sample, these devices require hand pressure or gravity; this may lead to messy and to potentially unsafe handling conditions during assembly and disassembly.

A new workflow has been developed to clarify, purify, and concentrate/buffer exchange a crude lentivirus sample. For added safety and improved handling, the clarification and purification steps are performed in a closed vacuum-based device. This purification produces high recovery of virus particles in about one hour with similar or improved results as compared to traditional methods. Purity was visualized by gel electrophoresis and confirmed by western blotting using an innovative vacuum-based immunodetection system that allows detection of the protein of interest in less than thirty minutes.


Here we show the results of the purification and the immunodetection of a Lentivirus-VSVG pseudotype that encodes green fluorescent protein (GFP). We demonstrate the viral titer, the percent recovery of infectious particles, and the purity of the virus sample.

Lentivirus Fast-Trap™ Virus Purification & Concentration Kit combined with SNAP i.d.™ Immunodetection enables results in less than three hours.

New Workflow: Virus Purification & Concentration




Result: Purified Virus in 1 to 2.5 hrs.



- Genetic modification
- RNAi
- Gene therapy
- Vaccine production

Immunodetection Workflow With the SNAP i.d. Protein Detection System



Protein Detection Method:

- Step 4: Electrophoresis-Samples denatured & reduced before SDS-PAGE
- Step 5: Transfer to Immobilon®-P membrane (Millipore #IPVH07850)
- Step 6: Immunodetection with SNAP i.d. Protein Detection System (Millipore #WBABVDBASE)
- Step 6A: Apply blocking reagent containing Tween® 20 surfactant and vacuum filter
- Step 6B: Add antibodies (primary, secondary) and wash between antibody additions
- Step 7: Add Immobilon® Western Chemiluminescent HRP substrate (Millipore #WBKLS0500) for chemiluminescence detection of protein of interest and visualize immunoreactive proteins

Performance of the Fast-Trap Lentivirus Purification & Concentration Kit and SNAP i.d. Protein Detection Workflow

Methods

Fast-Trap Purification

Virus: Crude lentivirus (University of North Carolina [UNC] Vector Core) was benzonase treated, clarified, purified, and concentrated using the new Fast-Trap kit workflow

SDS-PAGE

Samples: Virus samples normalized for protein concentration, denatured and reduced NuPAGE® 4 – 12% Bis-Tris gel (Invitrogen) 200 V, 35 minutes

Gel Markers: M = Mark12™ unstained standard or MagicMark™ XP Western standard (Invitrogen)

Stain: SYPRO® Ruby stain (Invitrogen) overnight, destain, and visualize, or membrane transfer for immunodetection

Membrane Transfer

Transfer: Semi-dry transfer method (BioRad Trans-Blot® SD Semi-dry Transfer cell) to Immobilon-P membrane at 10 V, 35 minutes with Tris-Glycine buffer

SNAP i.d. Immunodetection

Blocking: Add blocking buffer (0.5% non-fat dry milk in Tris Buffered Saline with 0.1% Tween-20 surfactant, TBST) and vacuum



1° Antibody: Probe with Mouse anti-HIV₁ p24 (Millipore #MAB8790) diluted 1:13,000 in blocking buffer, and incubate for 10 minutes

Wash: Wash with TBST buffer using constant vacuum

2° Antibody: Add Goat anti-Mouse IgG, HRP-conjugate (Millipore #AP124P) diluted 1:10,000 in blocking buffer, and incubate for 10 minutes

Wash: Wash with TBST buffer using constant vacuum

Detect: Incubate 5 minutes with Immobilon Western Chemiluminescent HRP substrate, and visualize by exposing to x-ray film



Ratio of Infectious to Non-Infectious Viral Particles

	Unpurified Lentivirus	Purified Lentivirus
Average Concentration by ELISA (vp/mL x 10 ¹⁰)	3.6	22.5
Infective particle Input (Ivp x 10 ⁶)	2.8	12.9
Ratio (Ivp/vp x 10 ⁴)	1.3	1.7

Table 1: The ratio for unpurified and purified lentivirus (using the Fast-Trap Purification and Concentration Kit) was determined by Enzyme-Linked ImmunoSorbent Assay (ELISA) Quick Titer™ Lentivirus Quantitation kit (Cell Biolabs). This ratio is maintained when using Fast-Trap Lentivirus Purification and Concentration kit.


Lentivirus Purification: Comparison of the Fast-Trap Kit to Other Formats/Methods (Syringe, Column, Precipitation, & Traditional Sucrose Gradient)

Processing Time & Virus Recovery


Purification Format	Key	Processing Volume (ml)	Processing Time (min.)	Recovery of Infectious Viral Particles, Ivp (%)
Vacuum, closed system	Fast-Trap Lentivirus	9	6	57%
Gravity syringe	S1	9	11	9%
Gravity syringe	S2	9	10	21%
Gravity syringe	S3	9	11	29%
Column	Col	9	300	38%
Precipitation/Centrifugation	Ppt	10	12 Hr. + 35	57%

Table 2: Crude lentivirus (Input 1.4×10^6 Ivp) was purified with the Fast-Trap Virus Purification and Concentration Kit, three types of gravity syringe chromatographic columns (S1, S2 & S3), one column (Col), and one precipitation (Ppt) method (following manufacturer's instructions). The table compares overall performance and processing time (for bind/wash/elute, or precipitation/centrifugation steps only). The Fast-Trap purification method has the highest recovery and the shortest processing time. The closed, vacuum-based Fast-Trap purification device was easy to handle without the messy assembly/disassembly process or flow regulation that some formats require.

Sample Purity & Capsid Protein P24 Detection



Fast-Trap Virus Purification and Concentration Kit versus Traditional Sucrose Gradient Method



Conclusions

- Results can be generated in less than three hours using the Fast-Trap Virus Purification and Concentration kit followed by SNAP i.d. immunodetection
- Fast-Trap Virus Purification and Concentration Kit:
 - Enables efficient recovery of high titer lentivirus
 - Maintains ratio of non-infective to infective viral particles
 - Outperforms other chromatographic and precipitation formats for processing time, percent recovery of infectious particles, purity, and handling
 - Provides an alternative to precipitation, column or chromatography-based, and traditional gradient purification methods
- Vacuum-based SNAP i.d. Protein Detection System:
 - Shortens immunodetection process to ~ 30 minutes
 - Offers a faster alternative to traditional methods without consumption of additional reagents or loss of sensitivity

Acknowledgements

➢ UNC Vector Core & Dr. R. Jude Samulski for providing lentivirus samples