

Technical Data Sheet

Dermatophytes Selective Agar (DTM) acc. to TAPLIN

Ordering number: 1.10896.0500

Culture medium proposed by TAPLIN et al. (1969, 1970) for the isolation and rapid differentiation of dermatophytes from specimens including those infected with other microorganisms.

The comparative studies of MERTZ et al. (1970) demonstrated that the selectivity of Dermatophytes selective agar (DTM) is superior to that of other media used for cultivating fungi.

According to ALLEN et al. (1970), this medium offers the advantage that the dermatophytes grow rapidly and produce an unmistakable color change.

IVD in vitro diagnosticum - For professional use only

Mode of Action

This culture medium contains the pH indicator phenol red and the selective inhibitors cycloheximide, gentamicin and chlorotetracycline, which partly suppress the growth of bacteria, yeasts and molds. When grown on DTM, most dermatophytes produce basic metabolites, which bring about an alkalization of the acidic culture medium, causing the phenol red to change its color from yellow to red. This color change may, however, occasionally be caused by other microorganisms, too. Many molds produce acidic metabolites, which do not change the color of the culture medium. According to the authors, it is thus possible to differentiate rapidly between dermatophytes and other fungi with a high degree of accuracy (approx. 97 %).

Typical Composition

Peptone from soymeal	10.0 g/l
D(+)-glucose	10.0 g/l
cycloheximide	0.5 g/l
gentamicin sulfate	0.1 g/l
chlorotetracycline	0.1 g/l
phenol red	0.2 g/l
agar-agar	17.0 g/l

Preparation and Storage, Cat. No. 1.10896. Dermatophytes Selective Agar (DTM) acc. to TAPLIN (500 g).

Preparation

Suspend 38 g/litre, autoclave under mild conditions (10 min at 121 °C), pour plates or prepare slant tubes.

pH: 5.5 ±0.2 at 25 °C.

The plates are clear and yellow-orange.

Specimen


e.g. Nails, hair, skin.

Clinical specimen collection, handling and processing, see general instructions of use.

Experimental Procedure and Evaluation

Inoculate the surface of the culture medium with specimens obtained by appropriate methods.

Incubation: 7 days, possibly up to 3 weeks at approx. 28 °C aerobically.



Candida albicans ATCC 10231



Geotrichum candidum DSMZ 1240

Storage

Usable up to the expiry date when stored dry and tightly closed at +15 to +25 °C. Protect from light.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed at +15 to +25 °C.



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Merck Millipore, the M logo and Merckoplate® are registered trademarks of Merck KGaA, Darmstadt, Germany. ATCC is a registered trademark of ATCC, Manassas, VA, USA. Lit. No. TN1840EN00

Quality Control

Control Strains	ATCC #	Incubation	Expected Results
Trichophyton mentagrophytes	18748	7 d at 28 °C	Growth poor / good; Colour change to red
<i>Trichophyton rubrum</i>	28188	7 d at 28 °C	Growth poor / good; Colour change to red
<i>Microsporum gallinae</i>	12108	7 d at 28 °C	Growth poor / good; Colour change to red
<i>Microsporum canis</i>	36299	7 d at 28 °C	Growth poor / good; Colour change to red
<i>Geotrichum candidum</i>	DSMZ 1240	7 d at 28 °C	Growth fair / good; Colour change to red ±
<i>Candida albicans</i>	10231	7 d at 28 °C	Growth good / very good; Colour change to red
<i>Aspergillus brasiliensis</i> formerly <i>Asp. niger</i>	16404	7 d at 28 °C	Growth none / poor
<i>Penicillium commune</i>	10428	7 d at 28 °C	Growth none / poor
<i>Saccharomyces cerevisiae</i>	9763	7 d at 28 °C	Growth none
<i>Bacillus cereus</i>	11778	7 d at 28 °C	Growth none
<i>Escherichia coli</i>	25922	7 d at 28 °C	Growth none
<i>Staphylococcus aureus</i>	25923	7 d at 28 °C	Growth none

Please refer to the actual batch related Certificate of Analysis.

Literature

ALLEN, A.M., DREWRY, R.A., a. WEAVER, R.E.: Evaluation of a new color indicator media for diagnosis of dermatophytosis. - Arch. Derm., 102; 68- 70 (1970).

MERTZ, W.G., BERGER, C.L., a. SILVA-HUTNER, M.: Media with pH-indicator for the isolation of dermatophytes. - Arch. Derm., 99; 203-209 (1969).

TAPLIN, D., ALLEN, A.M., a. MERTZ, P.M.: Experience with a new indicator medium (DTM) for the isolation of dermatophyte fungi, in "Proceedings of the International Symposium of Mycoses", scientific publication 205. Washington, D.C. Pan American Health Organization, 55-58 (1970).

Ordering Information

Product	Cat. No.	Pack size
Dermatophytes Selective Agar (DTM) acc. to TAPLIN	1.10896.0500	500 g

Merck KGaA, 64271 Darmstadt, Germany
 Fax: +49 (0) 61 51 / 72-60 80
 mibio@merckgroup.com
www.merckmillipore.com/biomonitoring

Find contact information for your country at:
www.merckmillipore.com/offices
 For Technical Service, please visit:
www.merckmillipore.com/techservice



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Merck Millipore, the M logo and Merckoplate® are registered trademarks of Merck KGaA, Darmstadt, Germany. ATCC is a registered trademark of ATCC, Manassas, VA, USA. Lit. No. TN1840EN00