1.16982.0001

Reflectoquant® Iron Test for RQflex® 20

Fe

1. Method

In acidic solution iron(II) ions react with a triazine derivative to form a red-violet complex that is determined reflectometrically.

2. Measuring range and number of determinations

Measuring range	Number of determinations
0.5 - 20.0 mg/l Fe ²⁺	50

3. Applications

This test measures only iron(II) ions. Samples must be reduced by ascorbic acid before iron(III) ions and total iron (sum of iron(II) and iron(III)) can be measured (see section 7).

Sample material:

Surface water
Seawater
Industrial water
Wastewater
Mineral water
Food and preserves (e.g. juices, wine, liquids for canned preserves)

4. Influence of foreign substances

This was checked individually in solutions with 10 and 0 mg/l Fe²⁺. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %			
Ag ⁺	100	Mn ²⁺	1000
Al ³⁺	1000	NH ₄ ⁺	1000
BO ₃ ³⁻	1000	Ni ²⁺	100
Ca ²⁺	1000	NO ₂ ⁻	100
Cd ²⁺	1000	NO ₃ ⁻	1000
Citrate	1000	Oxalate	1000
CN ⁻	10	Pb ²⁺	1000
Cu ²⁺	10	PO ₄ ³⁻	1000
Fe ³⁺	1000	S ₂ O ₃ ²⁻	1000
K ⁺	1000	Tartrate	1000
Mg ²⁺	1000	Zn ²⁺	1000
		Anionic surfactants ¹⁾	1000
		Cationic surfactants ²⁾	1000
		Nonionic surfactants ³⁾	1000
		NaCl	20%
		Na ₂ SO ₄	20%

¹⁾ tested with Na-dodecyl sulfate

²⁾ tested with N-cetylpyridinium chloride

³⁾ tested with polyvinylpyrrolidone

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test strips are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

Tube containing 50 test strips
1 bar-code strip

Other reagents:

MQuant® Iron Test, Cat. No. 1.10004, measuring range 3 - 500 mg/l Fe²⁺
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sodium acetate anhydrous for analysis EMSURE®, Cat. No. 1.06268
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
L(+)-Ascorbic acid for analysis EMSURE®, Cat. No. 1.00468
Ammonium iron(II) sulfate hexahydrate for analysis EMSURE®, Cat. No. 1.03792

6. Preparation

- Check the iron content with the MQuant® Iron Test.
Samples containing more than 20.0 mg/l Fe²⁺ must be diluted with distilled water.
- The pH must be within the range 1 - 4.**
If necessary, buffer the sample with sodium acetate or, respectively, adjust the pH with sulfuric acid.

7. Procedure

Important note: This test can only be used in conjunction with the reflectometer **RQflex® 20!** Observe the manual for the reflectometer RQflex® 20.

The following applies to the Iron Test:

Measurement procedure A

Stored reaction time: 15 sec

Press the START button of the reflectometer and - **this is imperative - at the same time** immerse **both reaction zones** of the test strip in the pretreated sample (**5 - 40 °C**) for **2 sec**.

Carefully allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel.

Immediately insert the strip all the way into the strip adapter with the reaction zones facing the display.

After the end of the reaction time, read off the result from the display in mg/l Fe²⁺. The result is automatically stored.

Determination of total iron:

To 10 ml of sample add approx. 10 mg of ascorbic acid and shake. After approx. 1 min perform the determination as described in section 7.

Calculation of the iron(III) content:

$$\text{mg/l Fe}^{3+} = \text{mg/l total iron} - \text{mg/l Fe}^{2+}$$

Notes on the measurement:

- If the measurement value exceeds the measuring range, repeat the measurement using **fresh**, diluted samples until a value of less than 20.0 mg/l Fe²⁺ is obtained.

Concerning the result of the analysis, the dilution (see also section 6) must be taken into account:

$$\text{Result of analysis} = \text{measurement value} \times \text{dilution factor}$$

- If the test strip is inserted into the adapter after the reaction time has expired, renewed depression of the START button may produce a false result.

8. Method control

To check test strips, measurement device, and handling (recommended before each measurement series):

Dissolve 0.7 g of ammonium iron(II) sulfate hexahydrate in distilled water, make up to 100 ml with distilled water, and mix. Fe²⁺ content: 1000 mg/l.

Dilute this standard solution with distilled water to 10.0 mg/l Fe²⁺ and analyze as described in section 7.

Additional notes see under

www.sigmaaldrich.com/qa-test-kits.

9. Notes

- Reclose the tube containing the test strips immediately after use.**
- At the end of each workday, cleanse the strip adapter thoroughly with distilled water or ethanol.

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