

sealed with an



The new solvent bottle
with advanced sealing
technology

Meet the most innovative packaging system for solvents: **the new M-Bottle**, featuring advanced sealing technology that eliminates leaks, assuring you that your solvent handling is safer and easier than ever before.

SUPERIOR SAFETY


- Advanced sealing technology effectively eliminates leaks
- Additional headspace ensures greater safety with solvents that have higher vapor pressure at ambient temperatures
- Thicker glass reduces risk of breakage
- Redesigned 38mm bottle neck

EASY HANDLING

- New cap with larger ribs allows easier removal and closure
- Redesigned pouring ring facilitates handling and dispensing

GREATER CONVENIENCE

- Markings on bottle shoulder indicate a nominal volume of 4L
- Designed to fit in standard storage spaces and bottle carriers



The Life Science business of Merck
operates as MilliporeSigma in
the U.S. and Canada.

Sigma-Aldrich®

Lab & Production Materials

Supelco®

Analytical Products

The new M-Bottle for solvents



New neck design
captures tamper-
evident ring

Innovative sealing technology
effectively prevents leaks

**Additional finger
clearance**
easier, safer pouring

Additional headspace
greater safety for
solvents with higher
vapor pressure at
ambient temperatures

Optimized shape
fits in all standard
bottle carriers and
storage spaces



To place an order or receive technical assistance

Order/Customer Service: SigmaAldrich.com/order

Technical Service: SigmaAldrich.com/techservice

Safety-related Information: SigmaAldrich.com/safetycenter

SigmaAldrich.com/m-bottle

We have built a unique collection of life science brands
with unrivalled experience in supporting you scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

