Merck

Pharma & Biopharma Raw Material Solutions

Liquid Solutions MyWay - Asia

Our portfolio of sterile filtered liquids essential for your processes – in stock and ready to ship

Keep your process on the right track with Liquid Solutions MyWay – a portfolio of high-quality sterile filtered liquids in pre-filled single-use bags that are in stock and ready-to-ship. It's that simple: You order, we ship.

The predefined portfolio includes commonly used solutions for your bioprocessing needs. Not only are they in stock, but they have a 12 month shelf life from the date of production. We are adding longer stability studies which will in the future give us a shelf life of 24 months.



Liquid Solutions MyWay: In stock and ready to ship portfolio of high quality sterile filtered liquids

Fast	Convenient	No Compromise on Quality
Shorter lead times compared to Bags come with stand	Order directly from our websiteBags come with standard Colder	Ready-to-use bags to mitigate risk and reduce footprint
	MPC/MPX connectors	• 12 month shelf life from the date of production (target: 24 months)
		EMPROVE® raw materials
		Produced in an ISO 9001:2015 certified facility



Robust Packaging Ensures Safe Delivery

Liquid Solutions MyWay single-use bags are readyto-use and available in sizes from 1 L to 190 L/200 L. All sizes are dual sourced to ensure better availability and supply chain robustness. To safeguard delivery, all shipping containers are ship-tested; hazardous products are packaged in stainless steel while nonhazardous solutions arrive in PET drums.

Exacting Standards Ensure the Highest Quality

All components meet Ph. Eur. requirements, with most of the raw materials also meeting USP. Certificates of analysis (CoA) and safety data sheets (SDS) for the finished goods are available and our manufacturing sites produce in an ISO 9001:2015 certified facility.

Ordering is Fast and Easy

For fast and easy online ordering, visit SigmaAldrich.com/liquid-solutions-my-way

Liquid Solutions MyWay Product Portfolio

Art. No.	Article	Bag Sizes	Typical Application
55120C	Phosphate Buffered Saline (PBS)	1 L, 20 L, 200 L	Buffer
55220C	Sodium chloride solution 1 M	1 L, 20 L, 200 L	Universal
58400C	Tris solution 1 M	1 L, 20 L, 200 L	Buffer Stock Solution
59910C	Packaged Water for Injection	50 L, 200 L	Universal
58300C	Sodium hydroxide solution 0.1 M	1 L, 20 L, 190 L	CIP
58310C	Sodium hydroxide solution 0.5 M	1 L, 20 L, 190 L	CIP
58320C	Sodium hydroxide solution 1 M	1 L, 20 L, 190 L	CIP
55256C	Acetic Acid concentrate 5 M (100 mM Acetic Acid, pH 3.2, at 50x dilution)	1 L, 20 L, 190 L	Universal
55258C	Citric Acid concentrate 1.5 M (50 mM Citric Acid, pH 3.2, at 30x dilution)	1 L, 20 L, 190 L	Elution, pH adjustment
58410C	Tris concentrate 2.5 M (50 mM Tris, pH 7.5, at 50x dilution)	1 L, 20 L, 200 L	Buffer Stock Solution, pH adjustment

To place an order or receive technical assistance:

Order/Customer Service: SigmaAldrich.com/order
Technical Service: SigmaAldrich.com/techservice
Safety-related Information: SigmaAldrich.com/safetycenter

Merck KGaA Frankfurter Strasse 250 64293 Darmstadt, Germany

SigmaAldrich.com



© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, SAFC, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.