

## 67652 *Acinetobacter* Millichrome™ plus Supplement

Growth promotion supplement for *Acinetobacter*.

### Composition

(1 vial is sufficient for 5 L of final media)

Growth and regulator factors

### Directions

Add 4.0 mL of the liquid *Acinetobacter* Millichrome™ plus Supplement into 1L *Acinetobacter* Chromocult® plus Agar base.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

