



QUALITY CONTROL RANGES

HUMAN AMYLIN ELISA CATALOG # EZHA-52K

STANDARD, LYOPHILIZED CATALOG #E8051-K

STANDARD	AMYLIN VALUE
Lot Number. 7905-20K	149 pM

The exact value of lyophilized standard is shown above.

Reconstitute with **1.0 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6051-K

Reconstitute each vial with **0.25 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	EA-120	10– 21	pM
QC 2	EA-220	46 – 96	pM

E6051-K: EA-120.220 (25 March 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.