

Quality Control Ranges

MILLIPLEX® Human Circadian/Stress Magnetic Bead Panel

Catalog Number HNCSMAG-35K

QC Catalog Number/ Lot Number HNCS-6035 Lot no 107.207

Analyte	QC Level	Lot #	Expected Range	Units
Cortisol	Control I	HNCS-107	5767 - 11977	pg/mL
	Control II	HNCS-207	23405 - 48610	pg/mL
Melatonin	Control I	HNCS-107	50 - 103	pg/mL
	Control II	HNCS-207	175 - 363	pg/mL

The life science business of Merck KGaA, Darmstadt, Germany operates as

MilliporeSigma in the U.S. and Canada.

The vibrant M, Milliplex, Millipore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Template 20300637, Ver 4.0

Cat#: HNCS-6035, Lot#: 107.207 (Rev 03 Dec 2025)

