



**QUALITY CONTROL RANGES**  
**MILLIPLEX<sup>®</sup> MAP Canine Pituitary Magnetic Panel**  
**Catalog # CANPIT-96K**  
**Control Catalog # CANP-6096 Lot# CANP-102 and**  
**CANP-7096 Lot# CANP-202**

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
ACTH	CANP-6096	CANP-102	32 - 66	pg/mL
	CANP-7096	CANP-202	248 - 515	pg/mL
FSH	CANP-6096	CANP-102	33 – 68	ng/mL
	CANP-7096	CANP-202	260 - 540	ng/mL
GH	CANP-6096	CANP-102	0.16 – 0.34	ng/mL
	CANP-7096	CANP-202	1.3 – 2.7	ng/mL
TSH	CANP-6096	CANP-102	51 – 107	pg/mL
	CANP-7096	CANP-202	400 – 831	pg/mL
LH	CANP-6096	CANP-102	1.2 – 2.4	ng/mL
	CANP-7096	CANP-202	10 - 20	ng/mL
Prolactin	CANP-6096	CANP-102	17 – 35	ng/mL
	CANP-7096	CANP-202	138 - 286	ng/mL

**Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.**

**Quality Control values in culture media are not tested.**

CANP-6096 Lot# CANP-102 and CANP-7096 Lot# CANP-202 (30 September 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

[emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.