MILLIPORE

CHEMICON now part of Millipore

RABBIT ANTI- CYP450 2E1 POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB10080	QUANTITY:	50 μL	
LOT NUMBER:		HOST/ISOTYPE:	Rabbit	
IMMUNOGEN:	Recombinant CYP450 2E1			
APPLICATIONS:	Western Blot Analysis: Wester at a dilution of 1:5,000.	n Blot Analysis of AB10080	1	2
	Lane 2 Poole	mbinant CYP450 2E1 ed Human Liver Microsomes tes expected band at 57 kDa	75 50 35	
	Optimal working dilutions must be determined by the end user.			
SPECIES REACTIVITY:	Human			
FORMAT:	Raw Rabbit Serum			
PRESENTATION:	Raw Rabbit Serum with 0.1%	Sodium Azide		
STORAGE/HANDLING:	Maintain at 2-8°C for up to 12 months from date of receipt.			
REFERENCES:	Viaggi, C, et al (2006). Cytochrome P450 and Parkinson's disease: protective role neuronal CYP 2E1 from MPTP toxicity. J Neural Transm Suppl : 173-6.			ctive role of
	Cannady, Ellen A, <i>et al</i> (2003). Expression and activity of cytochromes P450 2E1, 2A, and 2B in the mouse ovary: the effect of 4-vinylcyclohexene and its diepoxide metabolite. <i>Toxicol Sci</i> 73 : 423-30.			

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.