



Quality Control Ranges
Human Cardiovascular Disease Magnetic Bead Panel 4
Catalog No. HCVD4MAG-67K
QC Catalog No. HCVMG4-6067-4 Lot No. HCVM4-113.213

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
sE-Selectin	Control 1	74 - 154	ng/mL
	Control 2	292 - 607	ng/mL
Follistatin	Control 1	1.0 - 2.1	ng/mL
	Control 2	3.9 - 8.0	ng/mL
Pecam-1	Control 1	6.2 - 13	ng/mL
	Control 2	21 - 44	ng/mL
Pentraxin-3	Control 1	3.3 - 6.9	ng/mL
	Control 2	13 - 26	ng/mL
Tissue Factor	Control 1	3.7 - 7.7	ng/mL
	Control 2	14 - 29	ng/mL
Thrombomodulin	Control 1	7.8 - 16	ng/mL
	Control 2	31 - 64	ng/mL
Troponin T	Control 1	0.08 - 0.17	ng/mL
	Control 2	0.30 - 0.63	ng/mL

Note: The Quality Control Ranges are generated with serum matrix provided in the kit.

HCVMG4-6067-4: HCVM4-113.213 (22 September 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)