



QC RANGE SHEET
Human Cardiovascular Disease Magnetic Bead Panel 4
Catalog # HCVD4MAG-67K
QC Catalog #: HCVMG4-6067-4

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
sE-Selectin	Control 1	HCVM4-111	72 – 150	ng/mL
	Control 2	HCVM4-211	279 – 579	ng/mL
Follistatin	Control 1	HCVM4-111	1.0 – 2.0	ng/mL
	Control 2	HCVM4-211	3.9 – 8.1	ng/mL
Pecam-1	Control 1	HCVM4-111	6.0 – 13	ng/mL
	Control 2	HCVM4-211	18 – 37	ng/mL
Pentraxin-3	Control 1	HCVM4-111	3.3 – 6.8	ng/mL
	Control 2	HCVM4-211	12 – 26	ng/mL
Tissue Factor	Control 1	HCVM4-111	3.8 – 7.8	ng/mL
	Control 2	HCVM4-211	14 – 29	ng/mL
Thrombomodulin	Control 1	HCVM4-111	8.0 – 17	ng/mL
	Control 2	HCVM4-211	30 – 63	ng/mL
Troponin T	Control 1	HCVM4-111	0.09 – 0.18	ng/mL
	Control 2	HCVM4-211	0.32 – 0.67	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

HCVMG4-6067-4: HCVM4-111.211 (27 February 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany. © 2013 EMD Millipore Corporation, Billerica, MA 01821 USA