



**Quality Control Ranges**  
**Milliplex PLEXpedition Magnetic Bead Screening Panel**  
**Catalog# HPLX1-115SP, HPLX1-115SP-PX, HPLX1-115SP-PX80**  
**Control Catalog # HPLX1-6010-07 Lot # PLX17-101 and PLX17-201**

<u>Bead</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
HMMP1	Control 1	54 - 113	pg/mL
	Control 2	215 - 446	pg/mL
HMMP2	Control 1	1227 - 2549	pg/mL
	Control 2	5278 - 10962	pg/mL
HMMP7	Control 1	224 - 466	pg/mL
	Control 2	935 - 1942	pg/mL
HMMP9	Control 1	405 - 841	pg/mL
	Control 2	1625 - 3376	pg/mL
HMMP3	Control 1	225 - 467	pg/mL
	Control 2	934 - 1940	pg/mL
HMMP10	Control 1	416 - 864	pg/mL
	Control 2	1682 - 3493	pg/mL
HMMP12	Control 1	520 - 1080	pg/mL
	Control 2	2101 - 4363	pg/mL

**HPLX1-6010-07: PLX17-101.201 (04 September 2025)**

**[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)**

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.