



SIGMA-ALDRICH

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

PHENYLEPHRINE HYDROCHLORIDE SIGMA REFERENCE STANDARD

Product No. **P8155**

Lot No: 111K8814

Storage Temperature: Room Temperature

Formula: $C_9H_{13}NO_2 \cdot HCl$

Formula Weight: 203.67

CAS No: 61-76-7

Synonyms: Benzenemethanol, 3-hydroxy- α -[(methylamino)methyl]-, hydrochloride (R)-, (-) Hydroxy- α -[(methylamino)methyl]benzyl alcohol hydrochloride.

Product Description



This product meets USP specifications (Current through USP 25), and is traceable to USP Reference Standard lot K

Store in sealed vial at room temperature. Protect from light. Dry a portion at 105°C for 2 hours before using. Use promptly. Discard unused material.

TEST	USP 25 SPECIFICATIONS	Result Sigma R.S. Lot 111K8814
Infrared absorption spectrum KBr dispersion (197K)	compares to standard	conforms to USP RS lot K
Identification B	a solution (1 in 100) responds to the test for Chloride (191)	passes
Melting range (741)	between 140° and 145°	142-143 °C
Specific rotation, 50 mg per mL, in water (781S)	between -42° and -47°	-46.1°
Loss on Drying (731)	not more than 1.0%	passes
Residue on ignition (281)	not more than 0.2%	passes
Sulfate (221)	not more than 0.20 %	passes
Limit of ketones	less than 1 in 2000	passes
Chromatographic purity (621)	Related impurities not more than 1.0%, with no single impurity at more than 0.5%	passes
Chloride content	between 17% and 17.7%	17.5%
Assay	not less than 97.5% and not more than 102.5%	99.5%

See reverse side for IR spectrum of Sigma Reference Standard.

PHENYLEPHRINE HYDROCHLORIDE
SIGMA REFERENCE STANDARD
Product No. **P8155**
Lot 111K8814



11/2001

For R&D use only. Not for drug, household or other uses. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.