



**QUALITY CONTROL RANGES
CANINE C-PEPTIDE RIA
CATALOG # CCP-24HK**

**CALIBRATOR, LYOPHILIZED
CATALOG #8024-K**

CALIBRATOR	CANINE C-PEPTIDE VALUE
LOT NUMBER: CCP80-17K	22 ng/mL

The exact value of lyophilized calibrator is shown above.
Reconstitute with **2 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED
CATALOG #6000L-K**

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGES	UNITS
QC 1	UV-133	0.92 – 1.9	ng/mL
QC 2	UV-233	3.8 – 7.8	ng/mL