

02336 PALCAM Listeria Selective Supplement according to Van Netten et al.

An antimicrobial supplement recommended for the selective isolation and identification of *Listeria monocytogenes*.

Composition

(per vial, sufficient for 500 ml medium)

Ceftazidime	10.0 mg
Polymyxin B	50000 IU
Acriflavin Hydrochloride	2.5 mg

Directions

Rehydrate the contents of one vial aseptically with 5 ml sterile distilled water and aseptically add to 500 ml, molten, cooled (45-50 °C) PALCAM Listeria Selective Agar (Prod No 75977) / PALCAM Listeria Selective Enrichment Broth, Vegitone (Prod No 59859) / PALCAM Listeria Selective Agar, Vegitone (Prod No 15776). Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

