




## SIGMA QUALITY CONTROL TEST PROCEDURE

## Product Information

### Enzymatic Assay of MALIC ENZYME (E.C. 1.1.1.40)

#### PRINCIPLE:



Abbreviations used:

$\beta$ -NADP =  $\beta$ -Nicotinamide Adenine Dinucleotide Phosphate, Oxidized Form

$\beta$ -NADPH =  $\beta$ -Nicotinamide Adenine Dinucleotide Phosphate, Reduced Form

**CONDITIONS:** T = 25°C, pH = 7.4,  $A_{340\text{nm}}$ , Light path = 1 cm

**METHOD:** Continuous Spectrophotometric Rate Determination

#### REAGENTS:

- A. 100 mM Triethanolamine HCl Buffer, pH 7.4 at 25°C.  
(Prepare 100 ml in deionized water using Triethanolamine, Hydrochloride, Sigma Prod. No. T-1502. Adjust to pH 7.4 at 25°C with 1 M NaOH.)
- B. 100 mM L-Malic Acid Solution (Malic Acid)  
(Prepare 5 ml in deionized water using L(-)Malic Acid, Free Acid, Sigma Prod. No. M-1000.)
- C. 20 mM  $\beta$ -Nicotinamide Adenine Dinucleotide Phosphate, Oxidized Form, Solution (NADP)  
(Prepare 2 ml in deionized water using  $\beta$ -Nicotinamide Adenine Dinucleotide Phosphate, Sodium Salt, Sigma Prod. No. N-0505 or  $\beta$ -Nicotinamide Adenine Dinucleotide Phosphate, Sodium Salt, Sigma Stock No. 240-310.)
- D. 20 mM Manganese Chloride Solution ( $\text{MnCl}_2$ )  
(Prepare 25 ml in deionized water using Manganese Chloride, Tetrahydrate, Sigma Prod. No. M-3634.)
- E. Malic Enzyme Solution  
(Immediately before use, prepare a solution containing 0.25 - 0.50 unit/ml of Malic Enzyme in cold deionized water.)

**Enzymatic Assay of MALIC ENZYME  
(E.C. 1.1.1.40)**

**PROCEDURE:**

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	2.00	2.00
Reagent B (Malic Acid)	0.10	0.10
Reagent C (NADP)	0.05	0.05
Reagent D (MnCl <sub>2</sub> )	0.75	0.75

Mix by inversion and equilibrate to 25°C. Monitor the A<sub>340nm</sub> until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent E (Enzyme Solution)	0.10	-----
Deionized Water	-----	0.10

Immediately mix by inversion and monitor the increase in A<sub>340nm</sub> for approximately 5-10 minutes. Obtain the ΔA<sub>340nm</sub>/minute using the maximum linear rate for both the Test and Blank.

**CALCULATIONS:**

$$\text{Units/ml enzyme} = \frac{(\Delta A_{340\text{nm}}/\text{min Test} - \Delta A_{340\text{nm}}/\text{min Blank})(3)(\text{df})}{(6.22)(0.1)}$$

3 = Total volume (in milliliters) of assay

df = Dilution factor

6.22 = Millimolar extinction coefficient of β-NADPH at 340 nm

0.1 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

**UNIT DEFINITION:**

One unit will convert 1.0 μmole of L-malate and NADP to pyruvate, CO<sub>2</sub> and NADPH per minute at pH 7.4 at 25°C.

## Enzymatic Assay of MALIC ENZYME (E.C. 1.1.1.40)

### FINAL ASSAY CONCENTRATIONS:

In a 3.00 ml reaction mix, the final concentrations are 67 mM triethanolamine, 3.3 mM L-malic acid, 0.3 mM  $\beta$ -nicotinamide adenine dinucleotide phosphate, 5.0 mM manganese chloride and 0.025 - 0.050 unit malic enzyme.

### REFERENCE:

Geer, B.W., Krochko, D., Oliver, M.J., Walker, V.K. and Williamson, J.H. (1980) *Comp. Biochem. Physiol.* 65B, 25-34

### NOTES:

1. This assay is based on the cited reference.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

**Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.**