



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

10X CAPILLARY ELECTROPHORESIS RUNNING BUFFER

Product Number **B4930**
Storage at Room Temperature

Product Description

10X Capillary Electrophoresis Buffer for automated sequencing has been optimized for use with the ABI Prism® 3700 DNA Analyzer and Molecular Dynamics MegaBACE® capillary electrophoresis instruments. The product has been tested under stringent quality control standards to ensure proper conductivity and pH. It has also been application tested to provide consistency and quality performance in every lot.

Reagent

This product is provided as a 10X concentrated solution in plastic light-protective containers.

Precautions and Disclaimer

Capillary Electrophoresis Buffer is for laboratory use only. Not for drug, household, or other uses.

Preparation Instructions

Add 1 volume of the 10X capillary electrophoresis buffer solution to 9 volumes of deionized water. Cover the top of the container and mix the 1X solution by inversion.

Storage/Stability

Protect from light.
10X Capillary Electrophoresis Buffer may be stored at room temperature. 1X dilutions may be stored for up to 6 months at 2 °C to 8 °C.

3/2007 ALF/LSP

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.