


Eppendorf Certificate

Lithium Battery Test Summary / UN38.3

Article no.	2008561
1. Product Manufacturer	PMBL Limited
2. Battery Manufacturer	PMBL Limited
	Unit D Leach Road Chard Business Park Chard - Somerset TA 20 IFA
	Tel. +44 1460 66366 510488433@qq.com
3. UN 38.3 Test Lab	Pony Testing International Group Shenzhen Co. Ltd.
	1/F., Building A2, Jun Feng ZhongCheng ZhiZao Innovation Park, HePing Peace Community, FuHai Road
Address	BaoAn District, Shenzhen, Guangdong, China www.ponytest.com
4. Test Report Number	MNIYR67X94332721
5. Date of Test Report	19.08.2019
6. Description of Cell or Battery	Battery
Model	2008561
	Lithium Ion Battery Pack
Mass	48.286 g
Lithium content	
Nom. Voltage	3.7 V
Nom. Capacity	2.6 Ah
7. UN 38.3 Test Performed and Successfully passed	
T1 : Altitude simulation	Pass
T2 : Thermal Test	Pass
T3 : Vibration	Pass
T4 : Shock	Pass
T5 : External short circuit	Pass
T6 : Impact / Crush	Pass
T7 : Overcharge	Pass
T8 : Forced discharge	Pass
8. Assambled Battery Testing Requirements	
9. Edition of UN Manual of Tests an Criteria Used	ST/SG/AC.10/11/ Rev.6/Amend.1

Date:

21.04.2020



Simon Reiss
Product Life Cycle Manager

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf AG, Germany.
All rights reserved, incl. graphics and images. Copyright ©2020 by Eppendorf AG.

www.eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified