

17122 *Clostridium difficile* Supplement

An antibiotic supplement for the selective isolation of *Clostridium difficile*.

Composition

(per vial, sufficient for 500 ml medium)

D-Cycloserine 250 mg

Cefoxitin 8 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix vigorously (vortex) for 1 minute and aseptically add along with 7% (v/v) defibrinated horse blood to 500 ml sterile, molten *Clostridium Difficile* Agar Base (Cat. No. 17145). Mix well and pour into sterile petri plates. Sheep blood may be used in place of horse blood, but some strains of the organism will show a slightly reduced growth.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

