

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

3-Isobutyl-1-methylxanthine

BioUltra

Catalog Number 17018

CAS RN 28822-58-4

Synonyms: IBMX, 1-Methyl-3-isobutylxanthine, 3,7-Dihydro-1-methyl-3-(2-methylpropyl)-1H-purine-2,6dione, 3-Isobutyl-1-methyl-2,6(1H,3H)-purinedione

Molecular Formula: C₁₀H₁₄N₄O₂

Molecular Weight: 222.2

Product Description

Methylxanthines such as caffeine and theophylline inhibit adenosine 3',5'-cyclic monophosphate phosphodiesterase (cAMP PDE).1 3-Isobutyl-1methylxanthine (IBMX) has been shown to be a potent inhibitor of cAMP PDE, significantly more effective than theophylline.²⁻⁶ IBMX inhibits cyclic nucleotide PDE with subsequent inhibition of cyclic nucleotide hydrolysis, which results in accumulation of cyclic AMP and guanosine 3',5'-cyclic monophosphate. 7,8 In a study of cyclic AMP and insulin release by islets of Langerhans, IBMX at 1 mM caused a marked increase in the intracellular concentration of cyclic AMP in the presence of glucose.4

Caffeine and theophylline stimulate lipolysis in fat cells. IBMX was 20-fold more effective than theophylline, used at 0.05 mM.3 IBMX has been shown to promote the conversion of fibroblast cells into adipose cells. 9 As a PDE inhibitor, IBMX was shown to inhibit the growth of carcinoma cells both in vivo and in vitro in mice. 10

Various publications have cited use of this specific IBMX product in different systems, such as:

- Differentiation of cultured 3T3-L1 pre-adipocytes¹¹
- Isolating germinal vesical-stage oocytes¹
- Germinal vesicle transfer studies, in the micromanipulation medium¹³

Trace elemental analyses have been performed on the BioUltra 3-Isobutyl-1-methylxanthine. The Certificate of Analysis (CofA) provides lot-specific results. BioUltra 3-Isobutyl-1-methylxanthine is for applications which require tight control of elemental content.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses.

Preparation Instructions

IBMX is soluble in different organic solvents, as follows:

- Warm methanol at 50 mg/mL
- Ethanol at 10 mg/mL. Dissolves at 25 mg/mL only with sonication.
- Various publications report preparation in DMSO at 50 mg/mL¹³ and at 110 mg/mL.¹⁵

A 10 mM aqueous solution can be prepared by warming in a boiling water bath. 16

Storage/Stability

One publication reports storage of IBMX solutions in ethanol at 0.125 mg/mL at 4 °C for at least three months. 14 Aqueous solutions can be frozen in aliquots, then thawed for use by heating in a boiling water bath. These aliquots are stable for several months. 16 Solutions in DMSO may be stored at -20 °C. 15

References

- Chasin, M, and Harris, D.N., Adv. Cyclic Nucleotide Res., 7, 225-228 (1976).
- Data for Biochemical Research, 3rd ed. (Dawson, R.M.C. et al., eds.). Oxford University Press (Oxford, UK), pp. 326-327 (1986).
- 3. Beavo, J.A. et al., Mol. Pharmacol., 6(6), 597-603 (1970).
- 4. Montague, W., and Cook, J.R., *Biochem. J.*, **122(1)**, 115-120 (1971).
- Peytreman, A. et al., Endocrinology, 92(2), 525-530 (1973).
- 6. Ashcroft, S.J.H. et al., FEBS Lett., 20(3), 263-266 (1972).
- Klotz, U. et al., Naunyn-Schmiedebergs Archives Pharmacol., 296(2), 187-190 (1977).
- 8. Spaulding, S.W. and Burrow, G.N., Biochem. Biophys. Res. Commun., 59(1), 386-391 (1974).
- 9. Russell, T.R., Proc. Natl. Acad. Sci. USA, 76(9), 4451-4454 (1979).
- 10. Janik, P. et al., Cancer Res., 40(6), 1950-1954 (1980).
- 11. Ishibashi, K. et al., "Effect of Trans Fatty Acid on Insulin Responsiveness and Fatty Acid Composition of Lipid Species of 3T3-L1 Adipocytes", in Adipose Tissue (L. Szablewski, ed.). Intechopen, DOI: 10.5772/intechopen.76646 (May 30, 2018).

- 12. Gassler, J. et al., Methods Cell Biol., **144**, 389-407 (2018).
- 13. Liu, L., and Keefe, D.L., *Methods Mol. Biol.*, **371**, 191-207 (2007).
- 14. Schwertner, H.A. et al., Anal. Chem., **48(13)**, 1875-1878 (1976).
- 15. Eckel, J., *The Cellular Secretome and Organelle Crosstalk.* Academic Press/Elsevier, p. 167 (2018).
- 16. Salomon, Y., *Methods Enzymol.*, **195**, 22-28 (1991).

GCY,CKV,SM 05/20-1