



**QUALITY CONTROL RANGES
PYY (TOTAL) AND PYY (3-36) RIA KIT
CATALOG NO. PYYT-66HK AND PYY-67HK**

**STANDARD, LYOPHILIZED
CATALOG NO. 8066-K**

STANDARD	PYY VALUE
Lot Number. 19T66	1303 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **2 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED
CATALOG NO. 6066-K**

Reconstitute each vial with **1 mL** distilled or deionized water.

Kit	LOT NUMBER	EXPECTED RANGE	UNITS
PYY Total	HPYY-119	47 – 98	pg/mL
	HPYY-219	155 – 323	pg/mL
PYY 3-36	HPYY-119	53 – 111	pg/mL
	HPYY-219	154 – 320	pg/mL