

## Product Information

### 18284 Legionella Selective Supplement

An antimicrobial supplement recommended for the selective isolation of Legionella species

#### Composition

(per vial, sufficient for 500 ml medium):

Colistin sulphate	7500 U
Vancomycin	2.5 mg
Trimezhoprim	1.25 mg
Amphotericin B	1.25 mg

#### Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten Legionella Agar Base (e.g. Feeley Gorman Agar; Cat. No. 05598). Mix well and pour into sterile petri plates.