



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Human Neurodegenerative Disease Magnetic Panel 4 Catalog # HNDG4MAG-36K

Control Catalog # HNDG4-6036-4

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
S100B	Control I	HNDG4-109	58 - 121	pg/mL
	Control II	HNDG4-209	586 - 1218	pg/mL
A β 40	Control I	HNDG4-109	13 - 28	pg/mL
	Control II	HNDG4-209	165 - 343	pg/mL
A β 42	Control I	HNDG4-109	110 - 229	pg/mL
	Control II	HNDG4-209	1140 - 2367	pg/mL
sRAGE	Control I	HNDG4-109	81 - 168	pg/mL
	Control II	HNDG4-209	864 - 1795	pg/mL
GDNF	Control I	HNDG4-109	27 - 57	pg/mL
	Control II	HNDG4-209	284 - 590	pg/mL

HNDG4-6036-4: HNDG4-109.209 (18 November 2019)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.