



PURIFIED HUMAN APOLIPOPROTEIN A-II

CATALOG NUMBER:	ALP20
LOT NUMBER:	
QUANTITY:	500 µg
CONCENTRATION:	2.94 mg/mL
SOURCE:	Normal human plasma, fresh, non-frozen. Tested negative for HIV-1, HIV-2, HCV and HBc antibodies and HBsAg.
PURITY:	≥ 95% by SDS-PAGE. Single band on polyacrylamide electrophoresis.
MOLECULAR WEIGHT:	17,400
FORMAT:	Liquid in 10 mM NH ₄ HCO ₃ , pH 7.4. Non-sterile, no preservatives.
STORAGE/HANDLING:	Maintain at -20°C to -80°C at a concentration of ≥ 1 mg/mL for up to six months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com