

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

Diazoxide

Product Number **D9035** Store at Room Temperature

Product Description

Molecular Formula: C₈H₇CIN₂O₂S

Molecular Weight: 230.7 CAS Number: 364-98-7 Melting Point: 330-331 $^{\circ}$ C¹ λ_{max} : 268 nm (methanol)

 λ_{max} : 268 nm (methanol) Extinction coefficient: $E^{mM} = 11.3$ (methanol)¹

Synonyms: 7-chloro-3-methyl-

2H-1,2,4-benzothiadiazine 1,1-dioxide; 3-methyl-7-chloro-1,2,4-benzothiadiazine 1,1-dioxide¹

Diazoxide is a selective and ATP-sensitive K^+ channel activator in both vascular smooth muscle and pancreatic β -cells. It also inhibits insulin secretion by pancreatic β -cells. The use of diazoxide and other compounds to probe the ATP-sensitive K^+ channels in the TE671 human medulloblastoma cell line has been reviewed.

The effects of diazoxide on mitochondrial membrane potential, Ca^{2^+} transport, oxygen consumption and ATP generation in mouse pancreatic β -cells and rat liver mitochondria have been probed. The role of ADP in altering the responsiveness of cardiac potassium ATP channels to diazoxide (300 μ M) has been investigated. The protective effect of diazoxide in indomethacin-induced gastric injury in rats has been studied. A study of cultured human islets with diazoxide to probe insulin release has been published.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in dimethylformamide (50 mg/ml), yielding a clear, colorless to light yellow solution. It is also soluble in DMSO (10 mg/ml) and methanol (3.2 mg/ml). It is soluble in alcohol and alkaline solutions, and insoluble in water.¹

Storage/Stability

Solutions of this product are expected to be stable at 4 °C for 24-48 hours and may be frozen for longer storage.

References

- 1. The Merck Index, 12th ed., Entry# 3051.
- 2. Martindale The Extra Pharmacopoeia, 31st ed., Reynolds, J. E. F., ed., Royal Pharmaceutical Society (London, UK: 1996), pp. 849-850.
- Miller, T. R., et al., Pharmacological and molecular characterization of ATP-sensitive K⁺ channels in the TE671 human medulloblastoma cell line. Eur. J. Pharmacol., 370(2), 179-185 (1999).
- 4. Grimmsmann, T., and Rustenbeck, I., Direct effects of diazoxide on mitochondria in pancreatic B-cells and on isolated liver mitochondria. Br. J. Pharmacol., **123(5)**, 781-788 (1998).
- D'hahan, N., et al., Pharmacological plasticity of cardiac ATP-sensitive potassium channels toward diazoxide revealed by ADP. Proc. Nat. Acad. Sci. USA, 96, 12162-12167 (1999).
- Akar, F., et al., Protective effect of cromakalim and diazoxide, and proulcerogenic effect of glibenclamide on indomethacin-induced gastric injury. Eur. J. Pharmacol., 374, 461-470 (1999).
- Song, S. H., et al., Diazoxide attenuates glucoseinduced defects in first-phase insulin release and pulsatile insulin secretion in human islets. Endocrinology, 144(8), 3399-3405 (2003).

GCY/NSB 12/03