

67432 **Campylobacter Millichrome™ plus Supplement for 5L**

Selective supplement for detection of thermotolerant *Campylobacter*.

Composition:

(1 vial is sufficient for 5 L of final media)

Chromogenic and selective mix 1.05 g

Directions

- In a transparent vessel, add 210 mg of Campylobacter Millichrome™ plus Supplement in 10 ml of purified water.
- Swirl well until complete dissolution.
- Filter to sterilize at 0.45 µm.
- Add the 10 ml of the supplement solution to 1 liter of the melted Campylobacter Millichrome™ plus Agar base at 45-50 °C.
- Swirl or stir gently to homogenize.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

