

## Product Information

# 92301 *Listeria* mono Selective Supplement I

A selective supplement recommended by ISO Committee for the isolation of *Listeria* species.

### Composition

(per vial, sufficient for 500 ml medium)

Polymyxin B sulphate 38350 IU

### Directions

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add it to 460 ml of sterile, molten L. mono Differential Agar Base (Cat. No. 77408) along with sterile contents of L. mono Enrichment Supplement I (Cat. No. 03708) and sterile rehydrated contents of L. mono Selective Supplement II (Cat. No. 91603) or add in 470 ml of sterile, molten L. mono confirmatory Agar Base (Cat. No. 92302) along with sterile content of L. mono Enrichment Supplement II (Cat. No. 15895) and sterile rehydrated content of L. mono Selective Supplement II (Cat. No. 91603). Mix well and pour in sterile petri plates.