



QUALITY CONTROL RANGES
MILLIPLEx[®] MAP Human Adipokine Magnetic Bead
Panel 2
Kit Catalog # HADK2MAG-61K
Control Catalog # HADK2-6061-2

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
NGF	Control I	HADK2-111	74 - 154	pg/mL
	Control II	HADK2-211	317 - 659	pg/mL
IL-6	Control I	HADK2-111	102 - 212	pg/mL
	Control II	HADK2-211	495 - 1028	pg/mL
Insulin	Control I	HADK2-111	1035 - 2149	pg/mL
	Control II	HADK2-211	5078 - 10547	pg/mL
Leptin	Control I	HADK2-111	5373 - 11159	pg/mL
	Control II	HADK2-211	24089 - 50031	pg/mL
IL-8	Control I	HADK2-111	70 - 146	pg/mL
	Control II	HADK2-211	305 - 633	pg/mL
MCP-1	Control I	HADK2-111	148 - 308	pg/mL
	Control II	HADK2-211	679 - 1411	pg/mL
TNF α	Control I	HADK2-111	67 - 139	pg/mL
	Control II	HADK2-211	330 - 686	pg/mL
IL-1 β	Control I	HADK2-111	137 - 284	pg/mL
	Control II	HADK2-211	655 - 1361	pg/mL

HADK2-6061-2: HADK2-111.211 (27 April 2023)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA