

Quality Control Ranges

MILLIPLEX® Human Adipokine Magnetic Bead Panel 2

Catalog Number HADK2MAG-61K

QC Catalog Number/ Lot Number HADK2-6061-2/ HADK2-114.214

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
NGF	Control I	HADK2-114	73 - 151	pg/mL
	Control II	HADK2-214	327 - 678	pg/mL
IL-6	Control I	HADK2-114	87 - 181	pg/mL
	Control II	HADK2-214	433 - 898	pg/mL
Insulin	Control I	HADK2-114	933 - 1937	pg/mL
	Control II	HADK2-214	4670 - 9698	pg/mL
Leptin	Control I	HADK2-114	5028 - 10443	pg/mL
	Control II	HADK2-214	23102 - 47980	pg/mL
IL-8	Control I	HADK2-114	73 - 151	pg/mL
	Control II	HADK2-214	337 - 700	pg/mL
MCP-1	Control I	HADK2-114	139 - 288	pg/mL
	Control II	HADK2-214	617 - 1282	pg/mL
TNFα	Control I	HADK2-114	65 - 136	pg/mL
	Control II	HADK2-214	294 - 610	pg/mL
IL-1β	Control I	HADK2-114	128 - 266	pg/mL
	Control II	HADK2-214	624 - 1296	pg/mL

The life science business of Merck KGaA, Darmstadt, Germany operates as

MilliporeSigma in the U.S. and Canada.

The vibrant M, Milliplex, Millipore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Template 20300637, Ver 4.0

Cat#: HADK2-6061-2, Lot#: HADK2-6061-2 (Rev 03 FEB 2026)

