

Product Information

17775 Campylobacter Selective Supplement

An antibiotic supplement recommended for the selective isolation of *Campylobacter* species from human and animal derived clinical material and from contaminated food, water etc.. The *Campylobacter* selective supplement largely inhibits the accompanying microbial flora.

Composition

(per vial, sufficient for 1000 ml medium) Rifampicin 25.00 mg, Cefsulodin 6.25 mg, Polymyxin B sulfhate 20.000 Units; One package contains 5 vial.

Directions

Rehydrate the content of 1 vial aseptically with 10 ml of 50% ethanol. Mix well and aseptically add to 1000 ml of sterile, Wesley Broth (Cat. No. 17224) or *Campylobacter* Selective Agar (Cat. No. 21378). Mix well and dispense as desired.