

Product Information

CHLORAMPHENICOL SIGMA REFERENCE STANDARD

 Product No. **C6455**

Lot No: 55K8804

Storage Temperature: 2 – 8°

 Formula: C₁₁H₁₂Cl₂N₂O₅

Formula Weight: 323.13

CAS No: 56-75-7

 Synonyms: Acetamide, 2,2-dichloro-*N*-[2-hydroxy-1-(hydroxymethyl)-2-(4-nitrophenyl)ethyl]-, [*R*-(*R**,*R**)]-, *D-threo*-(-)-2,2-Dichloro-*N*-[β-hydroxy-α-(hydroxymethyl)-*p*-nitrophenethyl]acetamide.

Product Description


This product meets USP specifications (Current through USP 28, 1st Supplement), and is traceable to USP Reference Standard lot N1C074

Do not dry before using. Keep container tightly closed. Store in a cold dry place.

TEST	USP 28 SPECIFICATIONS	Result Sigma R.S. Lot 55K8804
Infrared absorption spectrum KBr dispersion (197K)	compares to standard	conforms to USP RS lot N1C074
Identification B high pressure liquid chromatography (621)	retention time of the major peak of the chromatogram in the <i>Assay preparation</i> corresponds to the major peak in the <i>Standard preparation</i> as obtained in the <i>Assay</i>	passes
Melting range (741)	between 149° and 153°C	149.5-151.5°C
Specific rotation (781S) 50 mg, undried, per mL, in dehydrated alcohol	between +17.0° and +20.0°	+19.7°
Crystallinity (695)	meets requirements	passes
pH (791)	between 4.5 and 7.5, in an aqueous suspension containing 25 mg	5.7
Chromatographic purity	no impurity greater than 1% and the sum of all impurities does not exceed 2%	passes
Assay high performance liquid chromatography (621)	not less than 97.0 percent and not more than 103.0. percent	100.0%


See reverse side for IR spectrum of Sigma Reference Standard.

CHLORAMPHENICOL
SIGMA REFERENCE STANDARD
Product No. **C6455**
Lot 55K8804



INFRARED SPECTRUM

KBr dispersion



4/2005

Sigma brand products are sold through Sigma-Aldrich, Inc.
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.
Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.
Please see reverse side of the invoice or packing slip.