

Product Information

Anti-CXCL11
antibody produced in chicken,
affinity isolated antibody

Catalog Number **GW22690**
Formerly listed as GenWay Catalog Number 15-288-22690,
Small inducible cytokine B11 Antibody.

Storage Temperature Store at -20 °C

Synonyms: Small inducible cytokine B11, CXCL11; Interferon-inducible T-cell alpha chemoattractant; I-TAC; Interferon-gamma-inducible protein 9; IP-9; H174; Beta-R1

Product Description

Chemotactic for interleukin-activated T-cells but not unstimulated T-cells. neutrophils or monocytes. Induces calcium release in activated T-cells. Binds to CXCR3. May play an important role in CNS diseases which involve T-cell recruitment. May play a role in skin immune responses.

NCBI Accession number: NP_005400.1
Swiss Prot Accession number: O14625

Gene Information: Human .. CXCL11 (6373)

Immunogen: Recombinant protein Small inducible cytokine B11

Immunogen Sequence: GI # 4885589, sequence 22 - 85

The product is a clear, colorless solution in phosphate buffered saline, pH 7.2, containing 0.02% sodium azide.

Species Reactivity: Human

Tested Applications: WB

Recommended Dilutions: Recommended starting dilution for Western blot analysis is 1:500, for tissue or cell staining 1:200.

Note: Optimal concentrations and conditions for each application should be determined by the user.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one week. For extended storage, store in -20 °C freezer in working aliquots. Repeated freezing and thawing, or storage in "frostfree" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

LPG,MAM12/09-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.