



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Rat Kidney Toxicity Magnetic Bead Panel 2
Kit Catalog # RKT2MAG-37K
Control Catalog # RKT2MG2-6037-2
Lot# RKT2-111 and RKT2-211

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Albumin	Control I	68 - 142	ng/mL
	Control II	343 - 712	ng/mL
AGP	Control I	8.7 - 18	ng/mL
	Control II	47 - 99	ng/mL
β-2-Microglobulin	Control I	3.0 – 6.2	ng/mL
	Control II	16 - 33	ng/mL
Cystatin C	Control I	0.06 - 0.13	ng/mL
	Control II	0.29 - 0.61	ng/mL
EGF	Control I	0.02 - 0.04	ng/mL
	Control II	0.08 - 0.17	ng/mL
Lipocalin-2/NGAL	Control I	0.08 - 0.17	ng/mL
	Control II	0.40 - 0.83	ng/mL

RKT2-6037-2: RKT2-111.211 (07 July 2023)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA