

**MOUSE ANTI-PARAINFLUENZA III
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB855-3

LOT NUMBER:

QUANTITY: 100 uL

SPECIFICITY: Specific for the HN protein, common to all parainfluenza III strains. This monoclonal exhibits no cross-reactivity to Parainfluenza I or II.

ISOTYPE: Blend of MAB855-1 and MAB855-2.

CLONE NAME: Blend of MAB855-1 and MAB855-2.

APPLICATIONS: General guidelines:
ELISA 1:800+
Indirect Immunofluorescence 1:100 - 200+
(fresh frozen tissue sections)
Final working dilutions must be determined by end user.

FORMAT: Ascites fluid, wet liquid, frozen.

STORAGE/HANDLING: Maintain at -20°C in aliquots for up to 12 months.
Short term storage (1-3 weeks) maintain at +4°C, the addition of 0.1% sodium azide is suggested.
Do not store in diluted format.
Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com